STORZ
KARL STORZ — ENDOSKOPE

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**Instructions
manual**

**Manual de
instrucciones**


**Endo-Koagulator
Modelle
26027 B/C/D**

**Endo-Coagulator
Models
26027 B/C/D**


**Endocoagulator
Modelos
26027 B/C/D**

- 1
- 2
- 3
- 4
- 5
- 6

- 7
- 8
- 9
- 10
- 11
- 12

- 13

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**Operating elements,
displays, connections
and their functions**

- ① Power switch, illuminates when on
- ② Jack for connecting instrument cord
- ③ Jack for connecting footswitch
- ④ Coagulation temperature display
- ⑤ Control to set coagulation temperature
- ⑥ Control to set coagulation time
- ⑦ Volume control for audible signal
- ⑧ Jack for connecting potential-equalization cable
- ⑨ Spare fuse compartment
- ⑩ Power fuse holders:
for 26027 B/D: 2 × T 1 A (SB)
for 26027 C: 2 × T 2 A (SB)
- ⑪ Power cord connector
- ⑫ Nameplate
- ⑬ Footswitch with connecting tube


Explanation of symbols:

-  Read the instructions carefully before
operating the unit!
-  Jack for potential-equalization cable

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Pictures of the instrument (A1)
Operating elements, displays, connections
and their functions (A2)

Installation and operating instructions

Unpacking (2)
Basic equipment (2)
Setting up (2)
Starting (3)
Function test (7)
How endocoagulation
works (8)
Range of applications of
endocoagulation (10)

Description of the instrument

Exchanging fuses (12)
Cleaning the instrument (13)
Maintenance (13)
Repairs (14)
Responsibility (14)
Warranty (14)

Technical description

Trouble shooting (15)
Specifications (16)
Classification (16)
Test certificate (16)
Technical information (17)

Spare parts, Recommended accessories

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Unpacking the instrument

Carefully remove the Endo-Coagulator and the accessories from the packing. Check the merchandise to make sure it is complete and that there are no damages. Notify the manufacturer or supplier if the delivered merchandise is incomplete or damaged in any way. If possible, keep the original packing material. It will be useful if you have to transport the instrument.

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
Installation

Important: if the Endo-Coagulator is operated in medical rooms, the electrical installation of the room must comply with the National Electrical Code 0107. The instrument is not designed for use in hazardous locations. If explosive anesthetic gases are used the instrument may not be operated in the danger range as displayed on the instrument.

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
Stand the instrument on a flat surface.

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
The instrument in accordance with DIN 42801 is equipped with a jack ⑧ for the potential-equalization cable.

Have the instrument grounded by a qualified technician.

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
Starting the Endo-Coagulator

Before plugging in the instrument, make sure that the voltage on the nameplate ⑫ corresponds to the voltage of the local power line!

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
Connect power cord. Insert the power cord into connector ⑩ as far as it will go.

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
Connecting the instrument to the power line should only be carried out outside potentially explosive locations.

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
Plug the connecting tube of the footswitch ⑬ into jack ③.

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
Plug the connecting cord of the coagulation instrument into jack ②.

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
Note: If the coagulation instrument has not been connected, a 2-tone warning signal will sound. The warning signal should stop as soon as the instrument has been connected. If not, this indicates a fault in the instrument or the connecting cord.

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
The coagulation temperature display ④ will now show the temperature of the instrument tip (normally room temperature).

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The unit is now ready for operation.
A beep sounds at 10-second intervals to indicate
that the unit is ready for operation.


The volume of the audible signal can be set with the
control knob ⑦ on the rear panel of the unit.

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
Set the desired coagulation temperature with the
control knob ⑤.

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
Set the desired coagulation time with the control
knob ⑥.

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
After the footswitch ⑬ has been activated, the
coagulation instrument heats up to the preset
temperature.

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The increase in temperature can be read on the temperature display (1).


In addition, the increase in temperature is accompanied by an audible signal which rises in pitch as the temperature climbs. This allows the surgeon to keep his eyes on the operative field.

Caution: To prevent tissue from accidentally being destroyed, the use of destructive heat should only be carried out under visual control.

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
After the preset coagulation time has been reached, the heating shuts off automatically. The pitch of the audible signal then decreases.

Note: Generally, the instrument should only be removed from the abdomen after the audible signal has stopped, thus indicating the temperature of the instrument tip has fallen to below 60° C.

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To interrupt the heating process before the preset coagulation time has been reached, release the footswitch (3).

The descending pitch of the audible signal indicates that the heating process has been shut off.

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
Caution: After you have switched off the unit, you must wait about 30 seconds before turning it back on.

If you detect a change in the sound of the audible signal, switch off the unit and have it inspected by a qualified technician.

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Function test


Note: Carry out a function test before each use.

Connect footswitch (3) and coagulation instrument.
Switch on power switch (1).
Make sure that power switch illuminates.
Set the desired temperature with control knob (5).

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
Set the coagulation time with control knob (6).

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Immerse the coagulation instrument in a fluid solution (make sure that it is sterile). Activate the footswitch and observe the temperature display (4).


Check that the temperature climbs to the value you have set and that it decreases again after the preset time has elapsed.

This should be accompanied by an audible signal whose pitch rises as the temperature increases and falls as the temperature decreases.

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Repeat procedure.

As soon as the preset temperature has been reached, release the footswitch.

This should interrupt the heating process, which is then shown by the temperature display (4) and by the falling pitch of the audible signal.

If the function test detects a malfunction, switch off the unit and have it inspected by a service technician.

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How endocoagulation works

After many years of experimental tests, we now know that the best method of thermal hemostasis is coagulation of protein at 100° C. This effect can be seen clearly when boiling an egg. The new endocoagulation technique is based upon this principle, however. The human body does not come into contact with the electric current. The tissue is heated passively by thermal radiation and conduction, in the same way as on an electric hot-plate.

With the endocoagulation method, it suffices to heat the instruments (crocodile forceps, point coagulator and/or myoma enucleator) to a temperature of between 100° C and 120° C in order to achieve adequate tissue coagulation. Temperatures higher than this must be avoided. If the tissue is heated too rapidly, this will result in desiccation and the tissue will stick to the instrument (similar to burnt food in a saucepan!) One should always endeavor to "cook" the tissue with as low a heating temperature as possible, since coagulation of the protein is sufficient. Desiccation of the protein is superfluous and will interfere with wound healing later on.

Why do endocoagulated wounds remain free of adhesions? Histochemical investigations have demonstrated that protein which is coagulated at a temperature of up to 120° C does not become sequestered. Protein which is coagulated in this way remains on the healthy tissue as a sort of protective dressing. It serves to promote immigration of histiocytes and fibrocytes. On the other hand, healthy peritoneum cannot grow firmly on this biologically dead, coagulated protein. This explains why even large coagulated areas, such as on the fundus of the uterus, remain free of adhesions after myoma enucleation.

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Medical range of applications of endocoagulation.

1. Hemostasis in the fallopian tubes for subsequent bloodless resection for laparoscopic sterilization
2. Coagulation of adhesions
3. Coagulation of endometriotic foci
4. Coagulation of the edges of wounds in the ovary
5. Hemostasis in the area of sactosalpinges for bloodless cleavage
6. Virtually bloodless enucleation of subserous or intramural myoma tubercles
7. Detaching large adhesions in the ovaries and fallopian tubes.

The endocoagulation technique should not be used for hemostasis of highly vascularized omental and intestinal adhesions. In such cases it is safer to use loop or endoligation methods.

Endocoagulation assists endoscopic abdominal surgery, but does not take the place of ligating and suturing in the omentum and intestinal area.

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**Range of applications
of endocoagulation**

The empirical values for time and temperature (e. g. 20 seconds and 110°C) are used for routine operations such as tubal coagulation and coagulation of endometritic foci, etc. When the pneumatic footswitch is activated (the footswitch does not operate with electric current!) the therapy instrument, e. g. the crocodile forceps, is heated to the set temperature, e. g. 110°C.

After the set coagulation time, the heat shuts off automatically. This can be recognized by the drop in signal frequency.

When temperature of the coagulation instrument drops below 60°C, the sound signal shuts off automatically. The coagulation instrument can then be removed from the abdomen.

The forceps can be cooled more rapidly by immersing it in the fluid of the Douglas pouch. Heating can be stopped at any time by releasing the footswitch.

If the tissue adheres slightly to the forceps, either coagulation has been carried out at too high a temperature or the tissue has been subjected to the high temperature for too long. Both result in the desiccation of the tissue adjoined to the forceps. At temperatures above 130°C the tissue is transformed to a sort of adhesive, a process similar to that of burnt milk sticking in a saucepan.

This effect can be avoided completely with experience and by using lower coagulation temperatures. Adhering tissue can be removed easily and at any time from the coagulation instrument with the trocar or endoscope. If bleeding occurs, coagulation is repeated using a lower temperature.

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When the endocoagulation technique is used, the human body does not come into contact with electric current. The heat of the coagulation instrument is transferred only by thermal conduction and thermal radiation. This means that only tissue which is in direct contact with the coagulation instrument can be heated. Transfer of the heat to distant points in the organism does not occur.

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For coagulation of endometriotic foci, the point coagulator or myoma enucleator must be lowered onto the endometriotic focus for a few seconds. This causes coagulation down to a depth of 1-2 mm. Healthy peritoneum turns white and an endometriotic focus turns brown (hemosiderin effect). This "thermo-color effect" is specific to endometriotic foci. Since thermal radiation and conduction occurs slowly, even endometriotic foci on the roof of the bladder can be coagulated at a temperature of 100-110°C.

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This method should not be used for endometriotic foci of the intestines. Such endometriotic foci frequently extend deep into the intestinal muscles. Coagulating them is thus very dangerous and must not be carried out.

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Large-area peritoneal wounds in the area of the parietal peritoneum (after extensive adhesiolysis) are coagulated by using the myoma enucleator or crocodile forceps in the fashion of a smoothing iron for adhesion prophylaxis.

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Vascularized adhesions can be severed neatly with hooked scissors after adequate coagulation (e. g. 110°C and 20 seconds). The following rule must always be observed:

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1. Grip the tissue to be coagulated.
2. Check that you have gripped only the tissue to be coagulated (beware: intestine etc.) by tugging the instrument.
3. Operate the pneumatic footswitch for coagulating the tissue and check that the sound signal rises in pitch and then falls after the coagulation is automatically stopped when you remove your foot from the footswitch.
4. Remove the crocodile forceps etc.
5. Inspect the coagulated tissue.
6. Sever the tissue neatly with hooked scissors.

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Cleaning the instrument

Caution: Always pull out power plug before cleaning!

The instrument surface should be cleaned with a mild cleaning agent or disinfectant. Do not allow any cleaning liquid to get into the instrument.

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Sterilization

The crocodile forceps, point coagulator and myoma enucleator must be cleaned with water and a soft brush. In order to prevent damage to the teflon coating of the crocodile forceps, do not use sharp or hard cleaning utensils! The instruments must then be sterilized in a steam sterilizer, gas sterilizer or by immersing them in appropriate disinfection solutions.

Avoid kinks in the green instrument cord. The cord can be sterilized in the same manner as the other instruments. Though steam sterilization is possible, it shortens the serviceable life of the cord.

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**Maintenance of the
Endo-Coagulator**


Protective maintenance is not absolutely necessary. Periodic maintenance however, can help detect a problem before it becomes too serious. This will increase the safety and the serviceable life of the unit. The manufacturer recommends having the endo-coagulator and instruments inspected once a year. Contact your local supplier or the manufacturer regarding maintenance service.

Information on maintenance contracts can be obtained from the manufacturer (see warranty card on the last page).

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
Exchanging fuses

Switch off the instrument. Disconnect the power plug.

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
Unscrew the power fuse holder ⑩ with a coin.
Open the spare fuse holder ⑨ by turning it to the left.

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Take out the new fuse with the appropriate value from the spare fuse holder ⑨ and mount it in the fuse holder.

Caution! Use only the stipulated fuse values.

Note: Do not forget to replace fuses you have taken out.

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
Secure the fuse holders ⑩.
Reconnect the power cord.
Carry out function test.

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Trouble shooting

Symptom:
Instrument does not go on.
(Pilot lamp does not illuminate).

Possible causes:
No power from the power line.
Defective power fuse.
Power cord is not properly connected to power cord
connector.

Remedy:
Check that there is electricity to the wall outlet.
Exchange fuses. Push cord firmly into power cord
connector.

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Symptom:
The dual-tone warning signal sounds.

Possible causes:
The coagulation instrument is not connected or is
defective, the instrument lead is defective.

Remedy:
Replace the instrument lead or the coagulation
instrument.

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Symptom:
The footswitch has been activated, but the
instrument does not heat up.

Possible causes:
The connecting tube of the footswitch either has not
been plugged in correctly, or it is kinked or
defective.
The footswitch is defective.

Remedy:
Check the connecting tube of the footswitch. If
necessary, replace the footswitch and connecting
tube.

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*Caution: Always unplug the instrument
beforehand!*

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Classification

According to MedGV:
this instrument belongs to Group 3
Type of moisture proofness:
protected against dripping water 
Type of protection against electric shocks:
Protection Class I
Degree of protection against electric shocks:
BF-type instrument 

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Test certificate

The type 26027 B/C/D instruments are manufac-
tured in accordance with DIN/IEC 601-1,
VDE 0750 Part 1.

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Technical information

Upon request the manufacturer will provide further information such as complete circuit diagrams, itemized lists of spare parts, adjustment instructions, and other material which can be of use to the qualified service technician when repairing the parts of the instrument the manufacturer has labeled repairable.

By making such information available the manufacturer does not authorize opening or repairs, with the exception of those manipulations described in this instructions manual. This also applies to trained service technicians.

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We reserve the right to make any changes in the construction of the instrument that will contribute to its further development and modification.

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