SERFAS

Energy System

Operation and Maintenance Manual

Rx ONLY

**English, Français, Deutsch, Italiano, Português,
Español, Nederlands, Dansk, Suomi, Norsk,
Svenska, Polski, Ελληνικά, 日本語, 中文, 한국어**

Contents

SERFAS Energy System.....	EN-1
SERFAS Système Energy.....	FR-37
SERFAS Energy System.....	DE-73
SERFAS Sistema Energy.....	IT-109
SERFAS Sistema Energy.....	PT-147
SERFAS Sistema Energy.....	ES-183
SERFAS Energy-systeem.....	NL-221
SERFAS Energy-system.....	DA-257
SERFAS Energy-järjestelmä	FI-293
SERFAS Energy-system.....	NO-329
SERFAS Energy System.....	SV-365
SERFAS Energy System.....	PL-401
Σύστημα SERFAS Energy	EL-437
SERFAS 能量系统	CHS-475
SERFAS Energy システム	JP-507
SERFAS 에너지 시스템	KO-543

SERFAS

Energy System

Operation and Maintenance Manual

Contents

Indications for Use	EN-2
Contraindications	EN-2
Warnings & Cautions.....	EN-2
Symbol Definitions.....	EN-6
Product Description and Intended Use	EN-7
Setting Up the Generator.....	EN-11
Setting Up the Probe.....	EN-12
Setting Up the Footswitch	EN-13
Powering On the System	EN-14
Measuring the Output Voltage of the SERFAS Energy Generator	EN-16
Cleaning.....	EN-21
Troubleshooting	EN-22
Technical Specifications	EN-24
Electromagnetic Compatibility.....	EN-28
Classifications and Approvals	EN-33
Warranty	EN-34
Service and Claims.....	EN-35

Indications for Use

The Stryker SERFAS Energy System is indicated for the resection, ablation, and coagulation of soft tissue and hemostasis of blood vessels in patients undergoing arthroscopic surgery of the knee, shoulder, ankle, elbow, hip and wrist.

Contraindications

The Stryker SERFAS Energy System is contraindicated for any non-arthroscopic procedure or in procedures where a non-conductive irrigant is used. The Stryker SERFAS Energy System is also contraindicated against use with patients for whom an arthroscopic procedure is contraindicated for any reason, and for use with patients having cardiac pacemakers or other electronic implants.

Warnings & Cautions

Please read this manual and follow its instructions carefully. The words **warning**, **caution**, and **note** carry special meanings and should be carefully reviewed:

Warning  **The personal safety of the patient or physician may be involved. Disregarding this information could result in injury to the patient or physician.**

Caution **Special service procedures or precautions must be followed to avoid damaging the instrument.**

Note Special information to make maintenance easier or important information more clear.

An exclamation mark within a triangle is intended to alert the user to the presence of important operating and maintenance instructions in the literature accompanying the product.

A lightning bolt within a triangle is intended to warn of the presence of hazardous voltage. Refer all service to authorized personnel.

Warning  **To avoid potential serious injury to the user and patient and/or damage to this device, the user must follow the following warnings:**

FIRE/EXPLOSION WARNINGS

1. DO NOT use this device in the presence of flammable anaesthetics, other flammable gases or objects, near flammable fluids such as skin prepping agents and tinctures, or oxidizing agents. Observe appropriate fire precautions at all times.
2. DO NOT use this device in oxygen enriched atmospheres, nitrous oxide (N₂O) atmospheres, or in the presence of other oxidizing agents, to prevent risk of explosion. Ensure that oxygen connections are not leaking.
3. Electrosurgical components, such as the Probe, may remain hot following activation. Keep all electrosurgical equipment away from flammable materials to avoid combustion.
4. To prevent the risk of fire, DO NOT replace Generator fuses. If it is suspected that fuses are damaged, return Generator to Stryker for repair.

PRIOR TO SURGERY

1. The operator of the SERFAS Energy System should be a qualified physician, having complete knowledge of the use of this equipment and awareness of the risks associated with surgical procedures.
2. The operator of the SERFAS Energy System should have current knowledge of technological advances in surgical products and techniques.
3. The operator of the SERFAS Energy System should read this manual thoroughly and be familiar with its contents prior to operating the equipment.
4. The operator of the SERFAS Energy System should be sure that the system functions as outlined in this manual prior to a surgical procedure. The SERFAS Energy System was fully tested at the factory before shipment.
5. SERFAS Energy components are designed to be used together as a system. Use only the appropriate SERFAS Energy Footswitch, Probe, and RF Generator.
6. The SERFAS and SERFAS Energy systems are independent systems. Their components are not compatible with each other.
7. Operators of the SERFAS Energy System should be experienced in arthroscopic and electrosurgical practices and techniques.
8. Carefully unpack the unit and ensure that all components are accounted for and remain undamaged from shipment. Inspect the probe cable for any damage to insulation. If damage to any component is detected, refer to the 'Service and Claims' section of this manual.
9. DO NOT use the SERFAS Energy System on patients with cardiac pacemakers or other electronic device implants. Doing so could lead to electromagnetic interference and possible death.

10. Ensure the proper connection of the primary power cord of the SERFAS Energy System to a grounded receptacle. DO NOT use extension cords or adapter plugs, to prevent risk of electric shock.
11. DO NOT wrap the SERFAS Energy Probe cable around metal objects or the induction of hazardous currents may result.
12. Position the cables to avoid contact with the patient, electrodes, cables, and any other electrical leads which provide paths for high frequency current.
13. Position the Generator so the fan directs the flow of air away from the patient.
14. When the SERFAS Energy System and physiological monitoring equipment are used simultaneously on a patient, position any monitoring electrodes as far as possible from the surgical electrodes. Monitoring equipment using high frequency, current-limiting devices is recommended. Needle monitoring electrodes are NOT recommended.
15. Smoke generated during electrosurgical procedures may be harmful to surgical personnel. Take appropriate precautions by wearing surgical masks or other means of protection.

DURING SURGERY

1. DO NOT use the SERFAS Energy System with non-conductive media (e.g. sterile water, dextrose, air, gas, glycine, etc.). Use only conductive irrigants such as saline or Ringer's lactate in order for system to function.
2. DO NOT allow the patient to come into contact with grounded metal objects or objects which have an appreciable capacitance to the earth, such as a surgical table frame or instrument table, etc., to prevent risk of shock. The use of antistatic sheeting is recommended for this purpose.
3. DO NOT activate the SERFAS Energy System for prolonged lengths of time when not in contact with tissue. Doing so may lead to unintentional damage to surrounding tissue.
4. When the SERFAS Energy System is activated, the conducted and radiated electrical fields may interfere with other electrical medical equipment. Provide as much possible distance between the RF Generator and other electronic medical equipment.
5. Select the lowest output power required to prevent patient injury.
6. Maintain the active electrode in the field of view at all times to avoid tissue damage.
7. Remove the probe from the surgical site and place away from metallic objects when not in use. Probes should be separated from other electrosurgical equipment to avoid inadvertent electrical coupling between devices. Inadvertent activation may cause user/patient injury and or product damage.

8. Keep the ends of the Probe cable connectors, Footswitch cable connectors, and Generator receptacles away from all fluids.
9. DO NOT activate the SERFAS Energy System until the probe is properly positioned in the patient.
10. Ensure that the probe tip, including the return electrode, is completely surrounded by irrigant solution during use.
11. Keep the activation indication lights and speaker in field of view and hearing at all times during activation. The light and sound are important safety features.
12. DO NOT contact metal objects, such as an endoscope or metal cannula, while activating the probe. Damage to the probe tip or other devices may result.
13. DO NOT obstruct the fan (located near the rear of the Generator).
14. Failure of the system may result in an unintended increase in output power.
15. When the SERFAS probe is activated, do not touch any part of your body to metallic or conductive objects that are in contact with the saline in the surgical site, as this will result in electrical shock.

AFTER SURGERY

1. DO NOT attempt to reuse or resterilize any product labeled “Single-Use”, as this may lead to equipment malfunction, patient/user injury, and/or cross contamination.
2. SERFAS Energy Probes will not activate if reprocessed or resterilized.
3. DO NOT use flammable agents for cleaning and disinfection of the SERFAS Energy Generator or Footswitch.
4. DO NOT remove the cover of the Generator to avoid electric shock and product damage.
5. Attempt no internal repairs or adjustments, unless specified otherwise in this manual. Units requiring repair should be returned to Stryker for repair.
6. Disconnect the SERFAS Energy Generator from the electrical output when inspecting fuses.

The warranty is void if any of these warnings are disregarded.

Warning Federal law (United States of America) restricts this device to use by or on order of a physician.

Stryker Endoscopy reserves the right to make improvements in the product(s) described herein. Product(s), therefore, may not agree in detail to the published design or specifications. All specifications are subject to change without notice. Please contact the local Stryker Endoscopy distributor listed in the Other Service section or call your local Stryker Endoscopy sales representative or agent for information on changes and new products.

Symbol Definitions

	Alternating Current
	Attention: consult accompanying documents
	Dangerous Voltage
	Defibrillation-proof, Type BF Applied Part
	Equipotential Ground
	Fuses
	Non-Ionizing Radiation
	Power On/Off
	Relative Humidity Range
	Temperature Range
	Protective Earth Ground
	Denotes compliance with UL 60601-1 and CSA C22.2 No. 601.1-M90.
	Denotes conformance to: 93/42/EEC Medical Device Directive
	Manufacturer
	Stryker European Representative
	This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact the manufacturer or other authorized disposal company to decommission your equipment.

Product Description and Intended Use

The Stryker SERFAS Energy System is a bipolar, radio-frequency (RF) electro-surgical system designed to deliver a high frequency output for use in arthroscopic and orthopedic procedures. The SERFAS Energy System consists of three main components:

1. The SERFAS Energy Generator
2. The disposable, single-use SERFAS Energy Probe
3. The SERFAS Energy Footswitch

The SERFAS Energy Generator

The SERFAS Energy Generator is the radio-frequency energy delivering unit, designed to deliver a high-frequency output to the treatment site to perform electro-surgical procedures. The front console panel features controls for adjusting the CUT mode power outputs and provides visual cues for CUT power levels, CUT and COAG activation, connection of Probe and Footswitch, and error code warnings.

The rear console panel provides ports for connecting the SERFAS Energy Generator to other Stryker equipment, including Sidne™ and future Stryker firewire compatible devices.

The features of the front and rear console panels are listed in Figures 1a and 1b.


Figure 1a: The SERFAS Energy Generator, front panel

1. Power Switch: Powers the Generator on and off
2. Error Indicator: A red indicator will illuminate to indicate an error. (See error code on Cut Set Point Display Window for determining type of error)
3. Coagulation Activation Indicator: A blue indicator will illuminate when RF energy is delivered due to the COAG pedal/button being pressed

4. Cut Activation Indicator: A yellow indicator will illuminate when RF energy is delivered due to the CUT pedal/button being pressed
5. Cut Set Point Display Window: Displays Cut power level setting and error codes
6. Cut Output Level Change: Increase or decrease the energy delivered during the Cut mode
7. Probe Indicator: Indicator will illuminate upon proper connection of Probe to Generator
8. Footswitch Indicator: Indicator will illuminate upon proper connection of Footswitch to Generator
9. Footswitch Connector: Connects to the SERFAS Energy Footswitch
10. Probe Connector: Connects to the SERFAS Energy Probes


Figure 1b: The SERFAS Energy Generator, rear panel

11. Side Port: Connects to the Side Console to enable voice operation
12. Hand Control Switch: Users who prefer not to use the hand control feature may deactivate this functionality by switching the hand control switch to the Off position
13. Firewire Connectors (Refer to iSwitch Operating and Maintenance Manual P/N 1000-400-700 for instructions for operating the Wireless Universal Footswitch with the SERFAS Energy Generator and the appropriate manuals for other Firewire devices)
14. Speaker Volume Control
15. Fuse Panel
16. AC Power Inlet
17. Equipotential Ground Plug

The SERFAS Energy Probe

The SERFAS Energy Probe is used to deliver the high-frequency energy to the treatment site within the patient. Probe appearance and functions will vary depending on the probe used.

The features of the probe are listed in Figure 2 below.


Figure 2: The SERFAS Energy Probe

1. Probe Tip: point of delivery of RF energy used for Cut and Coagulation
2. Probe Shaft: Insulated
3. Yellow CUT Button: Used for activation of Cut function
4. Blue COAG Button: Used for activation of Coag function
5. Cut Level Control Button: Used to adjust the Cut level
6. Suction Adapter: Used to connect probe to vacuum source (for SERFAS Energy Suction Probes only)
7. Cable Connector: Connects to front panel of SERFAS Energy Generator

The SERFAS Energy Footswitch

The SERFAS Energy Footswitch is an optional system feature that can be used to activate the SERFAS Energy Probes in both the CUT and COAG modes as well as to increase or decrease the CUT output level settings.

Figure 3 below lists the features of the Footswitch:

Figure 3: The SERFAS Energy Footswitch

1. Yellow Cut Pedal: Used for activation of CUT function
2. Blue Coag Pedal: Used for activation of COAG function
3. Minus (-) Pedal: Used to adjust CUT level down
4. Plus (+) Pedal: Used to adjust CUT level up
5. Footswitch Connector: Connects to front panel of SERFAS Energy Generator

Note Stryker considers instructional training, or in-service, as an integral part of the SERFAS Energy System. Your local Stryker sales representative will perform at least one in-service at your convenience to help set up your equipment and instruct you and your staff on its operation and maintenance. To schedule an in-service, contact your local Stryker representative after your equipment has arrived.

Note Use only the appropriate SERFAS Energy accessories with the SERFAS Energy Generator.

Setting Up the Generator

Initial Set-up of the Generator

Warning Prior to initial use, ensure that all package inserts, warnings, cautions, and Instructions for Use are read and understood.

Warning Before each use, check that all Generator indicator lights and audio signals are functional. Make sure that the power cable plug is properly connected to the Generator receptacle.

1. The SERFAS Energy Generator should be placed on a Stryker cart or on any sturdy table or platform. Refer to hospital procedures or local codes for detailed information.
2. Provide at least four to six inches of space around the sides of the SERFAS Energy Generator for convection cooling. It is normal for the top and rear panel to be warm when system is used continuously for extended periods of time.
3. Power cords use with the SERFAS Energy System should comply with appropriate electrical standards and be suitable for hospital use. Plug the SERFAS Energy Generator into a grounded receptacle. DO NOT use extension cords or three-prong to two-prong adapters. The mains power cord assembly should be checked periodically for damaged insulation or connectors.
4. Plug the power cord into back of Generator.
5. Plug the power cord into AC mains.

Setting Up the Probe

Warning Before each use examine all accessories and connections.

Ensure that the accessories function as intended. Improper connection may result in arcing or malfunction of the Probe or Generator which can result in an unintended surgical effect, injury, or product damage.

Warning Be sure that no liquid is present between connections to the RF Generator or Probe. Connection of wet accessories may lead to electric shock or electrical short.

1. Connect the probe connector on the probe cable to the gray probe receptacle on the front panel of the SERFAS Energy Generator (See Figure 4).

Figure 4: Connecting the Probe to the Generator

2. The Probe indicator, on the front panel, will illuminate following proper connection.

Warning Improper connection may lead to arcing, sparking or malfunction of the unit, which can result in injury to the patient, user, and/or damage to the system.

3. If probe is a suction style probe, connect vacuum hose to the suction adapter to utilize the suction feature.

Note To unplug the Probe from the Generator, grasp the knobbed portion of the connector and pull straight out.

Setting Up the Footswitch

If the optional footswitch is selected for use:

1. Connect Footswitch Cable to the Footswitch connector on the front panel of the Generator (See Figure 5).

Figure 5: Connecting the Footswitch to the Generator

2. The Footswitch indicator, on the front panel, will illuminate following proper connection.

Note To unplug the Footswitch from the Generator, grasp the knobbed portion of the connector and pull straight out.

Note Before operating the SERFAS Energy System, ensure that all components have been set up according to the instructions in the Setting up the SERFAS Energy System section of this manual. Be sure that all components are functioning properly, that all indicators are lit as appropriate, and that tone sounds during system self-test.

Powering On the System

1. Press the power switch to power the Generator on. A tone will sound during the system self-test. If a probe is connected during start-up, a beep will sound following the tone.
2. Be sure that the SERFAS Energy Generator and all components are properly connected. All components are properly connected when the following are present on the front panel of the SERFAS Energy Generator:
 - Probe indicator is lit
 - CUT level window displays a power level between 1 and 11
 - Footswitch indicator is lit (if optional Footswitch is used)

Note If probe is not properly attached, the probe indicator will not be lit.

Modes of Operation

CUT MODE

1. Cut level can be adjusted by using either the Up/Down arrow buttons on the front panel of the Generator, or by using the Plus/Minus pedals on the footswitch, or by using the Cut Level Control adjustment button on the probe.
2. Activation of the Cut mode can be achieved by pressing the yellow pedal on the Footswitch, or by pressing the yellow button on the Probe.
3. When Cut mode is activated, the yellow Cut indicator will illuminate and a higher-pitched tone will sound continuously for the duration of the activation.

COAG MODE

1. Power levels for the Coag mode are not adjustable and power settings for the Coag mode are dependant on the probe style used.
2. Activation of the Coag mode can be achieved by pressing the blue pedal on the footswitch, or by pressing the blue button on the Probe.
3. When Coag mode is activated, the blue Coag indicator will illuminate and a lower-pitched tone will sound continuously for the duration of the activation.

Hand Control Adjustment

The ability to use the hand control feature on the probe can be adjusted by using the Hand Control switch located on the back panel of the Generator. When in the On position, hand control is enabled. When in the Off position, hand control is disabled.

Volume Adjustment

Volume of tones can be adjusted by using the Speaker Volume Control knob located on the back panel of the Generator.

After Surgery

1. Withdraw Probe from the patient. Be sure that the Probe is NOT activated during withdrawal.
2. Unplug the Probe from the Generator
3. Properly dispose of single-use Probes.

Operating the SERFAS Energy System with the iSwitch:

Warning When the SERFAS Energy System is interconnected with other medical electrical equipment, leakage currents may be additive. To minimize total patient leakage current, any Type BF applied part should be used together with other Type BF applied parts. Ensure all systems are installed according to the requirements of IEC 60601-1-1.

Refer to iSwitch Operating and Maintenance Manual P/N 1000-400-700 for instructions for operating the iSwitch with the SERFAS Energy Generator.

Measuring the Output Voltage of the SERFAS Energy Generator

Below is the procedure to measure the power output of the presets on the SERFAS Energy Generator.

Equipment

- Electrosurgical analyzer (*such as the Fluke Model 454A tester or Fluke QA-ES Series II tester*)
- 2 test leads with banana clips
- SERFAS Energy Generator
- SERFAS Energy Probe
- SERFAS Energy Footswitch

Setup

1. Connect the electrosurgical analyzer to the power source.
2. Plug the two test leads in to the Active and Dispersive/Return input ports of the electrosurgical analyzer as shown in in the following picture. (*Refer to the electrosurgical analyzer user guide for location of Active and Dispersive/Return input ports.*)


Figure 1

3. Connect the SERFAS Energy Generator to the power source.

4. Connect the probe cable to the gray probe receptacle on the front panel of the SERFAS Energy Generator (see Figure 2).

Figure 2

5. Connect the footswitch cable to the footswitch connector on the front panel of the SERFAS Energy Generator (see Figure 3).

Figure 3

6. Identify the active and the dispersive output on the probe.

7. Connect the electrosurgical analyzer active input test lead to the active output of the SERFAS Energy probe (see Figure 4).


Figure 4

8. Connect the electrosurgical analyzer Dispersive/Return test lead to the SERFAS Energy probe (see Figure 5).

Figure 5

Test Procedure

1. Turn on the SERFAS Energy Generator.
2. Turn on the electrosurgical analyzer.
Set the electrosurgical analyzer to 250 ohm load. (*Refer to the electrosurgical analyzer user guide for instructions on setting the output power load.*)
3. On the SERFAS Energy Generator, select the Cut level by using either the Up/Down arrow buttons on the front panel, or the +/- pedals on the footswitch.
4. Activate the Cut mode at the selected level by pressing the yellow pedal on the footswitch.
5. Allow a few seconds for the measurement to stabilize.
6. Take note of the displayed power in watts and the generator cut level. The measurements should fall within the range specified in Table 1.

Table 1: SERFAS Generator Output Power at 250 Ohm Load

Cut Mode	Load	Output Power Range (Watts)	
		Minimum	Maximum
Level 1	250 ohms	5	8
Level 2	250 ohms	24	35
Level 3	250 ohms	42	63
Level 4	250 ohms	67	100
Level 5	250 ohms	86	129
Level 6	250 ohms	110	164
Level 7	250 ohms	132	199
Level 8	250 ohms	157	235
Level 9	250 ohms	180	270
Level 10	250 ohms	204	306
Level 11	250 ohms	233	350

7. Resume from Step 6 to test the other power output presets.

Note Note: Nominal impedance is 200 ohms. However, due to output modulation, the electrosurgical analyzer power signal is unstable at 200 ohms. Therefore, this test is conducted at 250 ohms, where modulation does not occur.

If this test reveals that power output levels that are outside of the suggested values, please return the SERFAS Energy generator for service, as described in Services and Claims.

Cleaning

SERFAS Energy Generator

The SERFAS Energy Generator cannot be sterilized. If the Generator needs to be cleaned, wipe it down with a damp cloth or sponge. Use only non-abrasive cleaning agents and DO NOT allow liquid to enter connectors on the Generator.

Footswitch

Caution Do not sterilize the SERFAS Energy Footswitch. DO NOT immerse the SERFAS Energy Footswitch in liquid. Use only non-abrasive cleaning agents and DO NOT allow liquid to enter connections. Doing so will result in damage.

Wipe the outside of the SERFAS Energy Footswitch with a mild detergent and water according to standard hospital practices.

Probes

Caution DO NOT attempt to clean, resterilize, or reprocess the Probes as this may lead to product failure or spread of infectious disease. SERFAS Energy Probes will not activate if reprocessed or resterilized.

Probes are shipped sterile and are intended to be used ONE TIME ONLY.

Equipment Disposal Guidelines

The SERFAS Energy Generator contains electronic printed circuit assemblies. Equipment should be disposed of at the end of its useful life in accordance with applicable national or institutional related policies relating to obsolete electronic equipment.

Dispose of any system accessories according to normal institutional practice relating to potentially contaminated items.

Troubleshooting

Note If a fault condition should occur, the display on the front panel will display an error code and an error beep will sound. The interpretation of the error code, as well as possible causes and solutions are listed in Table 1.

Note To recover system from Errors P1-P4, E1-E9, press any button or re-connect the Probe to the Generator. If errors persist, refer to Solution section in Table 1. If solutions do not prevent error occurrence, refer to the Service and Claims section of this manual to obtain repair service.

Recoverable Errors & Corresponding Codes		
Error Code	Reason for Error	Solutions
P1	Probe: Probe not ready	<ul style="list-style-type: none"> • Check Probe connection to Generator. Reconnect Probe.
P2	Probe: Probe expired	<ul style="list-style-type: none"> • Discard Probe and replace with new Probe
P3	Probe: Probe identification invalid	<ul style="list-style-type: none"> • Replace Probe
P4	Probe: Probe communication error	<ul style="list-style-type: none"> • Check Probe connection to Generator. Reconnect Probe. • Replace Probe
P5	Probe: Maximum time of use exceeded	<ul style="list-style-type: none"> • Discard probe and replace with new probe
E1	Unable to identify error	<ul style="list-style-type: none"> • Check Probe for damage • Replace Probe
E2	RF Power too high	<ul style="list-style-type: none"> • Check Probe for damage • Replace Probe
E3	RF Voltage too high	<ul style="list-style-type: none"> • Check Probe for damage • Replace Probe
E4	RF Current too high	<ul style="list-style-type: none"> • Check Probe for damage • Replace Probe
E5	RF delivery exceed countinous use limit	<ul style="list-style-type: none"> • Clear error and continue
E6	Generator temperature too high	<ul style="list-style-type: none"> • Ensure proper airflow on sides and rear of Generator
E7	Low impedance detected	<ul style="list-style-type: none"> • Check Probe for damage • Replace Probe
E8	Switch pressed upon Probe/Footswitch correction	<ul style="list-style-type: none"> • Check that switches on Probe, Footswitch, and Generator are not pressed

Non-Recoverable Errors & Corresponding Codes		
Error Code	Reason for Error	Solutions
F1	AC voltage incorrect	<ul style="list-style-type: none"> • Cycle power to the unit (turn the unit off and then on)
F4	COP Watchdog	<ul style="list-style-type: none"> • Same as above
F5	Oscillator fault	<ul style="list-style-type: none"> • Same as above
F6	Software fault	<ul style="list-style-type: none"> • Same as above
F7	Hardware fault	<ul style="list-style-type: none"> • Same as above
F8	Power-on self test error	<ul style="list-style-type: none"> • Same as above

Note If errors (F1-F8) persist, refer to the Service and Claims section of this manual to obtain repair service.

Solutions	
Problem	Possible Solutions
System does not power up after power switch is pressed	<ul style="list-style-type: none"> • Check that power cord is properly connected to Generator and into an appropriately grounded outlet.
Green footswitch indicator light does not illuminate	<ul style="list-style-type: none"> • Check that Footswitch is properly connected to the Generator and no damage exists to cable or connector. • Change Footswitch
Green Probe indicator light does not illuminate	<ul style="list-style-type: none"> • Check that Probe is properly connected to the Generator and no damage exists to cable or connector. • Change Footswitch

Note If disturbance occurs on the video monitor or any other electronic instrument, the user should:

Caution Ensure the Probe cable is not near any other instrument cables.

Caution Ensure the Endoscope has a non-metallic coupling ring (i.e. Stryker CE marked Endoscopes).

Technical Specifications

SERFAS Energy Probe

Overall Handpiece Length	17.5 cm
Overall Cable Length	3m
Working Length	10 - 17.5cm
Shaft Diameter	2.5mm - 5.0mm
Shaft Bend Angle	0-45°
Active Tip Orientation	0-90°
Supplied Sterile and For Single-use ONLY	
Sterilization Method	EtO

SERFAS Energy Generator

Dimensions	16.9" x 12.5" x 3.5" (42.9cm x 31.8 cm x 8.9cm)
Weight	5.6 kg (12.2 lb.)
Operating Temperature Range	10°C to 40°C
Operating Relative Humidity	15% to 80% (non-condensing)
Transport & Storage Temperature Range	-34°C (-29°F) to 65°C (149°F)
Transport & Storage Relative Humidity	10% to 85% (non-condensing)
Electrical Ratings	100-120V ~ 50/60Hz, 7A 220-240V ~ 50/60Hz, 4A
Inlet Fuse	8.0A, 250V
Speaker Volume	Adjustable between 45dB (minimum) and 65dB (max) at 1 meter

RF OUTPUT

Waveform	200 kHz + 1%, square wave, Crest factor <1.3 @200ohm
Maximum Power (CUT Mode)	400 Watts @ 200ohms
Maximum Power (COAG Mode)	90 Watts @ 200ohms

Generator Output Graphs

Output power at each set point with specified load resistance (per IEC 60601-2-2, sub clause 6.8.3) is given in the graphs below.


Figure 6: Output Power vs. Set Point at 200ohms Resistive Load


Figure 7: Output Power (CUT) vs. Load Resistance


Figure 8: Output Power (COAG) vs. Load Resistance


Figure 9: Maximum Open Circuit Voltage vs. Set Point


Figure 10: Output COAG Power vs. Set Point at 200ohms Resistive Load


Figure 11: Maximum Open Circuit Voltage vs. Set Point

Electromagnetic Compatibility

Like other electrical medical equipment, the SERFAS Energy System requires special precautions to ensure electromagnetic compatibility with other electrical medical devices. To ensure electromagnetic compatibility (EMC), the SERFAS Energy System must be installed and operated according to the EMC information provided in this manual.

Note The SERFAS Energy System has been designed and tested to comply with IEC 60601-1-2:2001 requirements for EMC with other devices.

Caution  The SERFAS Energy System may be interfered with by other equipment, including portable and mobile RF communication equipment, even if such equipment meets the applicable emissions requirements.

Warning  Do not use cables or accessories other than those provided with the SERFAS Energy System, as this may result in increased electromagnetic emissions or decreased immunity to such emissions.

Warning  If the SERFAS Energy System is used adjacent to or stacked with other equipment, observe and verify normal operation of the SERFAS Energy System in the configuration in which it will be used prior to using it in a surgical procedure. Consult the tables below for guidance in placing the SERFAS Energy System.

Guidance and Manufacturer's Declaration: Electromagnetic Emissions

SERFAS Energy Systemis intended for use in the electromagnetic environment specified below. The customer or the user of SERFAS Energy Systemshould ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic Environment - guidance
RF emissions CISPR11	Group 1	The SERFAS Energy Systemmust emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR11	Class B	SERFAS Energy Systemis suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	Class A	
Voltage Fluctuations/ flicker emissions IEC61000-3-3	Complies	

Guidance and Manufacturer's Declaration: Electromagnetic Immunity

SERFAS Energy Systemis intended for use in the electromagnetic environment specified below. The customer or the user of SERFAS Energy Systems should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electro-magnetic Environment: Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	±6kV contact ±8kV air	±2,4,6kV contact ±2,4,8kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	±1kV differential mode ±2kV common mode	±0,5, 1kV differential mode ±1, 2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec.	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of SERFAS Energy System requires continued operation during power mains interruptions, it is recommended that SERFAS Energy System be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	N/A	Power-frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: Ut is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration: Electromagnetic Immunity

SERFAS Energy Systems intended for use in the electromagnetic environment specified below. The customer or the user of SERFAS Energy Systems should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment: Guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80MHz to 2.5 GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the SERFAS Energy Systems system, including its cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p style="text-align: center;">Recommended Separation Distance $d = 1.17 \sqrt{P}$</p> <p style="text-align: center;">$d = 1.17 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.33 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^(a), should be less than the compliance level in each frequency range ^(b).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

(a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SERFAS Energy Systems system is used exceeds the applicable RF compliance level above, the SERFAS Energy Systems system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the SERFAS Energy System unit.

(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the SERFAS Energy System

The SERFAS Energy System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the SERFAS Energy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SERFAS Energy System as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power (W) of transmitter	Separation distance (m) according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Classifications and Approvals

Complies with medical safety standards:

- IEC 60601-1: 1998 + A1:1991 + A2:1995
- AS 3200.1.0: 1998
- IEC 60601-1-2: 2001
- IEC 60601-2-2: 1998
- UL 60601-1: 2003
- CSA C22.2 No. 601-1-M90

Classification

Type of protection against electric shock

- Class I Equipment

Degree of protection against electric shock

- Type BF Defibrillation-Proof Applied Part

Degree of protection against harmful ingress of water

- Generator: IEC 60601-2-2: Requirement per clause 44.3
- Probe: IEC 60601-2-2: Requirement per clause 44.6
- Footswitch: IEC60601-2-2: Requirement per clause 44.6,
IPX7 Water-tight Equipment

Mode of Operation

- Continuous with intermittent probe activation

Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.

Warranty

This warranty applies to customers in the United States only. Outside of the U.S.A., contact your Stryker sales representative or your local Stryker subsidiary. Stryker Endoscopy warrants the SERFAS Energy System against defects in both materials and workmanship to the registered owner at the time of purchase. All components of the system are covered by the warranty as described below.

This warranty does not apply to any unit which has been subject to misuse, neglect, improper installation or that which has been altered, adjusted, or tampered with by any person other than Stryker authorized personnel.

If, upon examination by authorized service personnel, it is determined that the malfunction is due to misuse or abuse, warranty provisions will not apply. An estimate of the cost of repair work will be given to the customer prior to servicing and repairing the unit.

The customer is responsible for returning the defective equipment to the factory at his or her own expense. Stryker or its representative will service the unit, repair or replace any defective parts thereof, and return the unit.

If, upon examination, it is determined that the fault has been caused by misuse or abnormal conditions of operation, the repairs will be billed to the customer as out-of-warranty repairs.

Instruments repaired under Stryker's standard repair program will be issued a thirty day warranty against defects in both materials and workmanship, provided the original warranty period had passed. Instruments submitted due to defects in materials and workmanship during the warranty period will be repaired at no charge to the customer.

The warranty set forth herein is exclusive and in lieu of all other warranties, remedies, obligations, and liabilities of Stryker, expressed or implied, including the implied warranties of merchantability and fitness for use and of consequential damages. These products are being sold only for the purpose described herein, and such warranty only runs to the purchaser. In no event shall Stryker be liable for any breach of warranty in any amount exceeding the purchase price of the product.

No agent, employee or representative of Stryker Endoscopy has the authority to bind the Company to any other warranty, affirmation, or representation concerning this instrument.

This warranty is valid only to the original purchaser of Stryker products directly from a Stryker authorized agent. The warranty cannot be transferred or assigned by the original purchaser.

The warranty periods for components of SERFAS Energy are as follows:

- SERFAS Energy Generator: 1 year from shipment date
- Footswitch: 90 days from shipment date
- Probes: Single Use Only

Caution Never open the SERFAS Energy Generator or attempt any service not described in this manual.

Service and Claims

If service is needed either during or after the warranty period:

1. Contact Stryker at 1-800-624-4422 or phone your local Stryker sales representative
2. Package all the components carefully in the original shipping container if possible
3. Ship SERFAS Energy, pre-paid and insured to:
Stryker Endoscopy Customer Service
Attn: Repair Department

Note Stryker reserves the right to make improvements in the product(s) described herein. Product(s) therefore may not agree in detail to the published design or specifications. All specifications are subject to change without notice. Please contact your local Stryker Distributor or call your local Stryker sales representative or agent for information on changes and new products.

For service in the USA, call your Stryker Endoscopy representative or Stryker Endoscopy Customer Service at 1-800-624-4422. Outside of the USA, please contact your Stryker Endoscopy distributor at one of the following locations.

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