

Electric Handpieces and Attachments

Sterilization & Cleaning

Sterilization

To obtain optimal performance and help prevent damage, perform one of the following sterilization procedures:

250°F Gravity:

- Gravity displacement sterilizer
- 250-254°F (121-123°C)
- Double wrapped in an instrument tray or fully perforated sterilization case
- 50-minute minimum exposure time
- 8-minute minimum drying time

270°F Gravity:

- Gravity displacement sterilizer
- 270-275°F (132-135°C)
- Double wrapped in an instrument tray or fully perforated sterilization case
- 35-minute minimum exposure time
- 8-minute minimum drying time

Dynamic Air Removal (PreVac or HiVac):

- Dynamic Air Removal sterilizer
- 270-275°F (132-135°C)
- Double wrapped in an instrument tray or fully perforated sterilization case
- 4-minute minimum exposure time
- 8-minute minimum drying time

NOTE: Stryker does not recommend "Flash" sterilization for routine sterilization of surgical instruments. "Flash" sterilization should only be used when individual surgical instruments require immediate sterilization and use.

"Flash" Autoclave:

- Gravity displacement sterilizer
- 270-275°F (132-135°C)
- Unwrapped in an instrument tray
- 10-minute minimum exposure time

"Flash" Dynamic Air Removal (PreVac or HiVac):

- Dynamic Air Removal sterilizer
- 270-275°F (132-135°C)
- Unwrapped in an instrument tray
- 4-minute minimum exposure time
- 1-minute or no dry time

*Sterilization Recommendations Validation based on AAMI protocol.

To Clean the Handpiece:

⚠ DO NOT immerse or place any component in a washer/sterilizer.

⚠ DO NOT allow water to run into the distal end of the handpiece or into electrical connections. Failure to comply may result in handpiece damage.



1. Wash exterior of all CORE components with Medical Grade Cleaner or soap solution.

⚠ CAUTION! DO NOT run water directly into open ends of handpiece.

2. Hold handpiece on an incline, distal tip pointing down, flush with running water to rinse away cleaning solutions.

3. Use only Stryker recommended cleaning brush to remove debris around bur and blade mounts.

4. Wipe dry with lint free cloth or force air dry.




To Clean a Handswitch:

1. Wipe clean with a mild detergent and water or alcohol.

2. Extend the finger rest to clean all surfaces and rinse with water.

3. Wipe dry with a lint free towel or force air dry.



To Clean an Attachment:



1. Clean an attachment with a brush or sponge, a hospital enzymatic cleaning solution and water.

2. Use Stryker recommended cleaning brush to clean cannula area. **DO NOT** use a brush in the nose or cannula of an Elite or Saber SD series attachment.

3. The distal tips of attachments may soak in a cup or beaker of detergent for 5-10 minutes.

4. Hold the attachment on an incline and flush with running water to rinse away the cleaning solution.

5. Use a lint free towel to dry or force air dry.



Lubrication Recommendations

⚠ Most attachments do not require lubrication.

To lubricate an Elite attachment, use Stryker TPS lubricant REF 5110-901-900.

To lubricate a Saber SD Series attachment use Stryker ContraClean REF 5110-901-920.

⚠ DO NOT lubricate after every use. Lubricate every 2 months or 40 sterilization cycles.