

INSTRUCTIONS



OLYMPUS ENDO-THERAPY

BIOPSY FORCEPS

FB-7U-1	FB-34C/K-1
FB-11K-1	FB-35C-1
FB-13E/K/Q/U-1	FB-36C/K-1
FB-15C/K-1	FB-37K/U-1
FB-19SX/C/K/N-1	FB-38W-1
FB-20C-1	FB-44D-1
FB-21SX/C/K-1	FB-50K/Q/U-1
FB-22C-1	FB-51K-1
FB-23K-1	FB-52C/K-1
FB-24E/K/Q/U-1	FB-53K/Q/U-1
FB-25F/K-1	FB-54K/Q/U-1
FB-26N-1	FB-55K/Q/U-1
FB-28R/U/Y-1	

AUTOCLAVABLE

USA: CAUTION: Federal law restricts this device to sale by or on the order of a physician.

CE

Contents

Symbols

Important Information-Please Read Before Use

- Intended Use
- Instruction Manual
- User Qualifications
- Instrument Compatibility
- Reprocessing and Storage
- Repair and Modification
- Signal Words

Chapter 1 Checking the Package Contents

- 1.1 Checking the Package Contents

Chapter 2 Instrument Nomenclature and Specifications

- 2.1 Nomenclature and Functions
- 2.2 Specifications

Chapter 3 Preparation, Inspection and Operation

- 3.1 Preparation
- 3.2 Inspection
- 3.3 Operation

Chapter 4 Reprocessing

- 4.1 General Policy
- 4.2 Required Reprocessing Equipment
- 4.3 Cleaning
- 4.4 Lubrication
- 4.5 Sterilization

Chapter 5 Storage

- 5.1 Inspection Before Storage
- 5.2 Storage

Symbols

The meaning(s) of the symbol(s) shown on the pack and/or this instrument are as follows:



Refer to instructions.

Faint, illegible text, likely bleed-through from the reverse side of the page.

Contents

Chapter 1: Choosing the Package

1.1. Choosing the Package

Chapter 2: Instrumentation

2.1. Instrumentation

2.2. Specifications

Chapter 3: Preparation, Inspection and Question

3.1. Preparation

3.2. Inspection

3.3. Question

Chapter 4: Processing

4.1. Processing

4.2. Instrumentation

4.3. Question

4.4. Instrumentation

4.5. Instrumentation

Chapter 5: Storage

5.1. Instrumentation

5.2. Instrumentation

Important Information-Please Read Before Use

Intended Use

Do not use these instruments for any purpose other than their intended uses.

- **FB-7U-1, FB-11K-1, FB-13E/K/Q/U-1, FB-23K-1, FB-25F/K-1, FB-26N-1, FB-28R/U/Y-1, FB-38W-1, FB-50K/Q/U-1, FB-51K-1, FB-53K/Q/U-1, FB-54K/Q/U-1, FB-55K/Q/U-1**

These instruments have been designed to be used with an Olympus endoscope to collect tissue within the digestive tract.

- **FB-15C/K-1, FB-19SX/C/N-1, FB-21SX/C/K-1, FB-34C/K-1, FB-44D-1, FB-52C/K-1**

These instruments have been designed to be used with an Olympus endoscope to collect tissue within the digestive tract, respiratory organs, female reproductive organs and urinary organs.

- **FB-20C-1, FB-22C-1, FB-24E/K/Q/U-1, FB-35C-1, FB-36C/K-1, FB-37K/U-1,**

These instruments have been designed to be used with an Olympus endoscope to collect tissue within the digestive tract and respiratory organs.

Instruction Manual

This instruction manual contains essential information using this instrument safely and effectively. Before thoroughly review this manual and the manuals of equipment which will be used during the procedure use the instruments as instructed.

Keep this and all related instruction manuals in a accessible location.

If you have any questions or comments about information in this manual, please contact Olympus.

User Qualifications

The operator of this instrument must be a physician medical personnel under the supervision of a physician and must have received sufficient training in clinical endoscopic technique. This manual, therefore, does explain or discuss clinical endoscopic procedures.

Instrument Compatibility

Refer to the Tables in Section 2.2, "Specifications confirm that this instrument is compatible with the ancillary equipment being used. Using incompatible equipment result in patient injury or equipment damage.

Reprocessing and Storage

This instrument was not sterilized before shipment. Before using this instrument for the first time, reprocess it according to the instructions in Chapter 4, "Reprocessing".

After using this instrument, reprocess and store it according to the instructions in Chapter 4, "Reprocessing" and Chapter 5, "Storage". Improper and/or incomplete reprocessing or storage can present an infection control risk, cause equipment damage or reduce performance.

Repair and Modification

This instrument does not contain user-serviceable parts. Do not disassemble, modify or attempt to repair it; patient or user injury and/or equipment damage can result.

Signal Words

The following signal words are used throughout this manual:

WARNING Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices or potential equipment damage.


NOTE Indicates additional helpful information.

Chapter 1 Checking the Package Contents

1.1 Checking the Package Contents



Match all items in the package with the components shown below. Inspect each item for damage. If the instrument is damaged, a component is missing or you have questions, do not use the instrument; immediately contact Olympus.

This instrument was not sterilized before shipment. Before using this instrument for the first time, reprocess according to the instructions in Chapter 4, "Reprocessing".



Chapter 2 Instrument Nomenclature and Specifications

2.1 Nomenclature and Functions



2.2 Specifications

The compatible Olympus endoscopes are listed in the Tables on the following pages. New endoscopes released after the introduction of this instrument may also be compatible for use in combination with this instrument. For further details, contact Olympus.


WARNING


Use this instrument only in combination with products recommended by Olympus. If combined with products not recommended by Olympus, patient or operator injury, malfunction or equipment damage may result.


Operating Environment


Ambient Temperature	10 to 40°C (50 to 104°F)
Relative Humidity	30 to 85%
Air Pressure	700 to 1060 hPa (0.71 to 1.08 kgf/cm ²) (10.1 to 15.4 psia)


Specifications


Model	FB-7U-1	FB-11K-
	Alligator Type	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	2300	1550
With Needle	NO	
Fenestrated Cups	NO	
	Working length less than 1850 mm;	Working length less than 1100 mm;
Compatible Olympus Endoscopes (All of these parameters should be met.)	EF, GIF, OGF, GF, GTF, JF, TJF, SIF (SIF-10 only), CF, PCF, OSF	EF, GIF, OC GF, GTF, C (Exclude I-, L-length), C
	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3 (Yellow)
Channel Inner Diameter (mm) (Color Code)	ø 3.7, ø 4.2, ø 6 (Orange)	ø 3.7, ø 4.2, ø 6 (Orange)
	ø 5.5 (Pink)	


Model	FB-13E-1	FB-13K-1	
	Standard Type		
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 3.4		
Working Length (mm)	1200	1550	
With Needle	YES		
Fenestrated Cups	NO		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 750 mm;	Working length less than 1100 mm; GIF
	Channel Inner Diameter (mm) (Color Code)	ø 3.7 (Orange)	ø 3.7, ø 6 (Orange)


Model	FB-13Q-1	FB-13U-	
	Standard Type		
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 3.4		
Working Length (mm)	1950	2300	
With Needle	YES		
Fenestrated Cups	NO		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1500 mm; GIF, TJF, CF (I-length only)	Working length less than 1850 r GIF, TJF, C
	Channel Inner Diameter (mm) (Color Code)	ø 3.7, ø 4.2, ø 6 (Orange)	ø 5.5 (Pink)


Model	FB-15C-1	FB-15K-1	
Alligator Jaws			
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 1.8		
Working Length (mm)	1050	1550	
With Needle	NO		
Fenestrated Cups	NO		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 600 mm; CHF, BF, ENF, CYF, HYF, OSF	Working length less than 1100 mm; EF, GIF, OGF, GF, GTF, CHF, CF (Exclude I-, L-length), OSF, BF, ENF, CYF, HYF
	Channel Inner Diameter (mm) (Color Code)	ø 2, ø 2.2 (Blue) ø 2.6 (Green) ø 2.8, ø 3.2 (Yellow)	ø 2, ø 2.2 (Blue) ø 2.6 (Green) ø 2.8, ø 3.2 (Yellow) ø 3.7, ø 6 (Orange)


Model	FB-19SX-1	FB-19K-1	
Standard Type			
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 1.8		
Working Length (mm)	700	1050	
With Needle	NO		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 400 mm; CHF, OSF (OSF-2-35 only), ENF, CYF, HYF	Working length less than 600 mm; CHF, OSF, ENF, CYF, HYF
	Channel Inner Diameter (mm) (Color Code)	ø 2, ø 2.2 (Blue) ø 2.6 (Green)	ø 2, ø 2.1 (Blue) ø 2.1 (Green) ø 2.8, ø 3.2 (Yellow)


Model	FB-19K-1	FB-19N-1
	Standard Type	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø1.8	
Working Length (mm)	1550	1800
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Working length less than 1100 mm; EF, GIF, OGF, GF, CHF, CF (Exclude I-, L-length), OSF, BF, ENF, CYF, HYF	Working length less than 1350 mm; EF, GIF, OGF, GF, GTF, JF, TJF, CHF, CF (Exclude L-length), OSF, BF, ENF, CYF, HYF
	ø 2, ø 2.2 (Blue) ø 2.6 (Green)	ø 2, ø 2.2 (Blue) ø 2.6 (Green)
Channel Inner Diameter (mm) (Color Code)	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow) ø 3.7, ø 4.2, ø 6 (Orange) ø 5.5 (Pink)


Model	FB-20C-1	FB-21S
	Standard Type	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.2	ø 1.8
Working Length (mm)	1050	700
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Working length less than 600 mm; CHF, OSF, BF	Working length less than 400 mm; CHF, OSI (OSF-2-3 only), EN, CYF, HYF
	ø 2.6 (Green) ø 2.8, ø 3.2 (Yellow)	ø 2, ø 2.2 (Blue) ø 2.6 (Green)
Channel Inner Diameter (mm) (Color Code)	ø 2.8, ø 3.2 (Yellow)	ø 2, ø 2.2 (Blue) ø 2.6 (Green)


Model	FB-21C-1	FB-21K-1
	Standard Type	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 1.8	
Working Length (mm)	1050	1550
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 600 mm; CHF, OSF, BF, ENF, CYF, HYF
	Channel Inner Diameter (mm)	Working length less than 1100 mm; EF, GIF, OGF, GF, GTF, CHF, CF (Exclude I-, L-length), OSF, BF, ENF, CYF, HYF
	Channel Inner Diameter (mm)	ø 2, ø 2.2 (Blue)
	Channel Inner Diameter (mm)	ø 2.6 (Green)
	Channel Inner Diameter (mm)	ø 2.8, ø 3.2 (Yellow)
	Channel Inner Diameter (mm)	ø 2.8, ø 3.2 (Yellow)
	Channel Inner Diameter (mm)	ø 3.7, ø 6 (Orange)


Model	FB-22C-1	FB-23K
	Standard Type	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.2	ø 2.4
Working Length (mm)	1050	1550
With Needle	YES	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 600 mm; CHF, OSF, BF
	Channel Inner Diameter (mm)	Working length less than 1100 mm; EF, GIF, C GF, GTF, CF (Exclu i-, L-length) OSF
	Channel Inner Diameter (mm)	ø 2.6 (Green)
	Channel Inner Diameter (mm)	ø 2.8, ø 3.2 (Yellow)
	Channel Inner Diameter (mm)	ø 2.8, ø 3.7, ø 6 (Orange)


Model	FB-24E-1	FB-24K-1	
	Standard Type		
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 2.4		
Working Length (mm)	1200	1550	
With Needle	YES		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 750 mm; EF, CF (S-length only), OSF, BF	Working length less than 1100 mm; EF, GIF, OGF, GF, GTF, CF (Exclude I-, L-length), OSF, BF
	Channel Inner Diameter (mm) (Color Code)	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow) ø 3.7, ø 6 (Orange)


Model	FB-24Q-1	FB-24U	
	Standard Type		
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 2.4		
Working Length (mm)	1950	2300	
With Needle	YES		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1500 mm; EF, GIF, OGF, GF, GTF, JF, TJF, CF (Exclude L-length), PCF (I-length only), OSF, BF	Working length less than 1850 mm; EF, GIF, CF, GF, GTF, TJF, SIF (SIF-10 or CF, PCF, BF
	Channel Inner Diameter (mm) (Color Code)	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow) ø 3.7, ø 4.2, ø 6 (Orange) ø 5.5 (Pink)


Model	FB-25F-1	FB-25K-1	
Standard Type			
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 2.4		
Working Length (mm)	1250	1550	
With Needle	NO		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 800 mm; EF, CF (S-length only), OSF	Working length less than 1100 mm; EF, GIF, OGF, GF, GTF, CF (Exclude I-, L-length), OSF
	Channel Inner Diameter (mm)	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow)
	(Color Code)		ø 3.7, ø 6 (Orange)


Model	FB-26N-1	FB-28R-	
Standard Type			
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 2.4		
Working Length (mm)	1800	2050	
With Needle	NO		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1350 mm; EF, GIF, OGF, GF, GTF, JF, TJF, CF (Exclude L-length), PCF (I-length only), OSF	Working length less than 1600 mm; EF, GIF, OGF, GTF, TJF, CF (Exclude L-length), (I-length only), OSF
	Channel Inner Diameter (mm)	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow)
	(Color Code)	ø 3.7, ø 4.2, ø 6 (Orange)	ø 5.5 (Pink)


Model	FB-28U-1	FB-28Y-1
	Standard Type	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	2300	2700
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1850 mm; EF, GIF, OGF, GF, GTF, JF, TJF, CF, PCF, SIF (SIF-10 only), OSF
	Channel Inner Diameter (mm) (Color Code)	Working length less than 2250 mm; EF, GIF, OGF, GF, GTF, JF, TJF, CF, PCF, SIF, OSF ø 2.8, ø 3.2 (Yellow) ø 3.7, ø 4.2, ø 6 (Orange) ø 5.5 (Pink)


Model	FB-34C-1	FB-34K-
	Standard Type	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 1.8	
Working Length (mm)	1050	1550
With Needle	YES	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 600 mm; CHF, OSF, BF, ENF, CYF, HYF
	Channel Inner Diameter (mm) (Color Code)	Working length less than 1100 mm; EF, GIF, OGF, GTF, JF, CF (Excluding I-, L-length), OSF, BF, E, CYF, HYF ø 2, ø 2.2 (Blue) ø 2.6 (Green) ø 2.8, ø 3.2 (Yellow) ø 3.7, ø 4.2 (Orange)
Other	Exclude BF-1T20/30, BF-XT20/30	


Model		FB-35C-1
Standard Type		
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	1050	
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 600 mm; OSF, BF
	Channel Inner Diameter (mm) (Color Code)	ø 2.8, ø 3.2 (Yellow)


Model		FB-36C-1	FB-36K
Alligator Jaws			
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 2.4		
Working Length (mm)	1050	1550	
With Needle	NO		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 600 mm; OSF, BF	Working length less than 1100 mm; EF, GIF, OLF, GF, GTF, CF (Excludes L-length OSF, BF)
	Channel Inner Diameter (mm) (Color Code)	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow) ø 3.7, ø 4.0 (Orange)


Model	FB-37K-1	FB-37U-1
	Rat Tooth	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	1550	2300
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1100 mm; EF, GIF, OGF, GF, GTF, CF (Exclude I-, L-length), OSF, BF
	Channel Inner Diameter (mm) (Color Code)	Working length less than 1850 mm; EF, GIF, OGF, GF, GTF, JF, TJF, SIF (SIF-10 only), CF, PCF, OSF, BF ø 2.8, ø 3.2 (Yellow) ø 3.7, ø 6 (Orange)


Model	FB-38W-1	
	Rat Tooth	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 1.5	
Working Length (mm)	2500	
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 2050 mm; CHF
	Channel Inner Diameter (mm) (Color Code)	ø 1.7 (Violet)


Model		FB-44D-1
Standard Type		
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 1.15	
Working Length (mm)	1150	
With Needle	NO	
Fenestrated Cups	NO	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 700 mm; CHF, BF, URF, HYF
	Channel Inner Diameter (mm) (Color Code)	ø 1.2, ø 1.5 (Not color coded)


Model		FB-50K-1	FB-50Q-1
Standard Type			
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 3.1		
Working Length (mm)	1550	1950	
With Needle	YES		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1100 mm; GIF	Working length less than 1500 mm; GIF, CF (l-length only)
	Channel Inner Diameter (mm) (Color Code)	ø 3.7, ø 6 (Orange)	ø 3.7, ø 4.2, ø 6 (Orange), ø 5.5 (Pink)


Model	FB-50U-1	FB-51K-1	
	Standard Type		
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 3.1	ø 5.1	
Working Length (mm)	2300	1550	
With Needle	YES		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1850 mm; GIF, TJF, CF	Working length less than 1100 mm; GIF-XT30 only
	Channel Inner Diameter (mm)	ø 3.7, ø 4.2, ø 6	ø 6
	(Color Code)	(Orange)	(Orange)
		ø 5.5 (Pink)	



Model	FB-52C-1	FB-52K-1	
	Rat tooth with alligator jaw (Swinging type)		
Shape of the Cups			
Maximum Insertion Portion Diameter (mm)	ø 1.8		
Working Length (mm)	1050	1550	
With Needle	NO		
Fenestrated Cups	YES		
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 600 mm; CHF, BF, ENF, CYF, HYF, OSF	Working length less than 1100 mm; EF, GIF, OC GF, GTF, C CF (Exclud I-, L-length) OSF, BF, E CYF, HYF
	Channel Inner Diameter (mm)	ø 2, ø 2.2 (Blue)	ø 2, ø 2.2 (Blue)
	(Color Code)	(Blue)	(Blue)
		ø 2.6 (Green)	ø 2.6 (Green)
	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3 (Yellow)	
		ø 3.7, ø 6 (Orange)	

Model	FB-53K-1	FB-53Q-1
	Rat tooth with alligator jaws (Swinging type)	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	1550	1950
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1100 mm; EF, GIF, OGF, GF, GTF, CF (Exclude I-, L-length), OSF
	Channel Inner Diameter (mm) (Color Code)	Working length less than 1500 mm; EF, GIF, OGF, GF, GTF, JF, TJF, CF (Exclude L-length), PCF(I-length only), OSF
	Channel Inner Diameter (mm) (Color Code)	Channel Inner Diameter (mm) (Color Code)
	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow)
	ø 3.7, ø 6 (Orange)	ø 3.7, ø 4.2, ø 6 (Orange)
		ø 5.5 (Pink)

Model	FB-53U-1	FB-54K-1
	Rat tooth with alligator jaws (Swinging type)	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	2300	1550
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1850 mm; EF, GIF, OGF, GF, GTF, JF, TJF, SIF (SIF-10 only), CF, PCF, OSF
	Channel Inner Diameter (mm) (Color Code)	Working length less than 1100 mm; EF, GIF, OGF, GF, GTF, CF (Exclude I-, L-length), O
	Channel Inner Diameter (mm) (Color Code)	Channel Inner Diameter (mm) (Color Code)
	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow)
	ø 3.7, ø 4.2, ø 6 (Orange)	ø 3.7, ø 6 (Orange)
	ø 5.5 (Pink)	

Model	FB-54Q-1	FB-54U-1
	Rat tooth with alligator jaws (Swinging type)	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	1950	2300
With Needle	NO	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1500 mm; EF, GIF, OGF, GF, GTF, JF, TJF, CF (Exclude L-length), PCF(I-length only), OSF
	Channel Inner Diameter (mm) (Color Code)	Working length less than 1850 mm; EF, GIF, OGF, GF, GTF, JF, TJF, SIF (SIF-10 only), CF, PCF, OSF
	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow)
	ø 3.7, ø 6 (Orange)	ø 3.7, ø 4.2, ø 6 (Orange)
		ø 5.5 (Pink)

Model	FB-55K-1	FB-55Q-1
	Rat tooth with alligator jaw (Swinging type)	
Shape of the Cups		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	1550	1950
With Needle	YES	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1100 mm; EF, GIF, OGF, GF, GTF, CF (Exclude I-, L-length), OSF
	Channel Inner Diameter (mm) (Color Code)	Working length less than 1500 mm; EF, GIF, OGF, GF, GTF, JF, TJF, CF (Exclude L-length), PCF(I-length only), OSF
	ø 2.8, ø 3.2 (Yellow)	ø 2.8, ø 3.2 (Yellow)
	ø 3.7, ø 6 (Orange)	ø 3.7, ø 4.2, ø 6 (Orange)
		ø 5.5 (Pink)

Model	FB-55U-1	
Shape of the Cups	Rat tooth with alligator jaws (Swinging type)	
		
Maximum Insertion Portion Diameter (mm)	ø 2.4	
Working Length (mm)	2300	
With Needle	YES	
Fenestrated Cups	YES	
Compatible Olympus Endoscopes (All of these parameters should be met.)	Model and Length	Working length less than 1850 mm; EF, GIF, OGF, GF, GTF, JF, TJF, SIF(SIF-10 only), CF, PCF, OSF
	Channel Inner Diameter (mm)	ø 2.8, ø 3.2 (Yellow)
	Channel Inner Diameter (mm)	ø 3.7, ø 4.2, ø 6 (Orange)
	Channel Inner Diameter (mm)	ø 5.5 (Pink)
Medical Device Directive		This device complies with the requirements of Directive 93/42/EEC concerning medical devices. Classification: Class I

Chapter 3 Preparation, Inspection and Operation

WARNING

- Before each case, prepare and inspect the instrument as instructed below. Inspect other equipment to be used with the instrument as instructed in their respective instruction manuals. Should the slightest irregularity be suspected, do not use the instrument. Contact Olympus. Damage or irregularity may compromise patient or user safety such as infection control risk, tissue irritation, punctures, hemorrhages, mucous membrane damage and may result in more-severe equipment damage.
- This instrument was not sterilized before shipment. Before using this instrument for the first time, reprocess it according to the instructions in Chapter 4, "Reprocessing". Do not use an instrument that has not been cleaned and sterilized. This poses an infection control risk or may cause tissue irritation.
- When using an instrument that has a Needle, be careful not to touch the Needle. Infectious substances attached to the Needle such as patient's blood or mucous, could pose an infection control risk and/or cause patient injury.

CAUTION

- Do not coil the Insertion Portion with a diameter of less than 15 cm. This could damage the Insertion Portion.
- Never use excessive force to open or close the Cups. This could damage the instrument.

3.1 Preparation

Equipment and Personal Protective Equipment

Prepare all equipment and personal protective equipment which will be used with the instrument in accordance with their respective instruction manuals. Appropriate personal protective equipment may include: Eye wear, a face mask, moisture-resistant clothing and chemical-resistant gloves.

Spare Instrument

Always have a spare instrument available.

Reprocessing Equipment

Prepare reprocessing equipment as described in Section 4.2, "Required Reprocessing Equipment" for immediate reprocessing after use.

3.2 Inspection

Wear the personal protective equipment as specified in Table on page 49.

Before each case, always inspect the instrument according to the following procedures.

If an abnormality in the instrument is detected, use a spare instrument, inspecting it thoroughly before use.

Inspection of the Sterile Package

Inspect the sterile package for tears, inadequate sealing or water damage. If the sterile package shows irregularities, the sterile condition of the instrument has been compromised. Use a spare instead.

Appearance Inspection

If any of the following steps reveals irregularities, do not use the instrument; use a spare instead.

1. While operating the Slider to open and close the Cups, confirm that the instrument is free from disconnection or looseness.
2. Make sure that the Cups close evenly and are properly aligned when the Slider is pulled.
3. When using an instrument with a Needle, push the Slider to open the Cups and confirm that the Needle is not detached or bent conspicuously.
4. Confirm that the Distal End of the instrument appears exactly as shown in the Tables in Section 2 "Specifications" and is not damaged.
5. Gently run your fingertips over the entire length of the Insertion Portion to check for any crushed areas, excessive bends, etc.

6. Make sure that there are no cracks on the Handle.

Inspection of Operation

If the Cups do not operate smoothly and as intended, do not use the instrument; use a spare instead.

1. Holding the instrument as shown in Figure 3.1, form a loop in the Insertion Portion approximately 20 cm in diameter.




Figure 3.1

2. Move the Slider and confirm that the Cups open and close smoothly.

3.3 Operation

The operator of the instrument must be a physician or medical personnel under the supervision of a physician and must have received sufficient training in clinical endoscopic technique.

This manual, therefore, does not explain or discuss clinical endoscopic procedures. It only describes basic operation and precautions related to the operation of this instrument.

WARNING

- When using the instrument, always wear appropriate personal protective equipment. Otherwise, blood, mucous and other potential infectious material from the patient could pose an infection control risk. Appropriate personal protective equipment may include: Eye wear, face mask, moisture-resistant clothing and chemical-resistant gloves that properly fit and are long enough so that your skin is not exposed.
- Do not insert the instrument into the endoscope unless you have a clear endoscopic field of view. If you cannot see the Distal End of the Insertion Portion in the endoscopic field of view or in X-ray images, do not use it. This could cause patient injury, such as punctures, hemorrhages or mucosal membrane damage. It may also damage the endoscope and/or instrument.
- Do not angulate the Bending Section of the endoscope (or operate the Forceps Elevator if applicable) abruptly while the Distal End of the Insertion Portion is extended from the Distal End of endoscope. This could cause patient injury, such as punctures, hemorrhages or mucosal membrane damage.

CAUTION

When using the instrument with a TV channel endoscope, never use electro-surgical accessories at the same time. This could cause patient, operator or assistant injury, such as thermal injury.

Inserting Into the Endoscope

WARNING

- Do not force the instrument if resistance to insertion is encountered. Reduce the angulation (or lower the Forceps Elevator if applicable) until the instrument passes smoothly. This could cause patient injury, such as punctures, hemorrhages or mucous membrane damage. It may also damage the endoscope and/or instrument.
- When inserting the instrument into the endoscope, hold the Slider firmly. Otherwise, the Cups may open and extend from the endoscope tip abruptly. This could cause patient injury, such as punctures, hemorrhages or mucous membrane damage. It may also damage the endoscope and/or instrument.
- Do not advance or extend the instrument abruptly. This could cause patient injury, such as punctures, hemorrhages or mucous membrane damage. It may also damage the endoscope and/or instrument.

CAUTION

- When inserting the instrument into the endoscope, hold it close to the Biopsy Valve or T-plug and keep it as straight as possible relative to the Biopsy Valve or T-plug. Otherwise, the Insertion Portion could be damaged.
- Before inserting the FB-51K-1 into the Biopsy Valve, open the Biopsy Valve. If you insert it without opening the Biopsy Valve, you may damage the Biopsy Valve and/or the instrument.

- When inserting the FB-51K-1 into the Biopsy Valve, hold it near the Distal End of the Insertion Portion to avoid damaging the instrument.

○ Combined With non-GF, non-GTF, non-TJF Endoscopes

1. Pull the Slider to close the Cups.
2. With the Cups closed, carefully insert the instrument into the Biopsy Valve or T-plug. (See Figure 3.2)




Figure 3.2

3. Advance the instrument until the Distal End of the Insertion Portion appears within the endoscopic field of view.

○ Combined With GF, GTF, JF or TJF Endoscopes

WARNING

When inserting the instrument into GF, GTF, JF or TJF endoscopes, raise the Forceps Elevator to its maximum height. If the Forceps Elevator is down, you will not be able to see the Distal End of the Insertion Portion in the endoscopic field of view. This could cause patient injury, such as punctures, hemorrhages or mucous membrane damage.

1. Raise the Forceps Elevator to its maximum height.
2. Pull the Slider to close the Cups.
3. With the Cups closed, carefully insert the instrument into the Biopsy Valve. (See Figure 3.2)
4. When the Distal End of the Insertion Portion contacts the Forceps Elevator, lower the Forceps Elevator.
5. Advance the instrument another 20 mm and raise the Forceps Elevator. You will see the Distal End of the instrument in the endoscopic field of view.

Collecting Tissue

WARNING

Do not force the Distal End of the Insertion Portion against body cavity tissue. This could cause patient injury, such as punctures, hemorrhages or mucous membrane damage.

1. To collect the target tissue, angulate the Bending Section or advance the instrument until it reaches the target site.
2. Push the Slider to open the Cups.
3. Press the open Cups against the target tissue.
4. Pull the Slider to collect the target tissue.

Withdrawing the Instrument From the Endoscope

WARNING

Do not withdraw the instrument from the endoscope quickly. This could scatter blood, mucous or other patient debris and pose an infection control risk.

CAUTION

- Do not withdraw the instrument from the endoscope while the Cups are open. This could damage the endoscope and/or instrument.
- If excessive resistance makes withdrawal difficult, adjust the angle of the endoscope until the instrument can be withdrawn smoothly. Forceful withdrawal could damage the instrument and/or endoscope.

1. If the endoscope is equipped with a Forceps Elevator, lower the Forceps Elevator.
2. Pull the Slider to close the Cups.
3. Withdraw the instrument from the endoscope.

Chapter 4 Reprocessing

WARNING

This instrument was not sterilized before shipment. Before using this instrument for the first time, reprocess it according to the instructions in this Chapter. Do not use an instrument that has not been cleaned and sterilized. This poses an infection control risk or can cause tissue irritation.

4.1 General Policy

- The medical literature reports incidents of patient cross contamination resulting from improper cleaning or sterilization. It is strongly recommended that reprocessing personnel have a thorough understanding of and follow all national and local hospital guidelines and policies. A specific individual or individuals in the endoscopy unit should be responsible for reprocessing endoscopic equipment. It is highly desirable that a trained backup be available should the primary reprocessing individual(s) be absent.
- All individuals responsible for reprocessing should thoroughly understand:
 - your institution's reprocessing procedures
 - occupational health and safety regulations
 - national and local hospital guidelines and policies
 - the instructions in this manual
 - the mechanical aspects of endoscopic equipment

- pertinent germicide labeling

Olympus Endo-Therapy Accessories are compatible with 2.0% to 3.2% glutaraldehyde solution. However, routine biological monitoring is not feasible with glutaraldehyde and, therefore, it should not be used to sterilize reusable medical devices that are compatible with other methods of sterilization that can be biologically monitored, such as steam sterilization.

WARNING

- Failure to properly clean and sterilize the instrument after each examination can compromise patient safety. During use, the instrument normally comes in contact with intact mucous membranes. To minimize the risk of transmitting diseases from one patient to another, after each examination the instrument must undergo thorough cleaning followed by sterilization.
- If the instrument is not cleaned meticulously, effective sterilization cannot be obtained. Clean the instrument thoroughly before sterilization to remove microorganisms or organic material which can limit the effectiveness of the sterilization process.

- Patient debris and reprocessing chemicals are hazardous. Wear personal protective equipment to guard against dangerous chemicals and infectious material. During cleaning and sterilization, always wear appropriate personal protective equipment, such as eye wear, a face mask, moisture-resistant clothing and chemical-resistant gloves that fit properly and are long enough so that your skin is not exposed. Always remove contaminated protective clothing before leaving the reprocessing area.
- The reprocessing procedures described in this manual should be completed the same day has been used. If reprocessing is delayed, residual organic debris will solidify and it may be difficult to effectively reprocess the instrument.

4.2 Required Reprocessing Equipment

Wear the personal protective equipment as specified in the Table on page 49.

1. Prepare the following equipment. The required amount of detergent solution, lubricant and other equipment depends on the number of instruments to be reprocessed.
2. Fill an immersion basin with detergent solution and fill a second immersion basin with lubricant at the temperatures and concentrations recommended by the manufacturers. Also fill the ultrasonic cleaner with a detergent solution appropriate for ultrasonic cleaning.

Equipment Needed for Reprocessing

To perform proper reprocessing, the equipment in the following Table is required. For details on preparation and directions for use of the following equipment, refer to the respective instruction manuals or contact the equipment manufacturer.

Contact Olympus for the names of specific brands of detergent solutions and lubricants.

Equipment Needed

Protective Equipment	Appropriate personal protective equipment may include: Eye wear, face mask, moisture-resistant clothing and chemical-resistant gloves.
Immersion Basin for Detergent Solution	Use a basin with a depth and diameter large enough to allow complete immersion of the instrument when the Insertion Portion is coiled with a diameter of not less than 15 cm.
Detergent Solution for Immersion	Use a neutral pH, low-foaming, medical grade detergent solution.
Ultrasonic Cleaner	Use a medical grade ultrasonic cleaner with a frequency range of 38 to 47 kHz, and with a depth and a diameter large enough to allow complete immersion of the instrument when the Insertion Portion is coiled with a diameter of not less than 15 cm. Compatible ultrasonic cleaners include OLYMPUS ULTRASONIC CLEANER KS-2
Detergent Solution for Ultrasonic Cleaning	Use a neutral pH, low-foaming, medical grade detergent solution with no abrasive.
Lubricant	Use a medical grade water soluble or low-viscosity emulsion type lubricant.

Immersion Basin for Lubricant	Use a basin with a depth and diameter large enough to allow complete immersion of the instrument when the Insertion Portion is coiled with a diameter of not less than 15 cm.
Lint-free Cloths	
Packages for Steam Sterilization	Use a packages compatible with steam sterilization (autoclaving). The packages should be large enough to accommodate the instrument when the Insertion Portion is coiled with a diameter of not less than 15 cm.
Sealing Device for packages	Sealing the packages may require a device such as a heat sealer. Prepare an appropriate searing device according to the packages to be used.
Autoclave	Use an autoclave that will operate at the conditions specified in Section 4.5, "Sterilization".

4.3 Cleaning

WARNING

When cleaning, avoid exposure to the reprocessing chemicals. It may pose an infection control risk or cause skin irritation.

CAUTION

- When reprocessing, do not coil the Insertion Portion with a diameter of less than 15 cm. This could damage the Insertion Portion.
- Never use excessive force to open or close the Cups. This could damage the instrument.

Immersion

WARNING

Immerse the instrument in detergent solution immediately after use. If the instrument is not cleaned immediately, it may be difficult to effectively reprocess and this could result in reduced performance.

1. When reprocessing the instrument for the first time after purchase, remove the Forceps Cap from the Cups and dispose of it.
2. Immerse the entire instrument in the detergent solution for the time specified in manufacturer instructions. If no time is specified, immerse between 5 minutes and 3 hours.
3. Remove the instrument from the detergent solution.

Ultrasonic Cleaning

1. Immerse the entire instrument in the ultrasonic cleaner containing detergent solution.
2. Clean ultrasonically for 30 minutes. For details of operation of the ultrasonic cleaner, refer to the instruction manual of the ultrasonic cleaner.
3. Remove the instrument from the detergent solution.

Rinsing

CAUTION

- After ultrasonic cleaning, rinse the instrument thoroughly to remove residual detergent. Residual detergent solution could cause tissue irritation in the next patient.

- Do not forcefully squeeze, wipe or scrub the instrument. This could cause damage to the instrument or result in reduced performance.

1. Rinse the instrument under clean running tap water.
2. Confirm that no debris is left on the surfaces of the instrument.
3. Wipe the exterior of the instrument with a clean, dry lint-free cloth.

4.4 Lubrication

WARNING

When lubricating, avoid exposure to the lubricant. It may pose an infection control risk or cause skin irritation.

CAUTION

- Do not coil the Insertion Portion with a diameter of less than 15 cm. This could damage the Insertion Portion.
- Never use excessive force to open or close the instrument. This could damage the instrument.

1. Immerse the Insertion Portion in the lubricant for 2 to 3 seconds.
2. Remove the instrument from the lubricant.
3. Operate the Slider to open and close the Cups two or three times.
4. Wipe the exterior of the instrument with a clean, dry lint-free cloth and allow the instrument to air dry.

4.5 Sterilization

Sealing the Package

WARNING

Before placing the instrument in package, always close the Cups. Otherwise they could tear the package during sterilization or storage and compromise its sterility, which could pose an infection control risk or cause tissue irritation.

CAUTION

- Do not coil the Insertion Portion with a diameter of less than 15 cm. This could damage the Insertion Portion.
- Never use excessive force to open or close the Cups. This could damage the instrument.

1. Before sterilization, the instrument must be thoroughly cleaned and dried. Residual moisture inhibits sterilization.
2. Coil the Insertion Portion and place the instrument in the package.
3. Seal the package. For details on sealing, refer to instruction manual of the package and the seal device.

Steam Sterilization (Autoclaving)

WARNING

- Use biological indicator recommended by your hospital policy and follow the manufacturer instructions, all national and local hospital guidelines, and policies.

- Always leave space between the packages in the autoclave. If the packages are placed too close together, effective sterilization will not be possible.
- Allow the packages to dry within the autoclave using the autoclave's drying cycle (if applicable) or by opening the door of the autoclave and allowing the packages to air dry. Handling a wet package can compromise its sterility.

1. Place the sealed package containing the instrument in the autoclave and sterilize in accordance with the conditions listed below. For details on operation of the autoclave, refer to the instruction manual for the autoclave or other manufacturer instructions.
2. After steam sterilization, let the instrument gradually cool down to room temperature. Sudden changes in temperature may damage the instrument.

NOTE Autoclavable products have a green reference label. Products that do not have green reference labels are not autoclavable.

	Temperature	Exposure Time
Prevacuum	132 to 134°C (270 to 274°F)	5 minutes

Table 4.1 Recommended Steam Sterilization(Autoclaving) Conditions

Chapter 5 Storage

WARNING

- Do not store the instrument in a sterile package that is damaged, wet, improperly sealed. Otherwise, the sterility of the instrument may be compromised and pose an infection control risk or cause tissue irritation.
- Do not store the sterile package containing the instrument in places where they will be damaged, wet, improperly sealed. Otherwise, the sterility of the instrument may be compromised and pose an infection control risk or cause tissue irritation.

CAUTION

Do not coil the Insertion Portion with a diameter of less than 15 cm. This could damage the Insertion Portion.

5.1 Inspection Before Storage

Prior to storage, inspect the sterile package as follows:

Confirm that the sterile package containing the instrument is free from tears, inadequate sealing or water damage. If tears, inadequate sealing or water damage is detected, repack and sterilize again as described in Section 4 "Sterilization".

5.2 Storage

Store the instrument in the sterile package at room temperature in a clean and dry environment. Do not store it in direct sunlight. Ensure that the packaged instrument is not crushed by surrounding objects during storage. Follow any additional storage instructions provided by the manufacturer of the sterile package.

OLYMPUS

OLYMPUS OPTICAL CO., LTD.

San-Ei Building, 22-2, Nishi Shinjuku 1-chome, Shinjuku-ku, Tokyo, Japan
Telex: 24209. Fax: (03)3340-2201/2. Telephone: (03)3340-2111.

OLYMPUS OPTICAL CO.(EUROPA) GMBH

(Premises/Goods delivery) Wendenstrasse 14-16, D-20097 Hamburg, Germany
(Letters) Postfach 10 49 08, D-20034 Hamburg, Germany. Telephone: (040)237730

OLYMPUS AMERICA INC.

Two Corporate Center Drive, Melville, N.Y.
11747-3157, U.S.A.
FAX: (516)844-5442. Telephone: (516)844-5000.

KEYMED LTD

KeyMed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, United Kingdom
Telex: 995283. Fax: (01702)465677. Telephone: (01702)616333.

OLYMPUS SINGAPORE PTE LTD

BLK 211, Henderson Road #13-03, Henderson Industrial Park, Singapore 159552
Telex: RS 35249 "OLYSIN". Fax: 2704704. Telephone: 2701351.

OLYMPUS BEIJING REPRESENTATIVE OFFICE

Room No. 3406 Beijing Jing Guang Center, Hu Jia Lou, Chao Yang Qu, Beijing, China
Fax: (10)6501-2085. Telephone: (10)6501-2084.

OLYMPUS MOSCOW CORPORATION

117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
Fax: (095)958-2277. Telephone: (095)958-2245.

OLYMPUS AUSTRALIA PTY. LTD.

1/104 Femtree Gully Road, Oakleigh, VIC 3166, Australia
Fax: (03)9543-1350. Telephone: (03)9265-5400.

