REF VC-10 (115VAC, 60 Hz)

Berkeley[®] VC[®]-10 Vacuum Curettage System

Instructions for Use and Technical Manual

ATTENTION: For your free copy of the VC-10 Setup Video, VHS or DVD format, please contact Customer Service at 1-888-524-7266 or 1-763-416-3000.

Gyrus ACMI, Inc.
136 Turnpike Road • Southborough, MA 01772-2104 • USA
1-888-524-7266 or 1-763-416-3000 • www.gyrusacmi.com

Table of Contents

Return/Repair/Warranty	3
Warnings and Cautions	3
Section 1.0 Unpacking and Initial Inspection	4
Section 2.0 General Description and Intended Use	4
Section 3.0 Description of Components	5
3.1 Accessory Components	5
3.2 Classifications	5
3.3 Technical Specifications	5
3.4 Nomenclature & Symbology	6
Section 4.0 Assemblies	6
4.1 Pump Assembly	6
4.2 Accessory Assembly	6
Section 5.0 Operating Instructions	7
5.1 Vacuum Operation	7
5.2 Vacuum Checks	7
5.3 Casters & Caster Locks	8
Section 6.0 Procedure	8
Section 7.0 Troubleshooting	9
Section 8.0 Environmental Protections (Health and Safety)	10
8.1 Disposal of Medical Waste Products	10
8.2 Disposal of Waste Electrical Equipment	10
Section 9.0 Maintenance	11
9.1 Fluid Overflow Protection/Safety Trap Collection System	11
9.2 Disposa-Filter	11
9.3 Fuse Replacement	11
Section 10.0 Cleaning & Disinfection	12
10.1 Cabinet	12
Section 11.0 Accessories	12

Return/Repair/Warranty

Return/Repair

All returns must have prior authorization. To comply with OSHA Bloodborne Pathogen Standard 29 CFR 1910.1030 regulations and U.S. Postal and Transportation law, all used medical devices returned for repair or replacement must be properly cleaned and decontaminated with a chemical germicide that has been cleared for use as a "Hospital Disinfectant." To ensure that the product has been properly decontaminated, a signed Decontamination Certificate should be enclosed in the package.

All medical devices returned to Gyrus ACMI for any reason must be shipped in accordance with Gyrus ACMI's return procedures (available upon request) and all applicable regulations. To obtain a return authorization number (RMA), return address, and instructions, please call Gyrus ACMI Customer Service at 1-888-524-7266 or 1-763-416-3000.

Please provide the product part number and lot/serial number. Product should be returned in its original packaging (when possible) and marked with the RMA number on the exterior of the package.

Limited Express Warranty

SHOULD THE PRODUCT BECOME INOPERABLE DURING NORMAL AND PROPER USE IN ACCORDANCE WITH APPLICABLE INSTRUCTIONS AND WITHIN THE TIME FRAME SPECIFIED BELOW FROM THE DATE OF SHIPMENT, GYRUS ACMI WILL REPAIR OR REPLACE THE PRODUCT, AT ITS SOLE OPTION, AT NO CHARGE. GYRUS ACMI MAKES NO OTHER WARRANTIES WITH RESPECT TO THE PRODUCT AND EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, AS TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER MATTER. IN NO EVENT SHALL GYRUS ACMI BE LIABLE FOR ANY CONSEQUENTIAL DAMAGES. IN NO EVENT SHALL GYRUS ACMI BE LIABLE FOR ANY BREACH OF WARRANTY IN ANY AMOUNT EXCEEDING THE PURCHASE PRICE OF THE PRODUCT.

This warranty applies only to the original purchaser and will be voided if the product(s) are serviced or repaired by anyone other than Gyrus ACMI or an organization duly authorized by Gyrus ACMI for such purpose.

Vacuum Curettage Systems..... 1 year

Warnings and Cautions

The following warnings and cautions apply to the use, care, and/or maintenance of the Gyrus ACMI® Berkeley VC -10 Vacuum Curettage System.

Warnings

(indicate that a danger to life or health can result from misusing the equipment)

1. Other than user/operator accessible and replaceable components identified and described within this document, there are NO serviceable parts or assemblies in the VC-10 Vacuum Curettage System. For technical support, contact the Gyrus ACMI Technical Support Hotline at 1-800-331-2264 or 1-763-416-3000.
2. Study this manual and other labeling thoroughly for safe handling and storage. Misuse of instruments can cause injury to the patient and could have an adverse effect on the procedure being performed. Do not drop instruments or allow them to be struck by other objects.
3. Procedures should be performed only by persons with adequate training and preparation. Consult the medical literature regarding techniques, complications, and hazards prior to any vacuum curettage procedure.
4. This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with nitrous oxide (N₂O).
5. This equipment is not suitable for use in hazardous atmospheric conditions.
6. Grounding reliability can only be achieved when the unit is connected to equivalent receptacle marked "Hospital Grade".
7. Do not use an instrument that fails to meet the criteria stated in the labeling, or that has been damaged.
8. Follow the labeling instructions regarding the disposal or reuse of accessories. Reuse of disposable accessories could compromise patient safety.
9. Use only those replacement parts specified in the labeling.

Cautions

(indicate that equipment or other property may be damaged or may malfunction by misuse)

1. To ensure satisfactory performance, perform the prescribed inspections and operational tests as recommended.
2. Do not allow fluids to enter the instrument cabinet.
3. To help avoid leaks after cleaning the fluid-overflow traps, before replacing the removable plastic bowls, confirm visually that the cork gasket in the underside of each bowl-mount flange is undamaged.
4. Do not move the device by pulling on the power cord or the foot pedal connecting hose.
5. The unit contains an internal thermal overheat mechanism to protect the motor from damage if the unit should overheat. Good ventilation and normal ambient temperatures (65–85 °F) during use should prevent the unit from overheating. If the unit is operated continuously for over 1 hour, the unit could overheat and shut down. If the unit shuts down, turn the unit OFF and allow the motor to cool. The motor will start again when the temperature has been reduced sufficiently.
6. Run the unit for approximately 5 minutes after each procedure (with the system open to atmosphere) to minimize moisture build-up that can lead to corrosion in the vacuum pump.
7. The VC-10 is NOT for field or transport use.

Section 1.0 Unpacking and Initial Inspection

Upon receipt, examine the instrument and associated accessories for damage. Do not use a damaged product or accessory. If there is any damaged product or accessory, contact Gyrus ACMI Customer Service at 1-888-524-7266 or 1-763-416-3000.

For technical support, contact the Gyrus ACMI Technical Support Hotline at 1-800-331-2264 or 1-763-416-3000.

The Berkeley VC-10 Vacuum Curettage System is shipped as an integral unit that includes the vacuum curettage pump, foot pedal, and power cord within its shipping carton. A starter kit of disposable accessory components is also provided in a separate carton (and may be placed within the unit storage compartment). Unpack the cartons carefully to avoid damage to the components and loss of any parts.

Section 2.0 General Description and Intended Use

The Berkeley VC-10 Vacuum Curettage System consists of a vacuum pump, collection canisters, connecting tubing and filters. A Vacurette® cannula (semi-rigid: straight and curved; or flexible) is attached to a collection set (tubing and handle) and attached to the patient port of the first collection canister lid.

The VC-10 system is equipped with dual collection and safety trap systems with automatic vacuum restriction. The system comes complete and ready to operate. The system features elastomeric polyurethane wheels (two of four with locking feature), a cord wrap on the back panel, and a storage compartment for supplies and accessories (accessible from the front door of the unit).

The VC-10 system is a high performance model that meets hospital operating room safety requirements. The device construction is a formed & welded steel structure with thermoformed plastic panels. All plastic panels and sheet metal components are painted with durable urethane paint. Polyurethane bumpers are added in four lower corners for impact protection.

The Berkeley VC-10 Vacuum Curettage System is designed for rapid transcervical aspiration of the contents of the uterus.

Section 3.0 Description of Components

3.1 Accessory Components

<u>Quantity</u>	<u>Description</u>
6	Berkeley SafeTouch® Collection Canisters, single-use disposable (for 3/8" tubing sets)
5	Berkeley SafeTouch Tissue Traps
10	Vacurette Cannulas, 8 mm, curved, sterile, single-use, disposable
5	18" Series Connecting Hose
5	3/8" Disposable Collection Tubing Set, sterile, single-use, disposable
1	Disposa-Filter Assembly
1	Berkeley VC-10 In-Service Video
1	Instructions for Use manual

See Section 11.0 for complete list of Berkeley Vacuum Curettage accessory components.

3.2 Classifications

Type Construction: Mobile Equipment

Protection Against Electric Shock Hazards—

Class of Protection: Class-I

Degree of Protection: Type-B

Rated Operation: Continuous

3.3 Technical Specifications

(for reference only—typical measurements/performance)

Operating Voltage—

VC-10: 115 VAC, 60 Hz

Power Consumption: 400VA

Fuse Ratings & Type—

VC-10: Qty:(2) (T6.3AH 250V)

Gyrus ACMI REF: P2700563-001

Earth/Chassis Leakage Current: Less than 150 μ a

Ground Resist: Less than 0.1 ohms case-to-ground connection

Horsepower Rating: 1/3 hp

Dimensions: 18.3" W x 25" D x 38" H
(46.4 cm W x 65.3 cm D x 96.5 cm H)

Weight: 110 lbs.

Vacuum Indicator Accuracy: Within \pm 5% of full-scale value

Environmental Specifications—

Operating Temperature: 5–40 °C, 70 kPa to 106 kPa atmospheric pressure

Storage Temperature: 0–50 °C, 60 kPa to 106 kPa atmospheric pressure

Relative Humidity: 10–90% non-condensing

Volume Reference: 1 cc (cubic centimeter) = 1 ml (milliliter)

1.	2.	3.	4.
	I O		
5.	6.	7.	
		Rx Only	

3.4 Nomenclature & Symbology

1. Type B Equipment
2. On/Off
3. Protective Earth Ground
4. Foot Pedal Actuator Control
5. Attention: Consult Manual
6. Quantity symbol: numeral in symbol (in place of #) indicates the quantity of units in package.
7. Caution: Federal (US) Law restricts this device to sale by or on the order of a physician.

Section 4.0 Assemblies

4.1 Pump Assembly

1. Verify that the line voltage rating shown on the back panel corresponds to the available power and that the power receptacle to be used is grounded:
 - a. VC-10 (115 VAC, 60 Hz)
2. Attach the appropriate power cord to the back of the unit (female end into the power entry module, male end into the power outlet) based on the available power as determined in #1 above.
3. If foot pedal control is desired, attach the foot pedal by fully inserting the foot pedal hose connector to its mating bulkhead receptacle on the front panel of the unit by pushing it straight in until one click sound is heard, confirming correct “full-engagement”. (It will be necessary to open the front door to gain access to the mating receptacle.) If the foot pedal connector cannot be inserted into the bulkhead receptacle, it may be necessary to fully depress the “thumb-slide” of the bulkhead connector downward, which will allow the insertion of the foot pedal connector into the bulkhead receptacle (see Figure 1).

Figure 1

4.2 Accessory Assembly

For 3/8" collection tubing set (REF 23116), see Figures 2 & 3:

1. Place a 3/8" primary and secondary collection canister and lids (REF 003987-901 or 003984-901) in the canister holders/receptacles formed into the top cover of the unit. Confirm that the canister lids are firmly in place.
2. Attach the short tubing end of the Disposa-Filter (B) (REF 54298) to the vacuum port on the pump (A) and the long tubing end to the “PUMP” port on the secondary collection canister lid (C).
3. Attach a series connecting hose (D) (REF 20714) to the “PATIENT” port on the second collection canister lid (E) and to the “PUMP” port on the primary collection canister (F).
4. Insert either a SafeTouch tissue trap (G) into the “PATIENT” port on the primary collection canister lid (H), or attach a permeable gauze sack (I) (REF 003997-901) to the underside of the “PATIENT” port on the primary collection canister.
5. To prevent the secondary canister from filling prematurely, an aspiration filter (J) (REF 003905-901) may also be attached to the underside of the “PUMP” port on the primary collection canister lid.
6. Attach a sterile 3/8" collection tubing set (K) to top of the tissue trap (L) or to the “PATIENT” port of primary collection canister lid (H).
7. Select and attach a sterile 7–12 mm Vacurette cannula (or 14 mm with adapter REF 003696) to the patient end of the collection tubing set (M).
8. Close off all open ports with sealing caps.

NOTE: For assembly of the 1/2" collection tubing set (REF BKT-506), refer to Gyrus ACMI document 181058.


Figure 2

Figure 3

Figure 4a – OFF Position

Figure 4b – ON Position


Figure 4c – Foot Pedal ON Position

Section 5.0 Operating Instructions

5.1 Vacuum Operation

The VC-10 power switch is a 3-position switch. The unit is shipped with the switch in the OFF position, which is the central position for this switch (see Figure 4a).

For vacuum operation to occur, the power switch must be placed in the ON position. The power switch is set in the ON position by fully depressing the left edge of the switch (see Figure 4b).

To turn the VC-10 device OFF at any time during the procedure or at the end of the procedure, return the power switch to the OFF position (central position for the 3-position switch).

For foot pedal operation, confirm that the foot pedal is properly attached. Note: Close the front door prior to using the VC-10 device. Fully depress the right edge of the power switch and then depress the foot pedal (refer to Section 3.4 Nomenclature and Symbology, and Figure 4c). Note: Releasing the foot pedal will discontinue the suction.

For intermittent suction, fully depress the foot pedal when suction is desired. If depressing the foot pedal does not result in suctioning, it may become necessary to vent the system by utilizing the “slide vent” of the sterile collection tubing set to allow the foot pedal actuation to continue.

Alternatively, if the Vacurette tip becomes clogged, clear clogged Vacurette tips as necessary and continue suctioning. If suctioning does not continue, see Section 7.0, #1a.

The maximum attainable vacuum with the vacuum-adjust valve in the fully clockwise position, at sea level, is approximately 71 cm Hg. There is a reduction of vacuum by 2.6 cm Hg per 1,000 feet (8.5 mm Hg per 100 meters) of elevation above sea level.

The vacuum-adjust knob is pre-set at the factory in the fully counterclockwise position for minimum vacuum. To increase vacuum, turn the vacuum-adjust knob clockwise (see Figure 5). To determine the maximum vacuum level available at any particular setting, turn the unit on and observe the vacuum reading on the gauge while completely occluding the intake opening of the collection bottle.

Vacuum Gauge Color Range (at outer dial “tick” marks)

0 to 7.4 cmHg:	White
7.5 to 50 cmHg:	Red
50 to 60 cmHg:	Yellow
60 to 75 cmHg:	Green

Vacuum-tight connections are assured when the collection system tapered fittings are properly connected and maintained. Proper sealing will maintain a consistent vacuum level throughout the system.

To minimize moisture build-up that can lead to corrosion in the vacuum pump, utilize the appropriate Disposa-Filter (see Figure 2) and monitor canister collections and related overflow protection/safety trap collection systems. Running the unit for approximately 5 minutes after each procedure (with the system open to atmosphere) will also minimize moisture build-up.

Figure 5

5.2 Vacuum Checks

1. Turn the VC-10 unit on, either by setting the power switch to the ON position (Figure 4b) or to the foot actuator position for foot pedal operation (Figure 4c).
2. Place a finger over the inlet at the end of the collection tubing handle. Observe the vacuum gauge until it stabilizes. The level shown on the gauge is the maximum vacuum level at the particular setting of the vacuum adjust knob.
3. Turn the vacuum adjust knob counterclockwise to decrease vacuum or clockwise to increase vacuum.
4. Verify that the vacuum gauge has stabilized after each knob adjustment before using the system in a procedure.
5. Slide the slip-ring on the rotating handle of the collection tubing assembly over the opening to maximize the suction to the tip of the Vacurette cannula. Slide the slip-ring away from the opening to reduce the suction at the tip of the Vacurette cannula.

Figure 6a – Unlocked

Figure 6b – Locked

5.3 Casters & Caster Locks

The VC-10 device has four casters (two locking casters on the front of the device, and two rolling casters without locking features on the rear of the device).

- a. Casters are “unlocked” when the locking lever is in the upward position (see Figure 6a).
- b. To lock casters, depress the locking lever of the caster with your foot (see Figure 6b).
- c. To unlock a locked caster, lift up the locking mechanism by foot to the unlocked position (see Figure 6a).

Section 6.0 Procedure

1. Assemble and connect appropriate accessories (see Section 4.0).
2. Dilate the cervix as needed and insert the cannula.
3. The sliding ring on the handle may be used to control the vacuum. If the cannula adheres to the uterine wall or becomes clogged, release the suction by moving the sliding ring from the closed position to the open position to decrease the vacuum prior to removal of the cannula.

WARNING: Do not apply excessive force while using the cannulas and F-sets. Excessive pushing, bending, or pulling could result in injury to the patient and/or breakage of the device.

4. At the end of the procedure discard the Vacurette cannula or F-set, disposable collection tubing set, and SafeTouch collection system components (including filters and tubing). These items are intended for single use only.

WARNING: Follow all applicable bloodborne pathogen procedures as indicated by OSHA and/or hospital/institutional requirements when disposing of single-use products. For further information, see Section 8.0 Environmental Protection (Health and Safety).

Section 7.0 Troubleshooting

Symptom/Possible Cause/Suggested Remedy

1. Unit does not start

- a) No power
 - Verify that the power cord is fully inserted into the VC-10 unit and also fully plugged into a properly working outlet.
 - Check fuses and replace as needed (see Section 3.3 for correct fuse type and part number).
- b) ON/OFF switch actuation or ON/OFF/Foot Pedal switch
 - Verify that the switch is in the ON position or verify that the switch is in the Foot Pedal position and that the foot pedal is properly connected to the unit and is depressed.
- c) Worn power cord
 - Check plug and cord for excessive wear. Replace as needed.
- d) Motor overheated
 - Allow motor to cool before continuing procedure.

2. Insufficient vacuum

- a) Improper setting
 - Adjust knob clockwise.
- b) Check vacuum
 - See Section 5.2.
- c) Disposa-Filter blocked
 - Replace as needed.
- d) Loose connections on disposable components
 - Tighten tubing connections, sealing caps, and lids on canisters.
- e) Tube couplings not connected or fully engaged
 - Confirm all tubing coupling connections behind rear panel.
- f) Collection canister bottle cracked
 - Replace as needed.
- g) Tubing blocked and/or leaking
 - Replace as needed.
- h) Slip ring on Vacurette handle open
 - Move slip ring over opening.
- i) Vacurette cannula blocked or leaking
 - Replace as needed.
- j) Safety trap(s) not properly installed or clogged
 - See Sections 4.2 and/or 9.1.

Section 8.0 Environmental Protections (Health and Safety)

8.1 Disposal of Medical Waste Products

The disposal of medical waste products (disposable devices and/or their contents) creates a risk similar to the risk associated with handling, storage, and removal of potentially contaminated medical waste products generated by a waste producing facility (e.g., hospital, laboratory, clinic, physicians office, or other similar facility).

It is the sole responsibility of the waste producing facility to properly obtain, understand, and follow all compliance guidelines, applicable policies and procedures, statutes, and regulations that govern regulated medical waste produced by their facility in conjunction with all applicable local, city, county, town, state, or country specific requirements for medical waste management, including determination and characterization of bio-hazardous medical waste (whether contaminated and/or infectious) which may include solid or liquid human waste (tissue, organs, products of conception, and body fluids) and proper instructions for handling, precautionary disinfection or decontamination (if required), containment and storage, marking or labeling, waste accumulation and duration of storage requirements, and waste removal and disposal procedures.

In addition, it is the sole responsibility of the waste producing facility to provide a bio-safe working environment, including training of staff members who come in contact with medical waste regarding appropriate policies and procedures, use of protective equipment, selection of storage containers and related supplies that prevent the release of waste prior to ultimate disposal, and any record keeping activities relative to the source, type of waste, and disposal history.

8.2 Disposal of Waste Electrical Equipment

It is the sole responsibility of the user/user facility to properly dispose of equipment in a safe manner consistent with all compliance guidelines, applicable policies and procedures, statutes, and regulations that govern waste electrical equipment discarded by their facility in conjunction with all applicable local, city, county, town, state, or country specific codes for waste management.

The safe disposal of waste electrical equipment should be undertaken only at the end of the equipment's useful life (i.e., the device is nonfunctional & nonrepairable), and should be undertaken in an environmentally conscious manner minimizing environmental impact.

It is the sole responsibility of the user/user facility to determine whether waste electrical equipment should be collected and disposed separately from or together with unsorted municipal waste. Public or private solid waste collection and disposal facilities that offer various forms of recovery (dismantling, reuse, and recycling, where possible) should be selected reducing the amount of waste being discarded and providing the highest level of protection for human health and the environment.

It is the sole responsibility of the user/user facility to determine the conditions under which the waste electrical equipment is to be collected, noting that items contaminated by the user may require confirmed decontamination by means of the best available treatment options.

Section 9.0 Maintenance

Figure 7

9.1 Fluid Overflow Protection/Safety Trap Collection System

The dual fluid overflow protection/safety trap collection system, located behind the rear panel of the cabinet, prevents fluids from entering, damaging, and/or contaminating the pump (see Figure 7).
NOTE: When the first overflow collection jar is full, the float ball blocks airflow from the vacuum pump and “suction” capabilities cease (while the vacuum pump may still be in the ON position).

The traps should be checked and cleaned periodically, and especially whenever any fluid reaches the Disposa-Filter. The unit should NOT be operated with any fluid in either of the fluid-overflow collection jars.

1. Hold a small pan underneath the fluid trap collection jar before unscrewing the jar to catch any fluid spillage due to jar removal.
2. Discard the contents per institutional procedures for disposing of bio-hazardous materials (refer to section 8.0 for general health and safety concerns).
3. Clean float ball, float ball mechanism, and collection jar thoroughly with an enzymatic detergent diluted according to the manufacturer’s recommendations.
4. Check the trap body for blockage of port openings.
5. Check the jar for cracks and leaks. Check gasket for wear. Replace worn or damaged parts; contact Technical Support or Customer Service for device repair or for part ordering information.
6. Dry all components well before reassembling, making sure to reinsert the float ball into the collection jar before reassembly.

Note: In the event that the second overflow collection jar (on the left when viewed from the rear of the unit) has fluid in it as well, repeat steps 1 through 6 for this jar.

9.2 Disposa-Filter

Replace the Disposa-Filter if it becomes soiled, clogged, or discolored/yellowed. To remove the Disposa-Filter, detach the connecting hose from the second collection canister and pull from the male connection/fitting on the top cover of the cabinet. Replace with new Disposa-Filter.

9.3 Fuse Replacement


Figure 8

1. Before replacing fuses, **disconnect the power cord from the power entry module** (lower rear panel of the unit).
2. To access the fuses, pry open and remove the small drawer/fuse holder at the bottom of the power entry module with a small, straight blade screwdriver or stylus. Remove fuses from the drawer/fuse holder (see Figure 8).
3. It is recommended to replace both fuses when fuse replacement is required. Replace fuses (Qty: 2) as necessary with correct fuse type, based on unit voltage (see Section 3.3 Technical Specifications for fuse ratings, type, and Gyrus Medical part number). Install fuses by placing them into the drawer/fuse holder (there is no specific orientation of the fuses required).
4. Fully reinstall the drawer/fuse holder into the power entry module. Reconnect the power cord into the power entry module. Confirm that the device is operational by plugging the unit into a functional wall outlet and turning the power switch to the ON position.

Section 10.0 Cleaning & Disinfection

The following cleaning and disinfection procedures are recommended by Gyrus ACMI. Use of any procedures not expressly recommended by Gyrus ACMI may adversely affect or damage Gyrus ACMI devices and may void Gyrus ACMI's Limited Express Warranty.

Follow all applicable bloodborne pathogen procedures as indicated by OSHA and/or your hospital requirements when cleaning, disinfecting, and/or sterilizing instruments.

10.1 Cabinet

The cabinet may be cleaned with a soft cloth or sponge dampened with a surface-cleaning solution.

Section 11.0 Accessories

Berkeley SafeTouch Collection System

The Berkeley SafeTouch Collection System is a completely disposable, closed containment system for use with the vacuum curettage system. The use of this system reduces the risk of exposure to bloodborne pathogens. It is available with or without a sealed tissue trap that eliminates handling the specimen. The collection system is sold in packages of 10 and comes complete with sealing caps.

REF	Description	Quantity
<u>Berkeley SafeTouch Collection System and Accessories</u>		
003984-901	3/8" Berkeley SafeTouch Collection System with Tissue Trap. Includes 10 ea. collection canisters, lids, tissue traps, and sealing cap sets.	1
003987-901	3/8" Berkeley SafeTouch Collection System without Tissue Trap. Includes 10 ea. collection canisters, lids, and sealing cap sets.	1
BKC-506	1/2" Berkeley SafeTouch Collection System. Includes 10 ea. 1/2" collection canisters, lids with preattached permeable gauze sacks with O-rings, and sealing cap sets.	1
003905-901	Berkeley SafeTouch Aspiration Filters	10/pkg
003853-902	Berkeley SafeTouch Tissue Traps	10/pkg
003997-901	Permeable Gauze Sacks with O-rings (for SafeTouch Bottles only)	10/pkg
101716-2	Large Blue Sealing Cap	1
101716-1	Small Blue Sealing Cap	1
003991-501	Blue Sealing Caps (3 large & 2 small)	1
<u>Disposable Collection Tubing Sets</u>		
23116‡	3/8" Disposable Collection Tubing Set	10/pkg
BKT-506‡	1/2" Disposable Collection Tubing Set	10/pkg
<u>VC Accessories</u>		
003696‡	14 mm Adapter	10/pkg
20714	18" Series Connecting Hose	10/pkg
54298	Disposa-Filter Assembly	10/pkg
54912	Disposa-Tops (for 3/8" tubing set)	10/pkg
56134-2	Disposa-Tops (for 1/2" tubing set)	10/pkg

‡ Supplied sterile, single-use, disposable

Berkeley Vacurette Cannulas

Berkeley Vacurette Cannulas are semi-rigid (straight and curved) or flexible (F-tip) vacuum aspiration currettes with rounded tips. They are used with the Berkeley Disposable Collection Tubing Sets (REF 23116 and REF BTK-506, sold separately) and the Berkeley Vacuum Curettage System.

Vacurette sizes 7–12 mm connect to 3/8" tubing (REF 23116) and Vacurette sizes 14 & 16 mm connect to 1/2" tubing (REF BTK-506). Vacurette size 14 mm can be used with 3/8" tubing with use of the 14 mm adapter (REF 003696)

REF	Style	Size	REF	Style	Size	REF	Style	Size
21852	straight	7 mm	21853	curved	7 mm	21663	F-tip	4 mm
21655	straight	8 mm	20317	curved	8 mm	21664	F-tip	5 mm
21413	straight	9 mm	21552	curved	9 mm	21665	F-tip	6 mm
21414	straight	10 mm	21553	curved	10 mm	21744	F-tip	7 mm
21415	straight	11 mm	21554	curved	11 mm	21745	F-tip	8 mm
21416	straight	12 mm	21555	curved	12 mm			
21464A	straight	14 mm	21593A	curved	14 mm			
21561	straight	16 mm	21594	curved	16 mm			

The F-Sets consist of an F-tip Vacurette cannula pre-attached to a collection tubing set.

REF	Style	Size	REF	Style	Size
23135	F-set	6 mm	23157	F-set	7 mm

The Berkeley Collection Tubing Sets, Vacurette Cannulas, and F-Sets are supplied sterile, single-use, disposable in individually sealed pouches, 10 per box.

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