

MicroChoice® Handpieces, Guards and Attachments



	Cat. No.
MicroChoice Elite™ High Speed Drill.....	5020-025
Electric	
Stainless Steel and Aluminum	
Max Standard Operating Speed 100,000rpm	
Torque 2in-oz	
Weight 7oz (196g)	
For use with Advantage®, E9000® and PowerPro® Consoles	
MicroChoice Medium Speed Drill	5020-021
Electric	
Stainless Steel and Aluminum	
Max Standard Operating Speed 25,000rpm	
Torque 6in-oz	
Weight 7oz (196g)	
For use with Advantage, E9000 and PowerPro, Consoles	
MicroChoice Medium Bur Guard	5020-060
Stainless Steel	
Use with 5091 Series Burs	
For use with 5020-025 and 5020-021 Handpieces	
MicroChoice Long Bur Guard	5020-061
Stainless Steel	
Use with 5092 Series Burs	
For use with 5020-025 and 5020-021 Handpieces	
MicroChoice Extra Long Bur Guard	5020-062
Stainless Steel	
Use with 5093 Series Burs	
For use with 5020-025 and 5020-021 Handpieces	
MicroChoice Laminectomy Bur Guard	5020-067
Stainless Steel	
Use with 5092-203 Bur Only	
For use with 5020-025 and 5020-021 Handpieces	
MicroChoice Tissue Retractor	5020-068
Stainless steel	
Use with 5092-236 Bur Only	
For use with 5020-025 and 5020-021 Handpieces	
MicroChoice 20° Angle Attachment.....	5020-063
Stainless Steel	
Use with 5092 Series Burs	
Use 1375-037 Spray Lubricant	
MicroChoice 20° X-Long Angle Attachment	5020-064
Stainless Steel	
Use with 5093 Series Burs	
Use 1375-037 Spray Lubricant	
MicroChoice Bur Changer	1375-003
Plastic	
Used for Inserting Wire Pass Burs into 70° and 90° Angle Attachments.	
MicroChoice 70° Contra Angle Attachment	5020-065
Stainless Steel	
Use with 5089 and 5090 Burs	
For use with 5020-025 and 5020-021 Handpieces	
Use 1375-037 Spray Lubricant	
MicroChoice 90° Angle Attachment.....	5020-066
Stainless Steel	
Use with 5089 and 5090 Burs	
For use with 5020-025 and 5020-021 Handpieces	
Use 1375-037 Spray Lubricant	
MicroChoice Dental 70° Contra Angle Attachment.....	5020-069
Stainless Steel	
For use with 5020-025 and 5020-021 Handpieces	
Use 1375-0374 Spray Lubricant	
Uses Friction Grip Burs Only. Not Available from Linvatec	