



# Spine Surgery

U Series Elite Attachments



Dedicated Saws



Integrated Irrigation



- Specialty-driven handpieces
- Increased power, precision burring
- Complete OR standardization

# Spine Surgery



## U Drill

- 75K RPM, Electronic Torque Feedback
- Tapered attachments, Telescoping burs for line of sight
- Graft Harvester



## Sagittal Saw

- New Rotating, Indexing Head
- 40% More Power than TPS Sagittal Saw
- New tighter retention blade mount



## Universal Driver

- New Trigger Oscillating Feature
- 40% More Torque than TPS
- 100% More Power than TPS



Rotating,  
Indexing Head



### CORE Spine Ordering Information

PART #	DESCRIPTION	PART #	DESCRIPTION
<b>5400-100</b>	<b>CORE U-Drill</b>	5100-10-334	Sagittal Saw Attachment
5400-9	CORE Handswitch		
5100-10-20	U Series O Straight Attachment	<b>5400-34</b>	<b>CORE Sagittal Saw</b>
5100-10-22	U Series O Angled Attachment	<b>5400-31</b>	<b>CORE Oscillating Saw</b>
5100-10-210	Irrigation Clip for O Attachment	<b>5400-37</b>	<b>CORE Reciprocating Saw</b>
5100-10-50	U Series M Straight Attachment	<b>5400-99</b>	<b>CORE Universal Driver</b>
5100-10-52	U Series M Angled Attachment	4100-125	Pin Collet
5100-10-200	Irrigation Clip for M/AM Attachments	4100-62	Wire Collet
5100-10-70	U Series AM Straight Attachment	4100-110	Synthes® Drill
5100-10-72	U Series AM Angled Attachment	4100-131	1/4" Jacobs® Drill
5100-10-80	U Series R Attachment	4100-132	5/32" Jacobs® Drill
5100-10-90	U Series RX Attachment	4100-133	1/4" Keyless Jacobs® Drill
		4100-135	Hudson® Modified Trinkle Drill
5100-10-450	U Series Elite Medium Straight	4100-160	Trinkle Reamer
5100-10-452	U Series Elite Medium Angled	4100-210	Synthes® Reamer
5100-10-270	Medium Irrigation Clip	4100-231	Jacobs® Reamer
5100-10-470	U Series Elite Long Straight	4100-235	Hudson® Modified Trinkle Reamer
5100-10-472	U Series Elite Long Angled	4100-260	Trinkle Reamer
5100-10-272	Long Irrigation Clip	4100-400	Sagittal Saw Attachment
5100-10-480	U Series Elite X-long Straight	4100-700	High Speed Bur Attachment
5100-10-482	U Series Elite X-long Angled		
5100-10-490	U Series Elite MIS		



For complete CORE System ordering information, contact your local Stryker Sales Representative.

**stryker®**  
MedSurgGroup

See Bur and Blade literature for complete listing of available cutting accessories

**Instruments Division**  
4100 East Milham Ave  
Kalamazoo, MI 49001  
(800) 253-3210

**Endoscopy Division**  
5900 Optical Court  
San Jose, CA 95138  
(800) 624-4422

**Leibinger Division**  
4100 East Milham Ave  
Kalamazoo, MI 49001  
(800) 962-6558