

# Venture™ 227

Economy Wheelchair

SEAT WIDTH AVAILABLE IN: 18"  
UPHOLSTERY AVAILABLE IN: BLACK



*The Standard Economy 227 wheelchair is the ideal choice for those who desire a quality wheelchair at a economical price. Equipped with ABS components, powder coated frame, precision seal bearings at all key locations for durability & smooth rolling, it still maintains unsurpassed workmanship and manufacturing standards. Quality of this model is never compromised with the low cost.*

**Tuffcare®**  
HIGHER STANDARD FOR HEALTH CARE™

# Venture™ 227

Economy Wheelchair

SEAT WIDTH AVAILABLE IN: 18"  
UPHOLSTERY AVAILABLE IN: BLACK

## STANDARD FEATURES

- \* Durable moisture resistant powder coated frame
- \* Fixed full ABS armrests
- \* 24" molded rear wheels with solid rubber tires
- \* Precision sealed bearings in rear wheels, caster wheels & caster stems
- \* High impact composite footplates
- \* High impact composite side panels
- \* High impact composite handrims
- \* 8" molded casters with solid rubber tires
- \* Cam-action front riggings: detachable swingaway footrests/elevating legrests
- \* Embossed, flame retardant, durable vinyl upholstery

## SPECIFICATIONS

MODEL	SEAT			BACK	FRONT RIGGINGS	OVERALL					WEIGHT (w/o Front Riggings)
227	18"	20"	16"	16 3/4"	16"- 21"	24 1/2"	36"	41"	31"	10 1/2"	36 lbs.
227E	18"	20"	16"	16 3/4"	16"- 21"	24 1/2"	36"	44 1/2"	31"	10 1/2"	36 lbs.
E	ELEVATING LEGRESTS										

## OPTIONAL ACCESSORIES

- \* Anti-theft bar
- \* Wheelchair tray
- \* Amputee adapter
- \* Headrest extension
- \* Wheel lock extension tubes
- \* One arm drive
- \* Anti-theft lock for front riggings
- \* Pneumatic tires
- \* I.V. rod attachment
- \* Anti-tipping device
- \* Anti-folding device
- \* Seat belt
- \* Projection handrim
- \* Leg straps
- \* Oxygen tank holder
- \* Pneumatic tires with airless inserts



Stylish composite handrim eliminates uncomfortable cold metal feel in winter time.



Swingaway detachable footrests or elevating legrests



Precision seal bearing in all key locations provides smooth noise free operation.