

Power  
Examination  
Table

# 405

---



**MIDMARK®**

# 405

## Power Examination Table

The Midmark 405 Power Examination Table is packed with special features for general practitioners, obstetricians, gynecologists, dermatologists and plastic surgeons. It features a tapered, powered back section for improved access to the head and neck area of patients. And it comes standard with an ergonomic hand control which effortlessly raises the back section from flat to a 75° chair. An 'Auto Return' feature brings the table to its lowest position to save valuable time; a 'Stop' button overrides the 'Auto Return'. For your convenience, a foot control accessory, which can be plugged in at either the head or foot end of the table, can also be added.

With a low height of 25", the 405 is more accessible to infirm, elderly and pregnant patients, as well as to those patients in wheelchairs. And once patients are on the table, it powers up to 37", which means no more stooping or bending to lift and examine patients.

For your convenience, the 405 comes standard with a pull out

drain pan, hidden stirrups and a concealed paper roll holder behind the backrest. To help you maintain infection control, the table comes with seamless, velcro-attached upholstery that's easy to remove and clean. And it comes with front and side drawers for all your storage needs.

The Midmark 405 Power Examination Table... it's the perfect choice for every facility!

### Features

*Power up, down and back, auto return, stop switch, infinitely adjustable stirrups that lock in one of three lateral positions, foot extension shelf, hospital grade electrical receptacle, concealed paper roll holder, hold down tear straps keep paper in place, easy to clean, removable upholstery, treatment pan, four high-impact easy to clean drawers, 5° pelvic tilt*

### Specifications

Minimum height: 25" (63.5 cm)  
Maximum height: 37" (94 cm)  
Length: 55" (140 cm)  
Length with foot extension: 73" (185 cm)  
Foot extension: 20.25" W x 17.75" L (51.4 x 45.1 cm)  
Pelvic tilt: Manual, 5°  
Foot extension tilt: Manual, 7°  
Table width (overhangs base): 27" (68.6 cm)  
Headrest taper: To 15.5" (39.4 cm)



### Paper roll holder

Rear of back section for 18" x 3.5" roll (45.7 x 8.9 cm)

### Treatment pan

Polystyrene 9.5" x 11.5" x 3" (24.1 x 29.2 x 7.6 cm)

### Drawers

Top 5" deep (12.7 cm), bottom 2.25" (5.7 cm)

Drawers are 15" W x 15.25" L (37.8 x 38.7 cm)

### Accessories

Knee crutches (fixed or articulating), stainless steel drain pan, urology drain pan, restraint straps, armboard, armboard bracket, procto rest, foot control and mounting brackets for Welch Allyn® exam light and Zeiss 99 Colposcope. All accessories, except procto rest and foot control may be ordered factory installed.

### Options

Factory installed drawers on either patient's right or left side and 60 watt drawer heater in top foot end drawer. Export 50 watt drawer heater.

### Electrical

115 VAC, 50/60 Hz, 12 AMPS  
220-240 VAC, 50/60 Hz, 10 AMPS (CE Mark)

Hospital grade duplex electrical receptacle on patient's left side  
UL Listed, with or without heater  
CUL listed, without drawer heater only

### Colors

Solid upholstery colors are Monaco, Soft Pearl, Wedgewood, Monet and Royal Carriage; coordinating Designer Patterns include English Gardens, Cathedral, Waterford, Tapestry and Cabernet. Side panels and painted steel surfaces are Pearl Grey.

CE



### ISO 9001 Certified

For more information or a demonstration, contact your Midmark dealer or Midmark Corporation, Versailles, Ohio.  
Call: 1-800-MIDMARK  
Fax: 1-800-365-8631  
*Outside the U.S.A.*  
Call: 1-937-526-3662  
Fax: 1-937-526-4604  
midmark.com



Because we care.

**MIDMARK®**