

stryker[®]
INSTRUMENTS



in the POWERFUL CHOICE
CORDLESS INNOVATION

SYSTEM 4 ROTARY HANDPIECE

STRYKER INSTRUMENTS, THE LEADER IN CORDLESS POWER,
BRINGS ITS 4TH GENERATION OF BATTERY POWERED INSTRUMENTATION...

...ADVANCING ORTHOPAEDICS WITH POWER, SPEED,
VERSATILITY, RELIABILITY AND PERFORMANCE.

performance
reliability
versatility
power
speed



4103 ROTARY HANDPIECE

- Up to 66% more torque at 4:1
- Speed/Torque ratios enhanced
- Modular design allows for increased versatility and reliability




TRIGGER ASSEMBLY


- 180° rotation allows for quick and easy selection of direction/safety options

SYSTEM 4 ROTARY ATTACHMENTS

HIGH SPEED DRILLING



4103-131 1:1 1/4" Jacob's Drill




4103-110 1:1 Synthes® Drill



4103-160 1:1 Standard Trinkle Drill



4103-82-131 Keyless Chuck Adaptor



4103-180 2:1 High Torque 1/4" Jacob's Drill

HIGH TORQUE REAMING




4103-231 3.25:1 1/4" Jacob's Reamer



4103-235 3.25:1 Hudson® Modified Trinkle Reamer



4103-210 3.25:1 Synthes® Reamer



4103-435 4:1 High Torque Hudson®/Modified Trinkle Reamer

SYSTEM 4 CHARGING AND DIAGNOSTICS



• Accepts System 2000 and Cordless Driver batteries



4110-120 CHARGING STATION

- New 4 bay design optimizes efficiency and storage space
- Battery gauge allows for capacity analysis and diagnosis for System 4 batteries
- Odometer reflects battery usage allowing for warranty and replacement tracking for System 4 batteries

4115 SYSTEM 4 BATTERY

- 30% Longer run time
- Additional sterilization options
- Increased durability



SYSTEM 4 SAWS



4108 SAGITTAL SAW 4106 RECIPROCATING SAW


- Reduced Decibel Level
- 30% Longer run time



4107 RECIP/STERNUM SAW

- Smaller
- Lighter
- Enhanced blade and guide assembly

SYSTEM 4 ACCESSORIES



Complete System 4 Ordering Information

SYSTEM 4 HANDPIECES

- | | |
|------|-------------------|
| 4103 | Rotary Handpiece |
| 4106 | Reciprocating Saw |
| 4108 | Sagittal Saw |

SYSTEM 4 ATTACHMENTS

- | | |
|-------------|---|
| 4103-131 | 1:1 1/4" Jacob's Drill |
| 4103-160 | 1:1 Standard Trinkle Drill |
| 4103-133 | 1:1 Keyless Drill |
| 4103-110 | 1:1 Synthes® Drill |
| 4103-135 | 1:1 Hudson® Modified Trinkle Drill |
| 4103-180 | 2:1 High Torque 1/4" Jacob's Drill |
| 4103-235 | 3.25:1 Hudson® Modified Trinkle Reamer |
| 4103-260 | 3.25:1 Standard Trinkle Reamer |
| 4103-210 | 3.25:1 Synthes® Reamer |
| 4103-213 | 3.25:1 Hudson® Reamer |
| 4103-231 | 3.25:1 1/4" Jacob's Reamer |
| 4103-435 | 4:1 High Torque Hudson®/Modified Trinkle Reamer |
| 4103-82-131 | Keyless Chuck Adaptor |

SYSTEM 4 ACCESSORIES

- | | |
|----------|-------------------------------|
| 4115 | Battery |
| 4110-120 | Battery Charger, 120V |
| 4110-215 | Charger Module (2115 Battery) |
| 4110-412 | Charger Module (4112 Battery) |
| 4110-415 | Charger Module (4115 Battery) |
| 4110-115 | Battery Adaptor |
| 4102-451 | Sterilization Case |
| 4102-452 | Mini Sterilization Case |
| 2106-152 | Sterilization Basket |

STERNUM APPLICATIONS

- | | |
|------------|---------------------|
| 4107 | Recip/Sternum Saw |
| 4107-008 | Sternum Blade Guard |
| 298-97-101 | Sternum Blade |