

Feel the advantage.

Cordless Driver 3 Accessories

Sterilization case

- Modular sterilization containers assist in the handling, storage, and sterilization of your Cordless Driver 3 components



4200-465-000
Sterilization Case

Ordering Information

Cordless Driver 3
4300-000-000 Cordless Driver 3 Handpiece

Drilling Attachments

4100-062-000 Wire Collet (.028" - .078" diameter)
4100-110-000 AO Small Drill
4100-125-000 Pin Collet (.078" - .125" diameter)
4100-126-000 Adjustable Pin Collet
4100-131-000 1/4" Keyed Drill Chuck
4100-132-000 5/32" Keyed Drill Chuck
4100-133-000 1/4" Keyless Drill Chuck
4100-134-000 1/8" Keyless Drill Chuck
4100-135-000 Hudson®/Modified Trinkle Drill
4100-160-000 Trinkle Drill
4100-225-000 Pin Collet (.126" - .157" diameter)
4100-600-000 1:1 High Speed Bur Attachment
4100-700-000 High Speed Bur Attachment

Reaming Attachments

4100-210-000 AO Large Reamer
4100-231-000 1/4" Keyed Reamer Chuck
4100-235-000 Hudson®/Modified Trinkle Reamer
4100-260-000 Trinkle Reamer
4100-260-000s4 DHS/DCS quick lock

Other Attachments

4100-120-000 T-Latch Attachment
4100-310-000 Right Angle Drive - AO Large Reamer
4100-335-000 Right Angle Drive, Hudson®/Modified
4100-355-000 Radiolucent Right Angle Drive
4100-400-000 Sagittal Saw Attachment
4100-410-000 Radiolucent Adapter

Accessories

4100-131-132 Universal Drill Chuck Key
4200-465-000 Sterilization Case
6212-000-000 Small System 6 Battery
6215-000-000 Large System 6 Battery
6110-120-000 System 6 Battery Charger
6126-000-000 System 6 Aseptic Battery Kit (Housing, Transfer Shield, and Battery)
6126-110-000 System 6 Non-sterile Battery
6126-120-000 System 6 Aseptic Housing
6126-130-000 System 6 Transfer Shield
6110-120-000 System 6 Battery Charger
6110-625-000 System 6 Battery Module (6215, 6212, and 6126-110 Batteries)
6102-455-000 System 6 Aseptic Sterilization Case

stryker®

Joint Replacements

Trauma, Extremities & Deformities

Craniomaxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Spine

Navigation

Endoscopy

Communications

Imaging

Patient Care & Handling Equipment

EMS Equipment

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 9100-001-124 Rev. A
UnDe/PS.

Copyright © 2009 Stryker
Printed in USA

Stryker Instruments
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 800 999 3811
toll free: 800 253 3210

www.stryker.com

stryker®

Surgical



The proof is in the performance.

Cordless Driver

3

POWER

VERSATILITY

COMPATIBILITY

Cordless Driver 3


Cordless Driver 3 extends the time-honored tradition of power and versatility in the Sports Medicine arena. With the improved Cordless Driver 3, Stryker continues to offer state-of-the-art battery powered equipment that you have come to expect from Stryker Instruments.



Battery Charging System

Battery Charger

- Quickly charges Cordless Driver 3 System 6 and non-sterile batteries.
- Modular design allows for 4 different batteries to charge simultaneously.
- Interactive screen allows user to assess the condition of the battery and read battery usage
- Battery charging station diagnoses battery levels and usage
- Conditioning cycle for discharging non-sterile batteries.



Battery Options



Aseptic Battery System

- Simple 3 piece aseptic delivery to the sterile field
- 90% more capacity than the System 5 aseptic battery
- **Powerful**
Significant increase in battery capacity
- **Simple**
Easy to handle and deliver in the sterile field
- **Versatile**
Modules accommodate all sterile and non-sterile System 6 batteries



6126-000-000 System 6 Aseptic Battery Kit

Drilling Attachments



Reaming Attachments



Sagittal Saw Attachment

- Variable speed for precision cutting
- Unique tool-free blade mount reduces blade whip & flex.

Cutting Accessories

- Stryker designs, produces and distributes the industry's highest quality cutting accessories.

