


Optical Esophagoscope

with Air-inflation device and Built-in Suction Channel.



12015 AA **HOPKINS® Straight Forward Telescope 0°, autoclavable,** diameter 5.5 mm, length 53 cm. Fiber optic light transmission incorporated. Color code: green.

12016 E **Esophagoscope Tube, size 11,** O. D. 12 mm, I. D. 11 mm. With suction channel and air inflation device including glass window plug with LUER-lock 12016 EW. With fiber optic light carrier 12016 EM.

12016 F **Esophagoscope Tube, size 12.** O. D. 13 mm, I. D. 12 mm. With suction channel and air inflation device including glass window plug with LUER-lock 12016 EW. With fiber optic light carrier 12016 FM.

12016 C **A. F. (anti-fog) Tube** with instrument channel, for use with Straight Forward Telescope 12015 AA.

12015 AA

12016 C

12016 E/F