

Ordering Information

A.T.S. 2000 System

| Catalog No. | Description |
|-------------|---|
| 60-2000-101 | A.T.S. 2000 Tourniquet System with Hoses |
| 60-2000-102 | A.T.S. 2000 Tourniquet System (with A.T.S. Cylindrical Cuffs) with Hoses |
| 60-2000-103 | A.T.S. 2000 Tourniquet System (with A.T.S. Cylindrical Cuffs) with Hoses, Tourniquet Stand and Basket |
| 60-2000-104 | A.T.S. 2000 Tourniquet System (without Cuffs) with Hoses, Tourniquet Stand and Basket |

Tourniquet Accessories

| Catalog No. | Description |
|-------------|---|
| 60-4022-001 | Tourniquet Stand (Includes Base and Pole) |
| 60-1908-001 | Tourniquet Accessory Basket |

A.T.S. Cylindrical Cuffs, Dual Port, Single Bladder

| Catalog No. | Length |
|-------------|----------------|
| 60-7500-001 | 8 in. (20cm) |
| 60-7500-002 | 12 in. (30cm) |
| 60-7500-003 | 18 in. (46cm) |
| 60-7500-004 | 24 in. (61cm) |
| 60-7500-005 | 30 in. (76cm) |
| 60-7500-006 | 34 in. (86cm) |
| 60-7500-007 | 42 in. (107cm) |

A.T.S. Contour Cuff, Dual Port, Single Bladder

| Catalog No. | Length |
|-------------|----------------|
| 60-7500-008 | 42 in. (107cm) |

A.T.S. Cylindrical Cuffs, Dual Port, Dual Bladder

| Catalog No. | Length |
|-------------|---------------|
| 60-7555-001 | 12 in. (30cm) |
| 60-7555-002 | 18 in. (46cm) |
| 60-7555-003 | 24 in. (61cm) |

Zimmer Sterile Disposable Tourniquet Cuffs with Positive-locking Connectors, Dual Port, Single Bladder (Sterile packed in a box of 10 with attached sterile tubing on each cuff)

| Catalog No. | Length |
|-------------|----------------|
| 60-7070-101 | 8 in. (20cm) |
| 60-7070-102 | 12 in. (30cm) |
| 60-7070-103 | 18 in. (46cm) |
| 60-7070-104 | 24 in. (61cm) |
| 60-7070-105 | 30 in. (76cm) |
| 60-7070-106 | 34 in. (86cm) |
| 60-7070-107 | 42 in. (107cm) |

Zimmer Sterile Disposable Tourniquet Cuffs with Positive-locking Connectors, Dual Port, Dual Bladder (Sterile packed in a box of 10 with attached sterile tubing on each cuff)

| Catalog No. | Length |
|-------------|---------------|
| 60-7080-101 | 12 in. (30cm) |
| 60-7080-102 | 18 in. (46cm) |
| 60-7080-103 | 24 in. (61cm) |

Zimmer Sterile Disposable Esmarch Bandages (Sterile packed in a box of 20)

| Standard No. | Latex-Free No. | Length |
|--------------|----------------|----------------------------|
| 60-2694-101 | 60-2694-011 | 3 in. x 9 ft. (8cm x 3m) |
| 60-2694-102 | 60-2694-012 | 4 in. x 9 ft. (10cm x 3m) |
| 60-2694-103 | 60-2694-013 | 6 in. x 9 ft. (15cm x 3m) |
| 60-2694-104 | 60-2694-014 | 3 in. x 12 ft. (8cm x 4m) |
| 60-2694-105 | 60-2694-015 | 4 in. x 12 ft. (10cm x 4m) |
| 60-2694-106 | 60-2694-016 | 6 in. x 12 ft. (15cm x 4m) |



Cuffs with Sleeves

(Sterile packed in a box of 10)

| Catalog No. | Description |
|-------------|--|
| 60-7070-051 | 8 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-052 | 12 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-053 | 18 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-054 | 24 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-055 | 30 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-056 | 34 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-057 | 42 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-151 | 8 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-152 | 12 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-153 | 18 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-154 | 24 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-155 | 30 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-156 | 34 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-157 | 42 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-051 | 8 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-052 | 12 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-053 | 18 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-054 | 24 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-055 | 30 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-056 | 34 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-057 | 42 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-151 | 8 in. SP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-152 | 12 in. SP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-153 | 18 in. SP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-154 | 24 in. SP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-155 | 30 in. SP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-156 | 34 in. SP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-157 | 42 in. SP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7080-051 | 12 in. DP, DB Cuff W/PLC* and Protective Sleeve |
| 60-7080-052 | 18 in. DP, DB Cuff W/PLC* and Protective Sleeve |
| 60-7080-053 | 24 in. DP, DB Cuff W/PLC* and Protective Sleeve |
| 60-7080-151 | 12 in. DP, DB Cuff W/PLC* and Protective Sleeve W/PLC* |
| 60-7080-152 | 18 in. DP, DB Cuff W/PLC* and Protective Sleeve W/PLC* |
| 60-7080-153 | 24 in. DP, DB Cuff W/PLC* and Protective Sleeve W/PLC* |
| 60-7085-051 | 12 in. SP, DB Cuff W/PLC* and Protective Sleeve |
| 60-7085-052 | 18 in. SP, DB Cuff W/PLC* and Protective Sleeve |
| 60-7085-053 | 24 in. SP, DB Cuff W/PLC* and Protective Sleeve |
| 60-7085-151 | 12 in. SP, DB Cuff W/PLC* and Protective Sleeve W/PLC* |
| 60-7085-152 | 18 in. SP, DB Cuff W/PLC* and Protective Sleeve W/PLC* |
| 60-7085-153 | 24 in. SP, DB Cuff W/PLC* and Protective Sleeve W/PLC* |

* Positive Locking Connector

Cuffs with Sleeves Each

(Sterile packed in a box of 1)

| Catalog No. | Description |
|----------------|---|
| 60-7070-053-01 | 18 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-056-01 | 34 in. DP, SB Cuff and Protective Sleeve |
| 60-7070-153-01 | 18 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-154-01 | 24 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-155-01 | 30 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7070-156-01 | 34 in. DP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-053-01 | 18 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-056-01 | 34 in. SP, SB Cuff and Protective Sleeve |
| 60-7075-153-01 | 18 in. SP, SB Cuff and Protective Sleeve W/PLC* |
| 60-7075-156-01 | 34 in. SP, SB Cuff and Protective Sleeve W/PLC* |

* Positive Locking Connector

Sleeves Only

(Sterile packed in a box of 10)

| Catalog No. | Description |
|-------------|--|
| 60-8000-008 | Protective Sleeve for 8 in. Disposable Cuff |
| 60-8000-012 | Protective Sleeve for 12 in. Disposable Cuff |
| 60-8000-018 | Protective Sleeve for 18 in. Disposable Cuff |
| 60-8000-024 | Protective Sleeve for 24 in. Disposable Cuff |
| 60-8000-030 | Protective Sleeve for 30 in. Disposable Cuff |
| 60-8000-034 | Protective Sleeve for 34 in. Disposable Cuff |
| 60-8000-042 | Protective Sleeve for 42 in. Disposable Cuff |

Contact your Zimmer Representative or visit us at www.zimmer.com.



ATS
2000
ATS

THE ATS 2000
TOURNIQUET SYSTEM



*The Choice for
Safety, Reliability
and Convenience*



Safety.

Zimmer's dual-port system

helps protect your patients by

monitoring and controlling actual

cuff/bladder pressures.

Reliability.

Advanced microprocessor technology

includes self-diagnostic tests and

continuous monitoring for consistently

dependable performance.

Convenience.

Fast inflation and deflation save time

for busy operating room personnel.



The safety you need.
The reliability and convenience you want.

Dual controls, dual cuffs, and dual readouts provide independent and precise monitoring and performance for bilateral or Bier Block procedures.



Accurate pressure readings are taken directly from each cuff.



- Low-profile cuffs stay in place without discomfort.
- Select reusable cuffs or disposable sterile cuffs to eliminate cross-contamination.
- Cuff sleeves help reduce wrinkles, pinching, and shearing of soft tissue.

Color-coded cuffs from 8 in. to 42 in. accommodate virtually all patients and limbs.



- Positive-locking connectors supply safer, leak-free inflation.
- Angled ports keep connectors clear of the surgical site.
- Visual and adjustable audio alarms readily detect pressure variations.

Easy-to-read display adjusts easily, too.



Simple setup includes user-defined default pressure and time settings.



- Lightweight portable unit can be mounted on an I.V. pole or placed on a table.
- One-hour battery backup – just in case.

The A.T.S. 2000 Tourniquet System

Combining safety, reliability and convenience so you can focus on surgery.

The built-in safety, reliability and convenience of the A.T.S. 2000 Tourniquet frees you to focus on surgery. Dual displays that monitor actual cuff pressure, precise microprocessor controls, and one-hour battery backup help assure patient safety. User-set alarms warn of unacceptable or hazardous conditions.

The dual-cuff display is easy to read, even from a distance, and touch pad controls are simple to read and use. Fast cuff inflation and deflation save time. A large selection of cuffs, disposable and reusable, fits virtually all patients. The lightweight unit can be conveniently positioned, and positive-locking connectors have angled ports to keep lines clear of the surgical site.

