




The OLYMPUS OES Operating hysteroscope (8 mm) provides maximum capability for operative hysteroscopy procedures.

This multi-media continuous flow instrument allows the inflowing distension media to pass through a special telescope irrigation channel in order to reach the distal end of the scope directly. This guarantees better viewing conditions even when reduced flow rates are used.



A 4605A
A 4606
A 4612
A 4680
A 4617

Hysteroscope – 8 mm System




Adapters



A 4680 Standard adapter for A 4612

A 0287 Sealing cap (not pictured)

Flexible Forceps (7 Fr.)



A 4616 Flexible biopsy forceps

A 4617 Flexible scissors

A 4618 Flexible foreign body forceps