

## IDEAL EYES HD Autoclavable C-Mount Arthroscopes



- » Improved Center To Edge Resolution
- » Resolution and brightness are enhanced across the entire field of view for optimal image quality
- » Enhanced outer radius line pairs per millimeter enhanced from 11.4 horizontal lines on current scopes to 12.8 horizontal lines on IDEAL EYES scopes
- » New Magnetic Drive Focus Ring Design
- » Optimal moisture resistance
- » Improved reliability during all sterilization methods
- » Precise focusing of the C-Mount scope with one hand

Part Number	Outer Diameter	View Direction	Field of View	Working Length	Light Post Direction	Locking Mechanism
502-704-010	4.0 mm	0°	105°	140 mm	Standard	Speed-Lock
502-704-030	4.0 mm	30°	105°	140 mm	Standard	Speed-Lock
502-704-070	4.0 mm	70°	105°	140 mm	Standard	Speed-Lock
502-704-130	4.0 mm	30°	105°	140 mm	Reverse Cant	Speed-Lock
502-804-030	4.0 mm	30°	105°	165 mm	Standard	Speed-Lock
502-804-070	4.0 mm	70°	105°	165 mm	Standard	Speed-Lock
502-804-130	4.0 mm	30°	105°	165 mm	Reverse Cant	Speed-Lock
502-727-010	2.7 mm	0°	105°	120 mm	Standard	Speed-Lock
502-727-030	2.7 mm	30°	105°	120 mm	Standard	Speed-Lock
502-727-070	2.7 mm	70°	105°	120 mm	Standard	Speed-Lock
502-627-010	2.7 mm	0°	80°	75 mm	Standard	J-Lock
502-627-030	2.7 mm	30°	80°	75 mm	Standard	J-Lock
502-627-070	2.7 mm	70°	80°	75 mm	Standard	J-Lock
502-627-130	2.7 mm	30°	80°	75 mm	Reverse Cant	J-Lock
502-723-010	2.3 mm	0°	80°	72 mm	Standard	J-Lock
502-723-030	2.3 mm	30°	80°	72 mm	Standard	J-Lock
502-723-070	2.3 mm	70°	80°	72 mm	Standard	J-Lock
502-723-130	2.3 mm	30°	80°	72 mm	Reverse Cant	J-Lock
502-719-010	1.9 mm	0°	65°	58 mm	Standard	J-Lock
502-719-030	1.9 mm	30°	65°	58 mm	Standard	J-Lock