

stryker[®]
INSTRUMENTS



is the POWERFUL CHOICE
CORDLESS INNOVATION

SYSTEM 4 ROTARY HANDPIECE

STRYKER INSTRUMENTS, THE LEADER IN CORDLESS POWER,
BRINGS ITS 4TH GENERATION OF BATTERY POWERED INSTRUMENTATION...

...ADVANCING ORTHOPAEDICS WITH POWER, SPEED,
VERSATILITY, RELIABILITY AND PERFORMANCE.

performance
reliability
versatility
power
speed



4103 ROTARY HANDPIECE

- Up to 66% more torque at 4:1
- Speed/Torque ratios enhanced
- Modular design allows for increased versatility and reliability



TRIGGER ASSEMBLY

- 180° rotation allows for quick and easy selection of direction/safety options

SYSTEM 4 ROTARY ATTACHMENTS

HIGH SPEED DRILLING



4103-131 1:1 1/4" Jacob's Drill




4103-110 1:1 Synthes® Drill



4103-160 1:1 Standard Trinkle Drill



4103-82-131 Keyless Chuck Adaptor



4103-180 2:1 High Torque 1/4" Jacob's Drill


HIGH TORQUE REAMING



4103-231 3.25:1 1/4" Jacob's Reamer



4103-235 3.25:1 Hudson® Modified Trinkle Reamer




4103-210 3.25:1 Synthes® Reamer



4103-435 4:1 High Torque Hudson®/Modified Trinkle Reamer

SYSTEM 4 CHARGING AND DIAGNOSTICS



- Accepts System 2000 and Cordless Driver batteries



4110-120 CHARGING STATION

- New 4 bay design optimizes efficiency and storage space
- Battery gauge allows for capacity analysis and diagnosis for System 4 batteries
- Odometer reflects battery usage allowing for warranty and replacement tracking for System 4 batteries

4115 SYSTEM 4 BATTERY

- 30% Longer run time
- Additional sterilization options
- Increased durability



SYSTEM 4 SAWS



4108 SAGITTAL SAW 4106 RECIPROCATING SAW


- Reduced Decibel Level
- 30% Longer run time



4107 RECIP/STERNUM SAW

- Smaller
- Lighter
- Enhanced blade and guide assembly

SYSTEM 4 ACCESSORIES



Complete System 4 Ordering Information

SYSTEM 4 HANDPIECES

4103	Rotary Handpiece
4106	Reciprocating Saw
4108	Sagittal Saw

SYSTEM 4 ATTACHMENTS

4103-131	1:1 1/4" Jacob's Drill
4103-160	1:1 Standard Trinkle Drill
4103-133	1:1 Keyless Drill
4103-110	1:1 Synthes® Drill
4103-135	1:1 Hudson® Modified Trinkle Drill
4103-180	2:1 High Torque 1/4" Jacob's Drill
4103-235	3.25:1 Hudson® Modified Trinkle Reamer
4103-260	3.25:1 Standard Trinkle Reamer
4103-210	3.25:1 Synthes® Reamer
4103-213	3.25:1 Hudson® Reamer
4103-231	3.25:1 1/4" Jacob's Reamer
4103-435	4:1 High Torque Hudson®/Modified Trinkle Reamer
4103-82-131	Keyless Chuck Adaptor

SYSTEM 4 ACCESSORIES

4115	Battery
4110-120	Battery Charger, 120V
4110-215	Charger Module (2115 Battery)
4110-412	Charger Module (4112 Battery)
4110-415	Charger Module (4115 Battery)
4110-115	Battery Adaptor
4102-451	Sterilization Case
4102-452	Mini Sterilization Case
2106-152	Sterilization Basket

STERNUM APPLICATIONS

4107	Recip/Sternum Saw
4107-008	Sternum Blade Guard
298-97-101	Sternum Blade