

Camera Control Unit	Product No.	Product Description/Name
	1088-010-000	1088 HD Camera Control Unit



Camera Head & Coupler	Product No.	Product Description/Name
	1088-210-122	1088 HD Camera Head w/ 24 mm Coupler, 3.15m cable



Camera Head Only	Product No.	Product Description/Name
	1088-210-105	1088 C-Mount Head, 3.15m cable
	1088-510-105	1088 C-Mount 90deg Urology camera Head, 3.15m cable
	1088-710-105	1088 C-Mount In-Line Head, 3.15m cable
	1088-910-105	1088 C-Mount Head, 6.3m cable



Coupler	Product No.	Product Description/Name
	1088-020-122	1088 HD, 24mm Coupler



Other High Definition Equipment	Product No.	Product Description/Name
	0220-190-000	X7000 Lightsource
	0240-050-888	SDC HD
	0240-030-900	Stryker Vision 19" Flat Panel Monitor



1088 HD HIGH DEFINITION CAMERA



A Revolution in Resolution

stryker
ENDOSCOPY

Stryker Endoscopy

5900 Optical Court
San Jose, CA 95138
Tel: 1.800.624.4422
Fax: 1.800.729.2917

www.strykerendo.com

1000-900-454 Rev A.

Change the way you perform your procedures with the stunning resolution and fidelity of the

Stryker HD Camera. The Hi Definition signal from this camera provides two thirds more information along with Progressive Scan technology resulting in true-to-life color reproduction, improved motion control, and the unbelievable visual precision to enable the most meticulous of surgical movements. Let us show you why the Stryker HD is A Revolution in Resolution.



Stryker's Hi Definition Camera

creates
incredible resolution
never before seen in the O.R.



Provides 66% greater focusing precision

HD Features:

- » Hi Definition (HD) Video – 1280 x 1024 native resolution
- » Incredibly small and ergonomic camera head – for superior comfort and control
- » Standard aspect ratio – allows you to maximize existing video equipment
- » Progressive scan technology – optimum image quality with superb resolution
- » Digital Outputs – DVI necessary to take advantage of digital Hi Definition fidelity
- » Four button camera head design – control of 6 functions from the camera head
- » Multi-Specialty Settings – user selectable specialty settings which customize video output for the unique needs of your cases

The Only Integrated Platform Built Around Superior HD Resolution

- » HD Camera – Hi Definition video means uncompromised clarity.
- » X7000 Lightsource – A lightsource designed around the unique optical needs of an HD system. Single-handed chord connections improve efficiency in the OR.
- » SDC – Capture, store, retrieve and annotate your high resolution images and video for clarity that will amaze your patients and colleagues, alike.
- » StrykerVision Monitor – A 19" flat panel custom designed with the color depth and enhancement requirements dictated by minimally invasive procedures.

