

finest surgical technology

Ackermann[®]

NEW DISPOSABLE TROCAR

This new system is consisting of a threaded transparent trocar sheath and a dilating („non-cutting“) obturator. *Dilating* means that the introduction of this obturator into the abdominal cavity is not performed by a strong push but by deliberate rotating motion of approximately 180°.

The benefits of such a dilating obturator are as follows:

- less port-site bleeding
- shorter wound scars
- decreased pain
- higher patient satisfaction
- cost effective
- available in custom packs
(2 cannulas one obturator)

Try our *antiSCAR* Technology



This new disposable trocar system is available in 2 models:

Art.No. 10-1088DD	5,5mm without stopcock
Art.No. 10-1088DDH	5,5mm with LL stopcock