

# M-Series Stretchers

**Power-washable design with  
700-pound weight capacity**

**Litter-mounted oxygen bottle holder  
helps eliminate caregiver bending**

**Glideaway™ siderails create a  
zero-transfer gap**

**Optional integrated scale system  
provides accurate, repeatable weights**



The advanced SM304 eliminates manual pushing through Stryker's Zoom® motorized technology.

## **Mobility, Simplicity, Versatility**

Fully-equipped and completely adaptable, the M-Series allows you to choose the mobility solution that best suits your needs. Supremely mobile and user-friendly, these stretchers greatly reduce the physical strain on clinicians without sacrificing patient care and comfort.




The SM104 uses a retractable fifth wheel to provide superior maneuverability and cornering.



The SM204's Big Wheel® reduces start-up force by 50% and steering effort by 60%.

## Accent Color Options



# M-Series Stretchers

**stryker**<sup>®</sup>



SM104 with Fifth Wheel



SM204 with Big Wheel



SM304 with Zoom Drive

## Medical

3800 E. Centre Ave.  
Portage, MI 49002  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

### Standard Features

- Thermoformed ABS base with extra capacity storage tray
- 700 lbs weight capacity
- 26" (66 cm) wide patient surface with 3" (8 cm) Enhanced Comfort Mattress
- Less than 21" low height
- Integrated utility tray
- Litter-mounted oxygen bottle holder
- Integrated transfer board
- Power-washable design (SM104 & SM204 only)
- Four-wheel, steel-ring brake system with dual-end activators
- Dual side-mounted foot control hydraulics with uni-lower pedal
- Glideaway™ siderails with integrated push/pull handles
- Four integrated IV receptacles
- Trend./reverse Trend.
- 8" (20 cm) Omni Surface™ casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Constant descent hydraulics
- Dual-end brake & steer controls

#### SM104: All Standard Features Plus CE

- Fifth Wheel steering

#### SM204: All Standard Features Plus CE

- Big Wheel steering

#### SM304: All Standard Features Plus C-US UL

- Zoom Motorized Drive System

### Specifications

	26" Litter		30" Litter	
• <b>Model Number</b>	<b>1007 (SM104)</b>		<b>1007 (SM104)</b>	
	<b>1015 (SM204)</b>		<b>1015 (SM204)</b>	
	<b>1025 (SM304)</b>		<b>1025 (SM304)</b>	
• <b>Overall Length</b>				
SM104 & SM204	84"	(213 cm)	84"	(213 cm)
SM304	85"	(216 cm)	85"	(216 cm)
• <b>Overall Width</b>				
Siderails Up	33.5"	(85 cm)	37"	(94 cm)
Siderails Down	30.25"	(77 cm)	30.75"	(78 cm)
• <b>Weight Capacity</b>	700 lbs	(317 kg)	700 lbs	(317 kg)
• <b>Height Range</b>				
High	34.5"	(88 cm)	34.5"	(88 cm)
Low (SM104 & SM204)	20.75"	(53 cm)	20.75"	(53 cm)
Low (SM304)	23"	(58 cm)	23"	(58 cm)
• <b>Litter Positioning</b>				
Backrest	0° - 90°		0° - 90°	
Knee Gatch	40°		40°	
Trend./reverse Trend. (SM104 & SM204)	±16°		±16°	
Trend./reverse Trend. (SM304)	±15°		±15°	
Maximum Angle	87.5°		87.5°	
• <b>Patient Surface</b>	26" x 75.25"	(66 x 191 cm)	30" x 75.25"	(76 x 191 cm)
• <b>Siderails</b>	14" x 57.75"	(36 x 147 cm)	14" x 57.75"	(36 x 147 cm)
• <b>Caster Diameter</b>	8"	(20 cm)	8"	(20 cm)

Stryker reserves the right to change specifications without notice.

### Warranty

- Two years parts, labor and travel

### Optional Features

- 30" (76 cm) wide patient surface
- 4" (10 cm) Enhanced Comfort Mattress
- 4" (10 cm) or 5" (13 cm) Ultra Comfort Mattress
- Four-sided brake & steer controls (SM104 & SM204 only)
- Pop-up push handles (SM104 & SM204 only)
- No-crank hydraulic knee gatch
- Exclusive backrest lift-assist design with drop seat
- Integrated scale system
- Three-sided hydraulic controls (SM104 & SM204 only)
- Defibrillator tray
- Defibrillator tray/foot extender/chart surface
- Dual-end siderail release
- Footboard/chartholder
- IV caddy
- IV poles
  - Permanent, folding two-or three-stage
- Serving tray
- Serving tray holder/footboard
- Siderail pads
- Upright oxygen bottle holder
- Restraint straps
- Retractable cord reel (SM304 only)