

Expanding Your Cordless Capabilities

**stryker**<sup>®</sup>  
**INSTRUMENTS**



***Cordless Driver 2***


more speed, more torque  
**Driven to Perfection...**

A History of Excellence...

...Driven to Perfection

**Cordless  
Driver 2**

Increased speed and torque are combined with proven cordless reliability to provide unsurpassed power and versatility.



**More Speed...  
More Torque...**

Proven power and reliability  
for applications such as:

- Sports Medicine
- Trauma
- Total Joint

**50% More Speed**


**20% More Torque**



## Drilling Attachments

Accommodates a variety of drilling demands.

Automatically adjustable wire & pin collets simplify the placement of wires and guide pins.



4100-133 1/4" Keyless Drill



4100-135 Hudson®/Modified Trinkle Drill



4100-134 1/8" Keyless Drill



4100-160 Trinkle Drill



4100-131 1/4" Jacobs Drill



4100-132 5/32" Jacobs Drill



4100-110 Synthes® Drill



Easy to use Keyless Chucks in 1/4" and 1/8" sizes.

Cordless Driver 2 attachments are color coded for easy selection and to reference quickly those used for drilling and reaming.

## Reaming Attachments

Power to tunnel with trephines and large diameter drill bits.

Reamer collets deliver high torque, excellent for reaming Intermedullary Canal and ACL reconstruction applications.



4100-235 Hudson Modified Trinkle Reamer



4100-260 Trinkle Reamer



4100-231 Jacobs Reamer



4100-210 Synthes Reamer



4100-335 Right Angle Drive Modified Trinkle Reamer



4100-310 Right Angle Drive Synthes Reamer



4100-355 Radiolucent Right Angle Drive

## Sagittal Attachment

Variable speed for precision cutting.


## Sagittal Blades

Stryker Precision Blades are available in various sizes and thicknesses to address the specific requirements of each procedure.


See bur & blade guides for a complete list of available cutting accessories.




Unique tool-free blade mount reduces blade whip & flex.




MEDIUM,  
NARROW,  
SAFEDGE™




AGGRESSIVE




ACL,  
10MM STOP




NARROW



WIDE



EXTRA SHORT  
NARROW



MEDIUM



## Sterilization & Storage/Battery Charging

Modular sterilization container assists in the handling, storage and sterilization of your Cordless Driver 2 components.

Battery charging station diagnoses battery levels and usage.

### 4200-465 Sterilization Case

- Holds handpiece and all accessories for complete sterilization.

- Convenient inner tray presents instruments on a sterile field.

- Modular design allows for 4 batteries to be charged at the same time.

- Indicator lights provide battery status feedback with "Charge", "Ready", and "Replace" indicators.

### 4110-120 Battery Charger

- Quickly charges Cordless Driver 2 and System 4 batteries.





Expanding Your Cordless Capabilities

**CD2 Handpiece**

PART #	DESCRIPTION
4200	Cordless Driver 2 handpiece (includes the wire collet 4100-62)

**Drilling Attachments**



PART #	DESCRIPTION
4100-62	Wire Collet (.28"-.78" diameter)
4100-110	Synthes® Drill
4100-125	Pin Collet (.78" - .125" diameter)
4100-126	Adjustable Pin Collet
4100-131	1/4" Jacobs® Drill
4100-132	5/32" Jacobs® Drill
4100-133	1/4" Keyless Drill Chuck
4100-134	1/8" Keyless Drill Chuck
4100-135	Hudson®/Modified Trinkle Drill
4100-160	Trinkle Drill

**Reaming Attachments**



PART #	DESCRIPTION
4100-210	Synthes® Reamer
4100-231	1/4" Jacobs® Reamer
4100-235	Hudson®/Modified Trinkle Reamer
4100-260	Trinkle Reamer



**Other Attachments**

PART #	DESCRIPTION
4100-310	Right Angle Drive - Synthes®
4100-335	Right Angle Drive - Hudson®/Modified
4100-355	Radiolucent Right Angle Drive
4200-355-32	Radiolucent Drive Drill Bit (3.2mm)
4200-355-35	Radiolucent Drive Drill Bit (3.5mm)
4200-355-40	Radiolucent Drive Drill Bit (4.0mm)
4200-355-45	Radiolucent Drive Drill Bit (4.5mm)
4100-400	Sagittal Saw Attachment (for use with the 2296 Micro Saw Blades)

**Cordless Driver 2 Accessories**

PART #	DESCRIPTION
4112	Battery
4110-120	Battery Charger, 120 VOLT
4110-112	Battery Adapter for use with 4112 Batteries and System 2000 Charger
4110-412	System 4 Charger Module for 4112 Battery
4100-131-132	Universal Jacobs® Key
4100-451	Sterilization Case and Tray
4200-465	Plastic Sterilization Case

See Bur and Blade literature for complete listing of available cutting accessories



For more information, contact your local Stryker Sales Representative.

4100 East Milham Ave  
Kalamazoo, MI 49001  
(800)253-3210  
(616)323-7700