

Power
Examination
Table

411



MIDMARK®

**Power
Examination
Table**

411

The Midmark 411 is the ultimate universal power table because it offers everything you need for your medical environment. With its ergonomic design, the 411 not only looks great... but it comes standard with both a programmable foot and hand control for added convenience. Program the table to your four favorite positions, and it will move quickly and precisely where you want it... every time.

A low height of 22½" makes the 411 more accessible to infirm, elderly and pregnant patients, as well as to those patients in wheelchairs. And a full height of 40" means no more stooping or bending to lift patients, so you and your staff will feel better at the end of each day.

Smart... efficient... dependable... and of course, the 411 coordinates with our other Midmark products to give your facility a look that's complete.

Features

Power up and down, tilt, back and foot, programmable hand and foot controls, infinitely adjustable stirrups that lock in one of three lateral positions, hospital grade electrical receptacles, paper roll holder hidden under the double arm headrest, auto return, stop switch, hold-down straps keep paper in place, accessory rails and stainless steel treatment pan



Specifications

- Height: 22.5" to 40" (57.15 to 101.6 cm)
- Table length: 70.5" (179.1 cm) extends to 81" (205.7 cm)
- Table width: 27" (68.6 cm)
- Trendelenburg: 30°
- Procto knee rest: Manual height adjustment of 11", 13", 15" and 17" to top of seat (28, 33.1, 38.1 and 43.2 cm)
- Stirrup angle: Straight leg, 15° and 30°
- Standard headrest: 11" H x 14" W at base tapering to 11" at top (27.9 x 35.6/27.9 cm)

Treatment pan

- Stainless steel 9" x 11.75" x 2.5" (22.9 x 29.8 x 6.4 cm)

Hand and Foot Controls

Hand and foot controls are standard

- 8-way power control buttons and four programmable positions on foot control and hand control
- 'Auto Return' button lowers table to low height position
- 'Stop' button overrides automatic return
- Easy to interpret diagrams
- Plug into either front or back of base pedestal
- Foot control helps to maintain sterile field
- Hand control stores conveniently on bracket

Accessories

Knee crutches (standard or articulating), chair arms, caster base, facial pad (fits standard and special procedures headrest), restraint straps, base rails, seat rails, IV pole w/clamp, vision block screen w/clamp, instrument tray w/clamp, armboard w/pad and clamp (fixed or articulating), hand procedure armboard w/pad, leg and clamp, special procedures headrest, wrist support headrest, side rails and Welch Allyn® light bracket

Options

- Traditional upholstery (*shown on cover*)
- Seamless upholstery (*shown on back*)

Electrical

115 VAC, 60 Hz, 12 AMPS
Two duplex hospital grade electrical receptacles, one on each side of pedestal base
UL Listed
CUL Listed

Colors

Solid upholstery colors are Monaco, Soft Pearl, Wedgewood, Monet and Royal Carriage; coordinating Designer Patterns include English Gardens, Cathedral, Waterford, Tapestry and Cabernet. Side panels and painted steel surfaces are Pearl Grey.



ISO 9001 Certified

For more information or a demonstration, contact your Midmark dealer or Midmark Corporation, Versailles, Ohio.
Call: 1-800-MIDMARK
Fax: 1-800-365-8631
midmark.com



Procto



OB/GYN



Flat table

Because we care.

MIDMARK®