

# CALCUTRIPT

Electrohydraulic System  
for treatment of stones.

**Special features:**

- Easy to operate.
- High safety for patient.
- High variability range of output energy.
- Test function key to check unit before use.
- Modular design with functional monitoring as servicing aid.

**Specifications:**

**Frequency selection:**  
4 steps each for intensity and pulse frequency.  
**Pulse voltage:** max. 1900 volt.  
**Max. discharge energy:** 1.5 ws.  
**Dimensions:** 520 mm × 160 mm × 380 mm (w × h × d).  
**Weight:** 16 kg.



27080 B **CALCUTRIPT Electrohydraulic System**,  
power supply: 230 VAC, 50/60 Hz,  
**including:**  
400 A **Mains Cord.**  
27085 F **Footswitch.**  
27080 KA **Probe Connection Cable.**  
27080 TA **Testprobe, 7 Fr.**  
27080 XA **Protection Sleeve** for probes  
(3 each). 3 Probes to be  
selected from the list below.

27080 C **Same**,  
power supply: 120 VAC, 50/60 Hz.  
27080 CJ **Same**,  
power supply: 110 VAC, 50/60 Hz.  
27080 D **Same**,  
power supply: 240 VAC, 50/60 Hz.

**Accessories:**

Probes	Diameter	Length	for use with
27080 FA	7 Fr.	50 cm	Cystoscope-Urethroscope, Stone Punch, Percutaneous Nephroscopes, Rigid Cholecystotomy/Cholecystectomy Set
27080 EA	4.5 Fr.	80 cm	Cystoscope-Urethroscope, Flexible Cysto-Urethro-Fiberscope 11272 AA, Rigid Uretero-Renoscope, Nephro-Fiberscope 11001 DD, Pediatric Percutaneous Nephroscope
27080 GA	3.0 Fr.	120 cm	Flexible Uretero-Renoscope, Ureterscope, Rigid Uretero-Renoscope, Flexible Cholangio-Fiberscope 11290 A/11292 B
27080 GL	3.0 Fr.	220 cm	Gastro-Intestinal Fiberscopes (Common duct stones)
27080 HA	1.6 Fr.	120 cm	Flexible Uretero-Renoscope, Rigid Uretero-Renoscope, Cholangio-Fiberscope 11290 A/11292 B

Package of 10.