


# Optical Forceps

for use with **HOPKINS**® Telescopes 10320 AA/BA.



- The jaws of the forceps and site of biopsy are well visible before and during the biopsy procedure.
- Foreign bodies can now be removed under precise visual control.

● With lock and A. F. (anti-fog) device for telescope.  
For use with **HOPKINS**® telescopes 10320 AA/BA and bronchoscope tubes sizes 8.5-7.5:

- 10350 K **Optical Forceps** for peanuts and soft foreign bodies.
- 10350 L **Optical Biopsy Forceps**, circular cup 4.3 mm.
- 10350 H **Optical Alligator Grasping Forceps**, large jaws.
- 10350 SB **Optical Scissors**.
- 10350 U **Optical Forceps, Universal**, for biopsy, removal

of foreign bodies and denaturated tissues.

For use with **HOPKINS**® telescope 10320 AA and bronchoscope tubes sizes 8.5-6.5 mm.

- 10350 LE **Optical Biopsy Forceps** with lock and A. F. (anti-fog) device for telescope, circular cup 2 mm.

The small jaws of the forceps 10350 LE are close to the telescope for easy access to narrow areas.

