

# Olympus OES Pro Telescopes

- High optical quality
- Outstanding durability
- Five-year Goldtip™ autoclavability guarantee
- Compatible with previous Olympus cystoscopes and resectoscopes

## OES PRO TELESCOPES

Telescopes are delivered with instrument tray WA05990A.

A22000A	Telescope, 4 mm, 0° direction of view, autoclavable
A22001A	Telescope, 4 mm, 12° direction of view, autoclavable
A22002A	Telescope, 4 mm, 30° direction of view, autoclavable
A22003A	Telescope, 4 mm, 70° direction of view, autoclavable
A22004A	Telescope, 4 mm, 110° direction of view, autoclavable
A22005A	Telescope, 4 mm, 30° direction of view, wide angle, autoclavable

## OES PRO LIGHT-GUIDE CABLE


WA03200A	Light-guide cable, size S, plug type, 3 m
A0460	Light-guide adapter for OES 4000 light-guide cables
O0332	Light-guide adapter for Wolf light-guide cables

## INSTRUMENT TRAY

WA05990A	Instrument tray for 2 telescopes. Part of delivery of A22000A-A22005A.
A0547	Protection tube for telescopes

## BRIDGE AND WORKING INSERTS

A22093A	Bridge for inspection of the bladder. Can be used with a 70° telescope A22003A in all resection sheaths
A22091A	Working insert, 6 Fr channel, for inner sheaths
A42091A	Working insert, 8 Fr channel for outer sheaths
A0298	Sealing cap, 0.1 mm, 6/pkg
O0131	Sealing cap, 0.7 mm, 10/pkg
A0282	Sealing cap, kink protection, 7 Fr, 3-ply lamella, 10/pkg



A22001A