

Ohmeda

Suction Products

Intermittent Suction Unit (ISU)



Ohmeda's Intermittent Suction Unit is widely regarded as the industry standard for quality and reliability. It is designed with an intermittent mode for gastric drainage as well as a continuous mode for tracheal, pharyngeal, or surgical suctioning.

- Dual modes — intermittent and continuous.
- Quiet, reliable operation.
- Dependable, time-tested gauge and controls.
- Provides FULL VAC regulation.
- Low maintenance, easy to service.
- Accepts Posi-Pulse® Drainage Assist Device, using adaptable Retrofit Kit.
- Glow-in-the-dark needle and increments provide optimum visibility.


Posi-Pulse Intermittent Suction Unit



Nasogastric drainage procedures are enhanced through the exclusive Posi-Pulse device and the Posi-Pulse ISU (see Posi-Pulse device description).

- Exclusive design for use with Posi-Pulse® Drainage Assist Device.
- Dual modes — intermittent and continuous.
- Quiet, reliable operation.
- Dependable, time-tested gauge and controls.
- Provides FULL VAC regulation.
- Low maintenance, easy to service.
- Glow-in-the-dark needle and increments provide optimum visibility.

Posi-Pulse® Drainage Assist Device



An innovation in suction therapy, Posi-Pulse offers convenience and increased patient safety in nasogastric drainage procedures. Problems of irrigation, gastric backflow, cross contamination and mucosal trauma are minimized.

- For use with Posi-Pulse ISU.
- Connects to ISU, using adaptable Retrofit Kit.
- Irrigates nasogastric tubes automatically.
- Eliminates vent lumen backflow.
- Reduces mucosal trauma.
- Easy to use.
- Disposable, single patient use.

6700-0000-900

Case of 25

6700-0110-900

Case of 10

6700-0101-800

Retrofit Kit