



Get a bigger view. A better view. A sharper view. Only with the new Welch Allyn MacroView™ Otoscope.



Get nearly a full view of the tympanic membrane

The new Welch Allyn MacroView Otoscope will forever change the way you view the ear — without changing the way you work. This next-generation otoscope features advanced optics that increase magnification while giving you nearly a full view of the tympanic membrane — virtually eliminating the need for panning. What's more, with its easy-to-focus eyepiece, the MacroView allows for enhanced visualisation to help you to see past wax, locate fluid behind the eardrum, and even perform exams without your glasses. The bottom line: MacroView gives you a bigger, better, sharper view of the ear — allowing you to perform everyday exams faster and more efficiently.

The Welch Allyn MacroView Otoscope gives you:

- Twice the field of view
- 30% greater magnification than a traditional otoscope
- Nearly a full view of the tympanic membrane
- Better clarity and definition of landmarks
- The ability to adjust focus for variable ear canal lengths or farsighted eyes
- Tip grip for secure fastening of ear specula, along with ejection feature for disposal

WelchAllyn®

Advancing Frontline Care™



The Welch Allyn MacroView™ Otoscope

Ordering Information

23820 Welch Allyn MacroView™ Otoscope with throat illuminator

23810 Welch Allyn MacroView™ Otoscope

The Welch Allyn Otoscope is also available in diagnostic sets, desk charger sets, and integrated wall systems.

Accessories

23804 Insufflator Bulb

06500-U Replacement Lamp

52432-U Universal KleenSpec® Paediatric Disposable Ear Specula

52434-U Universal KleenSpec® Adult Disposable Ear Specula

24400-U Set of Four Reusable Ear Specula (2, 3, 4, & 5 mm)



NEW! SofSeal™—fits on the end of Welch Allyn ear specula for improved seal and comfort during pneumatic otoscopy

24330 SofSeal, Medium—recommended for 3, 4, & 5 mm ear specula

24320 SofSeal, Small recommended for 2, 3, & 4 mm ear specula



NEW! Slotted Specula allows instrumentation under visualisation

52700 Instrumentation Tip (compatible with 23810 and 23820 MacroView Otoscope)

Find out how the new Welch Allyn MacroView Otoscope can help take your ENT exams to the next level. For more information go to www.welchallyn.com/macroview.

Welch Allyn European Customer Services, Navan Business Park, Dublin Road, Navan, Co. Meath, Republic of Ireland

Tel: +353 46 90 67790 Fax: +353 46 90 67754

Welch Allyn UK Ltd, Cublington Road, Aston Abbots, Buckinghamshire HP22 4ND, England

Tel: +44 (0) 207 365 6780 Fax: +44 (0) 207 365 9694

Welch Allyn SA, 394 Surrey Avenue, Ferndale 2194, P O Box 809, Fontainebleau 2032, South Africa

Tel: +27 11 777 7555 Fax: +27 11 777 7556

WelchAllyn®

Advancing Frontline Care™