

RICHARD  
WOLF




**UROLOGY  
LASER  
INSERT**

# UROLOGY LASER INSERT


8650.67

The Richard Wolf Laser Insert\* can be used cystoscopically in the bladder, or percutaneously in the renal pelvis, with a flexible fiber laser. The Laser Insert deflects the laser fiber ( $\leq 1.05\text{mm O.D.}$ )  $0^\circ\text{-}60^\circ$ , and couples to a laser connector with an M10x0.75 thread size.



## FOR LASER CYSTOSCOPY

Hofstetter Laser Insert . . . . .	8650.67
Lumina SL-Telescope (25°) . . . . .	8654.33
Lumina SL-Telescope (70°) . . . . .	8650.35
19.5 Fr. Cysto-Urethroscope Sheath	8650.031



## FOR LASER NEPHROSCOPY

Hofstetter Laser Insert . . . . .	8650.67	-
Lumina SL-Telescope (25°) . . . . .	8654.33	
Lumina SL-Telescope (70°) . . . . .	8650.35	-
24.5 Fr. Perc.Nephroscope Sheath .	8962.01	
Perc.Nephroscope Inner Sheath . .	8962.11	

\*developed by Dr. Hofstetter

For more information please contact:

**RICHARD WOLF MEDICAL INSTRUMENTS CORP.**

