

## Clinaseal Safety Centrifuge



**THE IDEAL PHYSICIAN'S OFFICE CENTRIFUGE!** Our quiet Clinaseal features exclusive Aero-Lok technology that helps protect personnel from aerosols while maintaining sample temperature stability. The motor chamber is isolated from the spinning chamber and the spinning chamber is in turn isolated from the outside environment. A brushless, permanently lubricated motor means no maintenance worries. The sturdy aluminum cabinet is covered with a durable, baked-on powder coating in an attractive tan color that compliments any environment. Private labeling is available.

- ⊕ **3150 RPM, 1132 x g R.C.F.\***
- ⊕ **Cool Running Brushless Motor**
- ⊕ **Accommodates 2ml-15ml specimen tubes**
- ⊕ **2 Year Warranty**
- ⊕ **9.5 in. wide x 10.0 in. deep x 9.0 in. high**
- ⊕ **Lid Safety Switch**
- ⊕ **Sealed design offers superior protection from biohazards**

### ORDERING INFORMATION

Model CS6C 6 Tube Capacity, 30 Min. Timer

Model CS8C 8 Tube Capacity, 30 Min. Timer

Model CS6X 6 Tube Capacity, On/Off Switch



## Variseal Variable Speed Centrifuge

**THE LOWEST PRICED VARIABLE SPEED CENTRIFUGE ON THE MARKET!** The Variseal combines all the features of our best-selling Clinaseal centrifuge with the added advantage of variable speed! The variable speed control conveniently indicates both RPM and R.C.F.\*. The brushless, permanently lubricated motor means no maintenance worries. This four-in-one instrument processes blood, urine, sperm, and fecal specimens, giving you unbeatable versatility in one compact, very affordable unit.

### ORDERING INFORMATION

Model VS6C 6 Tube Capacity, 30 Min. Timer

Model VS8C 8 Tube Capacity, 30 Min. Timer

- ⊕ **1000 - 3150 RPM, 114 - 1132 x g R.C.F.\***
- ⊕ **Cool Running Brushless Motor**
- ⊕ **Accommodates 2ml-15ml specimen tubes**
- ⊕ **2 Year Warranty**
- ⊕ **9.5 in. wide x 10.0 in. deep x 9.0 in. high**
- ⊕ **Lid Safety Switch**
- ⊕ **Sealed design offers superior protection from biohazards**

## Clinaspin Economy Centrifuge

**ROCK SOLID AND AFFORDABLE!** The Clinaspin gives you everything you want in a basic centrifuge. It separates blood samples with no frills. A durable polymer cover in an attractive tan color is mounted on a black, powder coated steel base. The Clinaspin is the ultimate in reliability. A brushless, permanently lubricated motor assures long life and sample temperature integrity.



- ⊕ **3200 RPM, 1170 x g R.C.F.\***
- ⊕ **Cool Running Brushless Motor**
- ⊕ **Accommodates 2ml-15ml specimen tubes**
- ⊕ **Locking Lid**
- ⊕ **2 Year Warranty**
- ⊕ **10 in. diameter x 8 in. high**

### ORDERING INFORMATION

Model 6C 6 Tube Capacity, 30 Min. Timer

Model 8C 8 Tube Capacity, 30 Min. Timer

Model 6X 6 Tube Capacity, On/Off Switch

\*Relative Centrifugal Force