

MIDAS REX PNEUMATIC TOOLS, INC.
 3001 RACE STREET FORT WORTH TEXAS 76111-4117
 800-433-7080 817-831-4181 FAX-817-838-2384

MIDAS REX INSTITUTE, INC.
 3000 RACE STREET FORT WORTH TEXAS 76111-4116
 800-433-7639 817-831-2604 FAX-817-834-4835

LAMINATED INSTRUMENTATION TEMPLATE



(This laminated sheet may be gas sterilized)

GENERAL INFORMATION:

ILLUSTRATIONS on this template are actual size.


ATTACHMENTS are designated by a letter name.

DISSECTING TOOLS are designated by attachment *letter* and dissecting tip design *number*. Examples: A-5, TU-14, AF-2, etc.





• Available 3/32" Stem Dissecting Tools in Ascending Order of Attachment Length


		ATTACHMENTS				
		X	B	K	A	L
TIP DESIGNS AVAILABLE	1	•	•	•	•	•
	2	•	•	•	•	•
	3	•	•	•	•	
	3D		•	•	•	
	4	•	•	•	•	
	5	•	•	•	•	
	6	•	•		•	
	8D			•	•	
	9			•	•	
	10	•			•	
	13			•	•	•



Available only in standard (100) diamond grit.




Coarse (60) diamond grit
Standard (100) diamond grit



Available only in coarse (60) diamond grit.

• Available 1/8" Stem Dissecting Tools in Ascending Order of Attachment Length




TIP DESIGNS AVAILABLE	ATTACHMENT						
	O	M	AM	TU	R	RX	
1	•	•	•	•	•	•	•
2	•	•	•	•	•	•	•
3	•	•	•	•	•	•	•
4	•	•	•	•	•	•	•
5	•	•	•	•	•	•	•
6	•	•	•	•	•	•	•
7	•	•	•	•	•	•	•
8	•	•	•	•	•	•	•
8D							
8DC							
8H							
9							
10	•	•	•	•	•	•	•
11							
12	•	•	•	•	•	•	•
13							
14							
15							
16							
17							
17H							
18							
18H							
19							
19H							
23							
24							
25							
31	•	•	•	•	•	•	•
31D							
32	•	•	•	•	•	•	•
32D							
33	•	•	•	•	•	•	•
33D							
35	•	•	•	•	•	•	•
35D							
51							
52							





Metal Dissection - All WH dissecting tools and M-51, M-52, AM-51, and TU-51.

Bioceramic dissection - WH-1, M-51, AM-51, and TU-51

Extended collet (Available for AM, TU and R lengths)



The above family of AF dissection tools are for use in the AF, BF, and MF attachments.



Dissecting Tool Guides available for AC, Q, and A-L Series.

