



# Steerable Optical Biopsy Instruments

For use with Universal-Bronchoscopes 10318 B/C size 8.5 and 7.5 and **HOPKINS**® Telescopes 10320 AA/BA/DA.



- 10357 E **Deflecting Device.**
- 10358 A **Flexible Biopsy Forceps.**
- 10358 B **Same, smaller size, for transbronchial lung biopsy.**
- 10359 **Flexible Brush for cytology.**



10358 A 10357 E      10358 B 10357 E      10468 B 10357 E

The sharp cutting jaws can be pushed mechanically against the intact mucosa membrane. During this process specimen can even be obtained from segments where the usual type of forceps would slip off.