

MICROAIRE



Pneumatic Large Power Sterilization Case
9700-000

Features:

- Polymer & Aluminum construction provides durability
- Base holds 3 Handpieces & Hose
- Insert Tray is molded with cavities for 11 couplers & Pin Driver
- 10" x 20" x 4" (25cm x 50cm x 10cm)

Air Hoses & Accessories

MicroAire offers a wide variety of Air Hoses and other Accessories for use with your Series 7000 and other Pneumatic Instrument Systems

For more information, please see your MicroAire Catalog, contact your MicroAire Sales Representative, or visit us online at www.microaire.com

INSTRUMENTS & ACCESSORIES

CATALOG NUMBER	DESCRIPTION
7100-200	Drill/Reamer with Hall® Connector
7200-200	Oscillating Saw with Hall® Connector
7400-200	Reciprocating Saw with Hall® Connector
	Handpieces are also available with the following air hose connections: Linvatec/3M® "Maxi", A.O. Synthes®, Stryker®.
7100-001	Zimmer Reamer Coupler, 250 RPM
7100-002	Trinkle Coupler, 750 RPM
7100-003	Jacobs Coupler, 750 RPM
7100-004	Keyless Coupler 1/4", 750 RPM
7100-005	Keyless Coupler 1/8", 750 RPM
7100-006	A.O. Synthes Coupler, 750 RPM
7100-008	Hudson Reamer Coupler, 250 RPM
7100-010	Jacobs Coupler, 250 RPM
7100-011	Harris/Aesculap Reamer, 250 RPM
7100-111	A.O. Synthes Reamer Coupler, 250 RPM
7100-112	A.O. Synthes Non-Reversing Reamer Coupler, 250 RPM
7100-045	Wire/Pin Driver Coupler, 1.6mm- 2.8mm
7100-050	Pin Driver Coupler, 3.2mm- 4.0mm
7100-051	Zimmer Coupler, 750 RPM
7100-058	Hudson Coupler, 750 RPM
7100-474	Pulse Lavage Coupler
4100-021	Automatic Screwdriver Bit Holder for 4100-024/025/026
4100-023	Automatic Screwdriver Bit Holder for 4100-027
4100-024	Screwdriver Bit, Single Slot
4100-025	Screwdriver Bit, Phillips
4100-026	Screwdriver Bit, Cruciate
4100-027	Screwdriver Bit A.S.I.F. Trinkle Shank
7013-000	Large Power "Quiet" Air Hose with Hall® Connector - 10 Foot Length
9013-000	Standard Air Hose with Hall® Connector - 10 Foot Length
9700-000	Sterilization Case for Series 7000 Pneumatic System



SERIES 7000
PNEUMATIC LARGE POWER INSTRUMENT SYSTEM



1641 Edlich Drive Charlottesville, VA 22911
TOLL FREE (800) 722-0822 (434) 975-8000
FAX (800) 648-4309 (434) 975-4131

Visit us online at:
www.microaire.com



A member of The Marmon Group of companies

MicroAire
SURGICAL INSTRUMENTS



7100 DRILL/REAMER

DRILL
 SPEED 750 RPM
 TORQUE 50 in. lbs
REAMER
 SPEED 250 RPM
 TORQUE 150 in. lbs



7200 OSCILLATING SAW

SAW
 SPEED 15,000 CPM
 ARC 5°
 12 POSITION
 ROTATING HEAD



7400 RECIPROCATING SAW

SAW
 SPEED 15,000 CPM
 360°
 BLADE ROTATION

DRILLING



Jacobs Coupler 750 RPM 7100-003
 NAME SPEED CATALOG #



Zimmer Coupler 750 RPM 7100-051
 NAME SPEED CATALOG #



Keyless Coupler 1/4" 750 RPM 7100-004
 NAME SPEED CATALOG #



Hudson Coupler 750 RPM 7100-058
 NAME SPEED CATALOG #



A.O. Synthes Coupler 750 RPM 7100-006
 NAME SPEED CATALOG #



Keyless Coupler 1/8" 750 RPM 7100-005
 NAME SPEED CATALOG #



Trinkle Coupler 750 RPM 7100-002
 NAME SPEED CATALOG #

REAMING



Harris/Aesculap Reamer Coupler 250 RPM 7100-011
 NAME SPEED CATALOG #



Synthes Non-Reversing Reamer Coupler 250 RPM 7100-112
 NAME SPEED CATALOG #



Zimmer Reamer Coupler 250 RPM 7100-001
 NAME SPEED CATALOG #



Hudson Reamer Coupler 250 RPM 7100-008
 NAME SPEED CATALOG #



Synthes Reamer Coupler 250 RPM 7100-111
 NAME SPEED CATALOG #



Jacobs Coupler 250 RPM 7100-010
 NAME SPEED CATALOG #

PIN DRIVING



Jacobs Coupler 750 RPM 7100-003
 NAME SPEED CATALOG #



Wire/Pin Driver Coupler 750 RPM 7100-045
 NAME SPEED CATALOG #
 NOTE: ACCEPTS WIRES AND PINS FROM 1.6 MM TO 2.8MM



Pin Driver Coupler 750 RPM 7100-050
 NAME SPEED CATALOG #
 NOTE: ACCEPTS PINS FROM 3.2MM TO 4.0MM

K-Wires & Steinmann Pins SIZE FROM 1.6MM TO 4.0MM

**SCREWDRIVING COUPLERS
 BITS & HOLDERS**



Trinkle Coupler 750 RPM 7100-002
 NAME SPEED CATALOG #



Screwdriver Bit Holder 4100-021
 NAME CATALOG #



Screwdriver Bit Holder 4100-023
 NAME CATALOG #

Screwdriver Bit- Single Slot

Trinkle Shank 4100-024
 NAME CATALOG #

Screwdriver Bit- Cruciate

Trinkle Shank 4100-026
 NAME CATALOG #

Screwdriver Bit-Phillips

Trinkle Shank 4100-025
 NAME CATALOG #

Screwdriver Bit-A.S.I.F.

Trinkle Shank 4100-027
 NAME CATALOG #

NOTE: The 4100-023 bit holder fits only the 4100-027 A.S.I.F. screwdriver bit.

PULSE LAVAGE



Pulse Lavage Coupler 7100-474
 NAME CATALOG #

NOTE: Uses 4740 tubing sets

