

# MediCap<sup>®</sup> MVR Lite

industry's most advanced and compact recordings.





#### **HD Medical Video Recorder**

Record Full HD video from your endoscope, arthroscope, surgical camera, ultrasound, etc. New state-of-the art recording platform records Full 1080p at 60 frames per second.



## Compact and robust

This rugged MVR features a stainless steel case and fluid resistant front panel. The compact size allows limitless mounting and installation options.



### **Intuitive Use**

It's so easy, no training is needed. Just plug it in and go. You'll be recording HD videos within minutes of opening the box.

Medical Imaging Made Easy™

# **Simple Features**



Medical Grade Designed, manufactured and certified for use in the OR / medical environment.



Remote control The included remote allows surgeons to record video, stills, and power on/ off at the touch of a button.



Photos Capture and store HD photos during your video recording.



E>

Foot Pedal Capture high quality HD video and HD photos to a USB drive using an endoscope trigger or foot pedal (optional).

# MediCap<sup>®</sup> MVR Lite Specs<sup>1</sup>

| Device Type                | <b>HD Medical Video Recorder</b><br>Records high-definition video and still<br>images from medical video sources to<br>external USB drive.   |
|----------------------------|--|
| Case/Buttons               | Stainless steel, fluid resistant   |
| User interface             | Overlay on external monitor<br>(record, pause, storage information)  |
| Video Inputs               | HDMI, DVI (with adapter)   |
| Video Outputs              | HDMI, DVI (with adapter)   |
| Video Formats              | H.264 / MPEG4 (1080p60 / 1080p50)  |
| Input video<br>resolutions | 1080p50/60, 1080i50/60, 720p50/60,<br>480p60, 576p50, 1280x1024,<br>1920x1200p50/60  |
| Image Formats              | JPG  |
| xternal Storage            | USB flash drive, USB hard drive  |
| USB Support                | USB3.0, high speed,<br>1 front connector, 2 rear connectors  |
| Footswitch/<br>Trigger     | 3.5mm mini jack. Compatible with standard footswitches (not included) and closed-contact camera triggers.                                    |
| Mechanical                 | Size: 160 x 140 x 40mm / 6.3 x 5.5 x 1.6"<br>Weight: 0.76kg / 1lb 11oz   |
| Environmental              | Storage: -40° to +85° C<br>Operating: -20° to +35° C   |
| Power                      | AC 110-240V, 50/60Hz, 0.4A,<br>3-prong AC jack   |
| Certifications             | CE 93/42/EEC, IEC 60601-1:2005,<br>ANSI/AAMI ES60601-1:2005<br>CAN/CSA C22.2 No.60601, 1:2008,<br>E361228, FCC Part 15, FDA Class 1,<br>RoHS |
| Kit Contents               | MediCap MVR Lite, USB flash drive,<br>power cable, DVI-to-HDMI adapter,<br>remote control, user guide  |
|                            |  |





Back Panel Connectors

# Also available...



#### MediCap® USB300

All in one full HD/3D recorder with built in monitor, network integration and patient information features

MediCapture, Inc. -North America

2250 Hickory Road, Suite 200 Plymouth Meeting, PA 19462, USA

(610) 238-0700





© 2015 MediCapture, Inc. MVRLITE-DS01-EN-150813-MC

## MediCapture Europe B.V.


Kraanpoort 10 6041 EG Roermond The Netherlands





#### MediCapture APAC

8FI- 5, No. 1, Jing Shan North Road Taipei City Chong Chen District 100 Taiwan, Taipei, Republic Of China



Medical Imaging Made Easy™

