



Rigid Optical Forceps


USA Elite System Rigid Optical Forceps are designed for use with the USA Elite System M3™ 12° or 30° Telescopes and all USA Elite System II Cystourethroscope Sheaths 21–25 Fr in size. These forceps also connect to all USA Elite System Resectoscope Sheaths without adapters.



E8213
Cup Tearing Biopsy Forceps obtain a large tissue sample.



E8216
Cup Cutting Biopsy Forceps obtain a tissue sample with a clean cut.



E8217
Stone-Crushing Forceps facilitate the removal of stones and stone fragments from the bladder.



E8215
Grasping Forceps' serrated jaws facilitate the removal of foreign bodies, stone fragments, and stents.

Features

- Snap-In Snap-Out™ locking mechanisms
- Complete line of rigid urology forceps
- Compatible with USA Elite System 21–25 Fr cystourethrosopes and resectoscopes
- Non obstructive jaws


Benefits

- Rapid, secure assembly/disassembly; perfect alignment of components
- Ideal for multiple applications
- System versatility
- Unobstructed field of view; easy to target location

Rigid Optical Forceps

- E8213** USA Elite System Rigid Optical Cup Tearing Biopsy Forceps
- E8215** USA Elite System Rigid Optical Grasping Forceps
- E8216** USA Elite System Rigid Optical Cup Cutting Biopsy Forceps
- E8217** USA Elite System Rigid Optical Stone-Crushing Forceps

USA Elite System™ Cystourethroscope Kits

	Kit Catalog No.	Color	Size	M3 Telescope	Sheaths	Albarran Bridge	Short Bridge	Visual Obturators	Forceps	Light Source & Light Guide	Accessories	Shipping Case
Master Kit	E-CUK	---	All	M3-30 (30°) M3-70 (70°)	E117, E121, E123, E125	E151, E152	E153, E159	E161, E163, E165	E8213, E8215	---	96E	99-0337
Office Kit	E-CUOK	---	All	M3-30 (30°) M3-70 (70°)	E117, E121	E151	E153, E154		E8213, 332F8	ALU-2B G93	96E	99-0338

™ Trademark of Circon Corporation or its Subsidiaries in the U.S. and/or other countries of the world