



USA Elite System™ Gynecologic Resectoscope Accessories

- GYNE-BR** USA Elite System Operative Hysteroscope Bridge, allows insertion of accessories up to 2.5 mm (7.5 Fr). Adapts to all USA Series Continuous Flow Resectoscope Inner Sheaths and USA Series Standard Flow Sheaths.
- GYA-732** Semi-Rigid Gynecologic Alligator Forceps, 7 Fr, 32 cm
- GYP-732** Semi-Rigid Gynecologic Punch Biopsy Forceps, 7 Fr, 32 cm
- GYS-732** Semi-Rigid Gynecologic Scissors, 7 Fr, 32 cm
- GE8216** USA Elite System Rigid Optical Cup Cutting Biopsy Forceps. Snaps directly into USA Elite System resectoscope sheaths. Compatible with USA Elite System M3G Telescopes
- GE8218** USA Elite System Rigid Optical Hook Scissors, double-action, hook jaws. Snaps directly into USA Elite System resectoscope sheaths. Compatible with USA Elite System M3G Telescopes



USA Elite System™ Gynecologic Resectoscope Kits and Sheath Systems


	Kit Catalog No.	Color Code	Size	M3G Telescope	Iglesias Working Element	Visual Obturators	Sheaths	Accessories Included
Rotating CFR Kit, 12°	GYR-525-12	Gold	8.5 mm (25.6 Fr)	M3G-12 (12°)	GEIWE	GERVO-CF25	GERIS-CF25, GEROS-CF25	Shipping Case, 75U, GRE, GRB
Rotating CFR Kit, 30°	GYR-525-30	Gold	8.5 mm (25.6 Fr)	M3G-30 (30°)	GEIWE	GERVO-CF25	GERIS-CF25, GEROS-CF25	Shipping Case, 75U, GRE, GRB
Rotating CFR Kit, 12°	GYR-527-12	Gray	9.2 mm (27.6 Fr)	M3G-12 (12°)	GEIWE	GERVO-CF27	GERIS-CF27, GEROS-CF27	Shipping Case, 75U, GRE, GRB
Rotating CFR Kit, 30°	GYR-527-30	Gray	9.2 mm (27.6 Fr)	M3G-30 (30°)	GEIWE	GERVO-CF27	GERIS-CF27, GEROS-CF27	Shipping Case, 75U, GRE, GRB
Rotating CFR System	GERSV-CF25	Gold	8.5 mm (27.6 Fr)	---	---	GERVO-CF25	GERIS-CF25, GEROS-CF25	---
Rotating CFR System	GERSV-CF27	Gray	8.5 mm (27.6 Fr)	---	---	GERVO-CF27	GERIS-CF27, GEROS-CF27	---
CFR Kit, 12°	GYE-525-12	Yellow	8.5 mm (25.6 Fr)	M3G-12 (12°)	GEIWE	GEVO-CF25	GEIS-HCF25, GEOS-HCF25	Shipping Case, 75U, GRE, GRB
CFR Kit, 30°	GYE-525-30	Yellow	8.5 mm (25.6 Fr)	M3G-30 (30°)	GEIWE	GEVO-CF25	GEIS-HCF25, GEOS-HCF25	Shipping Case, 75U, GRE, GRB
CFR Kit, 12°	GYE-527-12	White	9.2 mm (27.6 Fr)	M3G-12 (12°)	GEIWE	GEVO-CF27	GEIS-HCF27, GEOS-HCF27	Shipping Case, 75U, GRE, GRB
CFR Kit, 30°	GYE-527-30	White	9.2 mm (27.6 Fr)	M3G-30 (30°)	GEIWE	GEVO-CF27	GEIS-HCF27, GEOS-HCF27	Shipping Case, 75U, GRE, GRB
CFR System, Round Tip	GEST-HCF25	Yellow	8.5 mm (25.6 Fr)	---	---	GEVO-CF25	GEIS-HCF25, GEOS-HCF25	---
CFR System, Round Tip	GEST-HCF27	White	9.2 mm (27.6 Fr)	---	---	GEVO-CF27	GEIS-HCF27, GEOS-HCF27	---
Standard Kit, 12°	GYE-324-12	Yellow	8 mm (24 Fr)	M3G-12 (12°)	GEIWE	EVOR-24	EMRS-24B	Shipping Case, 75U, GRE, GRB
Standard Kit, 30°	GYE-324-30	Yellow	8 mm (24 Fr)	M3G-30 (30°)	GEIWE	EVOR-24	EMRS-24B	Shipping Case, 75U, GRE, GRB

™ Trademark of Circon Corporation or its Subsidiaries in the U.S. and/or other countries of the world


Arth
 Gastro
 Gyn
 Lap
 Urology
 Video
 Illum
 Ordering
 Index

USA Series™ Gynecologic Resectoscope Accessories

- GYN-BR** USA Series Operative Hysteroscope Bridge, allows insertion of accessories up to 2.5 mm (7.5 Fr). Adapts to all USA Series Continuous Flow Resectoscope Inner Sheaths and USA Series Standard Flow Sheaths.
- GYA-732** Semi-Rigid Gynecologic Alligator Forceps, 7 Fr, 32 cm
- GYP-732** Semi-Rigid Gynecologic Punch Biopsy Forceps, 7 Fr, 32 cm
- GYS-732** Semi-Rigid Gynecologic Scissors, 7 Fr, 32 cm
- GY8216** USA Series Rigid Optical Cup Cutting Biopsy Forceps. Snaps directly into USA Series resectoscope sheaths. Compatible with USA Elite System M3G Telescopes.
- GY8218** USA Series Rigid Optical Hook Scissors, double-action, hook jaws. Snaps directly into USA Elite System resectoscope sheaths. Compatible with USA Elite System M3G Telescopes.



USA Series™ Gynecologic Continuous Flow Resectoscope Kits and Sheath Systems*

	Kit Catalog No.	Color	Size	M3G Telescope	Iglesias Working Element	Visual Obturators	Sheaths	Accessories Included
CFR Kit, 12°	GYH-525-12	Yellow	8.5 mm / 25.6 Fr	M3G-12 (12°)	GIWE	25-VO-CFR	G25-IS-HCF, G25-OS-HCF	Shipping Case, 75U, GRE, GRB
CFR Kit, 30°	GYH-525-30	Yellow	8.5 mm / 25.6 Fr	M3G-30 (30°)	GIWE	25-VO-CFR	G25-IS-HCF, G25-OS-HCF	Shipping Case, 75U, GRE, GRB
CFR Kit, 12°	GYH-527-12	White	9.2 mm / 27.6 Fr	M3G-12 (12°)	GIWE	27-VO-CFR	G27-IS-HCF, G27-OS-HCF	Shipping Case, 75U, GRE, GRB
CFR Kit, 30°	GYH-527-30	White	9.2 mm / 27.6 Fr	M3G-30 (30°)	GIWE	27-VO-CFR	G27-IS-HCF, G27-OS-HCF	Shipping Case, 75U, GRE, GRB
CFR System, Round Tip	G25-ST-HCF	Yellow	8.5 mm / 25.6 Fr	---	---	25-VO-CFR	G25-IS-HCF, G25-OS-HCF	---
CFR System, Round Tip	G27-ST-HCF	White	9.2 mm / 27.6 Fr	---	---	27-VO-CFR	G27-IS-HCF, G27-OS-HCF	---


* USA Series Gynecologic Standard Flow Resectoscope components are sold individually.

™ Trademark of Circon Corporation or its Subsidiaries in the U.S. and/or other countries of the world

USA Series™ Continuous Flow Operative Hysteroscopes¹
Designed with Multi-Functional Features for Operative Hysteroscopy

SNAP-IN SNAP-OUT™

A true Continuous Flow Operative Hysteroscopy System with advanced M3G™ optic design and superior flow characteristics.



GY7-CFH Single Channel System with M3G Telescope

GY8-CFH Double Channel System with M3G Telescope

Features

- True continuous flow
- Non-communicating channels
- 7 or 8 mm outer diameter with one or two 2.3 mm channels
- Snap-In Snap-Out™ locking mechanism
- Compatible with USA Elite System™ M3G™ Telescope optics
- Available in 0°, 12°, and 30° direction of view
- 1 cm markings on sheath


Benefits

- Maintains view clear of blood/debris and sustains constant pressure in uterus
- Separate optics, working/distension channels enable greater control of operative procedure
- Allows simultaneous use of instruments/laser fibers
- Quick, easy, secure assembly
- Excellent clarity, depth of field, resolution, and video compatibility
- Choice of telescope's direction of view
- Used to sound the uterine cavity

- GY7-CFH USA Series Continuous Flow Operative Hysteroscope System (7 mm) includes:**
- GYB-7** USA Series Continuous Flow Operative Bridge/Obturator, single 2.3 mm (7 Fr) channel
 - GYC-7** USA Series Continuous Flow Outer Sheath, 7 mm
- GY8-CFH USA Series Continuous Flow Operative Hysteroscope System (8 mm) includes:**
- GYB-8** USA Series Continuous Flow Operative Bridge/Obturator, dual 2.3 mm (7 Fr) channels
 - GYC-8** USA Series Continuous Flow Outer Sheath, 8 mm

COMPANION PRODUCTS

CS-B612	Accessory Port Seals, non-latex
GYA-732	Semi-Rigid Gynecologic Alligator Forceps, 7 Fr, 32 cm
GYP-732	Semi-Rigid Gynecologic Punch Biopsy Forceps, 7 Fr, 32 cm
GYS-732	Semi-Rigid Gynecologic Scissors, 7 Fr, 32 cm
M3G-0	USA Elite System M3G Direct Telescope (0°)
M3G-12	USA Elite System M3G Operative Telescope (12°)
M3G-30	USA Elite System M3G Foroblique Telescope (30°)
M3G-0A	USA Elite System M3G Autoclavable Direct Telescope (0°)
M3G-12A	USA Elite System M3G Autoclavable Operative Telescope (12°)
M3G-30A	USA Elite System M3G Autoclavable Foroblique Telescope (30°)



Excellent continuous flow through isolated single or double working channels.

For more information on the USA Elite System M3 and M3G Autoclavable Telescopes, please refer to p. GY-15.

¹ U.S. Patent No. 5,320,091
 ™ Trademark of Circon Corporation or its Subsidiaries in the U.S. and/or other countries of the world