



# STERILE DISPOSABLE BLADES

## STERILE DISPOSABLE BLADES




Used with Dyonics® EP-1® Motor Drives  
(REF #s 7205354 and 7205355)








Used with the PS3500EP® Motor Drive  
(REF # 3476)

### Orbit® Rotatable Curved Disposable Blades




Packaged 6 per box/sterile


	EP-1 REF #	PS3500EP REF #	Description	Color
	7205320*	7204700	4.5 mm Orbit Incisor®, curve length 17 mm	Lime Green
	7205321	7204927	4.5 mm Orbit Synovator®, curve length 17 mm	Forest Green
	7205322	7204929	4.5 mm Orbit RazorCut®, curve length 17 mm	Pink
	7205323	7204928	4.5 mm Orbit Full Radius, curve length 17 mm	Yellow
	7205687	N/A	4.5 mm Orbit Incisor Plus, curve length 17 mm	Violet

\* Item Shown

STERILE DISPOSABLE BLADES






Used with Dyonics® EP-1® Motor Drives  
(REF #s 7205354 and 7205355)




Used with the PS3500EP® Motor Drive  
(REF # 3476)


Curved Disposable Blades

Packaged 6 per box/sterile

	EP-1® REF #	PS3500EP® REF #	Description	Color
	7205332	4136	4.5 mm Curved Incisor®, curve length 17 mm	Lime Green
	7205334	4138	4.5 mm Concave Synovator®, curve length 17 mm	Forest Green
	N/A	4349	4.5 mm Concave Full Radius, curve length 17 mm	Yellow



Used with Dyonics® EP-1® Motor Drives  
(REF #s 7205354 and 7205355)



Used with the PS3500EP® Motor Drive  
(REF # 3476)


**Straight Disposable Blades**

Packaged 6 per box/sterile


	EP-1 REF #	PS3500EP REF #	Description	Color
	7205316	3446	4.5 mm TurboWhisker®	Navy Blue
	N/A	3449	4.0 mm TurboCutter®	Turquoise
	7205312	4223	3.5 mm Incisor®	Medium Grey
	7205313	3810	4.5 mm Incisor	Lime Green
	7205314	4191	5.5 mm Incisor	Mocha
	7205310	3826	4.5 mm Synovator®	Forest Green
	7205311	4190	5.5 mm Synovator	Canary Yellow
	7205317	4222	3.5 mm RazorCut®	Purple
	7205318	4280	4.5 mm RazorCut	Pink
	7205319	4284	5.5 mm RazorCut	Light Pink
	7205308	3439	3.5 mm Cutter	Red
	7205309	3440	4.5 mm Cutter	Blue
	N/A	3441	4.5 mm Trimmer	Green
	7205305	3442	3.5 mm Full Radius	Beige
	7205306	3443	4.5 mm Full Radius	Yellow
	7205307	3444	5.5 mm Full Radius	Orange
	7205315	3529	4.5 mm TurboTrimmer®	Powder Blue
	7205345	7205398	4.5 mm Incisor Plus	Violet
	7205459	N/A	5.5 mm Incisor Plus	Teal

D4

STERILE DISPOSABLE BURRS








Used with Dyonics® EP-1® Motor Drives  
(REF #s 7205354 and 7205355)




Used with the PS3500EP® Motor Drive  
(REF # 3476)

Straight Disposable Burrs


Packaged 6 per box/sterile


	EP-1 REF # (Wide Mouth)	PS3500EP REF # Wide Mouth	Slotted	Description	Color
	7205324	7205660	3450	4.0 mm Abrader	Aqua
	7205325	7205661	3451	5.5 mm Abrader	Black
	7205326	7205662	3452	4.0 mm Acromionizer	Mauve
	7205327	7205663	3453	5.5 mm Acromionizer	Brown
	7205330	7205666	4304	4.0 mm StoneCutter™	Maroon
	7205331	7205667	4305	5.5 mm StoneCutter	Olive
	7205328	7205664	4306	4.0 mm NotchBlaster®	Lilac
	7205329	7205665	4307	5.5 mm NotchBlaster	Peach
	7205668	N/A	N/A	4.0 mm AcromioBlaster™	Sage Green
	7205669	N/A	N/A	5.5 mm AcromioBlaster	Brick Red

Wide Mouth Burr



Slotted Style Burr





**Helical Disposable Burr**

Designed with an inner core featuring helical flutes, a spherical sharpened end, and enclosed by an outer sleeve with a slanted end. Quickly abrades bone, cuts soft tissue, evacuates the delivery site, reduces clogging and—as a result—saves valuable O.R. time. Packaged 6 per box/sterile

	<b>EP-1° REF #</b>	<b>Description</b>	<b>Color</b>
	<b>7205727</b>	<b>4.5 mm HeliCut™</b>	<b>Slate</b>

## LONG DISPOSABLE BLADES

**Dyonics® EP-1® Long Shaver Blades**

18 cm working length. Packaged 3 per box/sterile


REF #	Description	Color
7205335	4.5 mm curved, concave, Full Radius	Yellow
7205336	4.5 mm curved, convex, Full Radius	Yellow
7205341	4.5 mm straight, Full Radius	Yellow
7205343	4.5 mm straight, Incisor®	Lime Green

**PS3500EP® Long Shaver Blades**

18 cm working length. Packaged 3 per box/sterile

REF #	Description	Color
7205358	4.5 mm curved, concave, Full Radius	Yellow
7205359	4.5 mm curved, convex, Full Radius	Yellow
7205413	4.5 mm straight, Full Radius	Yellow
7205414	4.5 mm straight, Incisor	Lime Green

STERILE DISPOSABLE MINI BLADES/BURRS



**Sterile Disposable Mini Blades and Burrs**

Used with the Dyonics® EP-1® Mini Motor Drive (REF # 7205357), PS3500EP® Mini Motor Drive Unit (REF # 3391)

Indications: Small joint arthroscopy and Functional Endoscopic Sinus Surgery (FESS)


Packaged 6 per box/sterile

	REF #	Description	Color
	3606**	2.9 mm Cutter, 7 cm length	Beige
	3407**	2.0 mm TurboWhisker®, 7 cm length	Yellow
	3418**	2.9 mm TurboWhisker®, 7 cm length	Green
	3410	2.0 mm Full Radius, 7 cm length	Blue
	3419	2.9 mm Full Radius, 7 cm length	Red
	3553	2.9 mm Barrel Abrader®, 7 cm length	Purple
	3530	2.9 mm Abrader, 7 cm length	Orange
	3819	3.5 mm Abrader, 7 cm length	Aqua

D8

\*\* Contraindicated for the FESS procedure.

## STERILE DISPOSABLE MINI BLADES—LONG AND CURVED







Used with Dyonics® EP-1® Mini Motor Drives (REF # 7205357) and PS3500EP® Mini Motor Drive Unit (REF # 3391)

- Curved blade allows greater access to anatomy
- Blade Window Lock feature sets blade to lock at any position in window, as selected by the surgeon, protecting tissue from blade edges and controlling suction when not in cutting mode
- Long and curved disposable mini blades are identical to the current mini blades, with one clear difference: they are one centimeter longer (8 cm)

### Long and Curved Disposable Mini Blades

Indications: Small joint arthroscopy and Functional Endoscopic Sinus Surgery (FESS)  
Packaged 6 per box/sterile

	REF #	Description	Color
	7205474	3.5 mm Full Radius, 8 cm length	Plum
	7205477	3.5 mm RazorCut®, 8 cm length	Pumpkin (orange)
	7205558	3.5 mm Curved Concave, Incisor®, 8 cm length	Brick Red
	7205476	3.5 mm Incisor® Plus, 8 cm length	Sage Green