


Optical Urethrotome

with **HOPKINS**® Straight Forward Telescope 27015 A/27005 A/AA,
for treatment of strictures, size 21 Fr.

STORZ
KARL STORZ — ENDOSKOPE



- 27068 D **Urethrotome Sheath**, 21 Fr.,
with channel for Filiform-Bougies and 2 LUER-lock
adaptors.
- 27068 DO **Obturator** for urethrotome sheath 21 Fr., 27068 D.
- 27068 CD **Telescope Bridge** with channel for instruments 5 Fr.
- 27068 F **Supplementary Sheath**, half-round,
to insert a ballon catheter,
fits urethrotome sheath 27068 D.
- 27068 S **Supplementary Sheath**,
for continuous irrigation and suction,
fits urethrotome sheath 27068 D.