


Stone Punch

with **HOPKINS III**® Straight Forward Telescope 27005 AR, enlarged view
or **HOPKINS III**® Forward-Oblique Telescope 27005 B/BA, enlarged view.



MAUERMAYER Stone Punch, 25 Fr.,
consisting of:

27005 AR **HOPKINS III**® Straight Forward Telescope 0°,
enlarged view, diameter 4 mm, **autoclavable.**
Fiber optic light transmission incorporated.
Color code: green.

or

27005 B **HOPKINS III**® Forward-Oblique Telescope 30°,
Enlarged View, diameter 4 mm.
Fiber optic light transmission incorporated.
Color code: red.

or

27005 BA **HOPKINS III**® Forward-Oblique Telescope 30°,
Enlarged View, diameter 4 mm, **autoclavable.**
Fiber optic light transmission incorporated.
Color code: red.

27077 A **Punch Working Element.**

27077 B **Punch Sheath** with Central Valve, including
connecting tubes for suction and irrigation, 25 Fr.,
straight beak with obturator 27077 DO.

27077 G **Insert Tube** with channel for flexible instruments, 7 Fr.,
with atraumatic beak for urethroscopy.

27077 H **Insert Tube** with channel for flexible probes
for electrohydraulic lithotripsy, 27080 FA.

Larger stones can be broken up until they are in the size for removal
through the punch only.

The grasping of stones in smaller diverticula is possible because the
length of the instrument is convenient.

Evacuators see page URO-ACC 5.

CALCUTRIPT Electrohydraulic System and Probes

see page URO-Units 3 A.