


# Resectoscope Sheaths with rotatable inner tube

for continuous flow irrigation and suction.

NEW

**Special features:**

- Ceramic insulating beak insert for burn resistance.
- Ceramic insulation on inner tube for easy service.
- For use with standard working elements 27040 and 27050.
- Rotatable inner tube.



**26 Charr.**

27050 SL **Resectoscope Sheath** including connecting tubing for in- and outflow, 26 Fr., oblique beak, **rotatable** inner tube with ceramic insulation, for use with working elements 27040 C/D/E and 27050 C/D/E/EN. Color code: yellow.

27050 XA **Inner Tube**, rotatable, with ceramic insulation, for use with resectoscope sheath 27050 SL, 26 Fr.

**28 Charr.**

27050 SM **Resectoscope Sheath** including connecting tubing for in- and outflow, 28 Fr., oblique beak, **rotatable** inner tube with ceramic insulation, for use with working elements 27040 C/D/E and 27050 C/D/E/EN. Color code: black.


27050 XB **Inner Tube**, rotatable, with ceramic insulation, for use with resectoscope sheath 27050 SM, 28 Fr.

27094 BY **Adaptor** for use with optical forceps 27072 BL/BS/MC and with resectoscope sheaths 27040 SL/SD/SM/SE, 27050 SL/SM.

Standard Obturators and Deflecting Obturators see page RES-OB 1 A.

# Accessories

Obturators for use with resectoscope sheaths 24–28 Fr.




**Standard Obturators**

- 27040 OA For use with 27/28 Fr. sheaths 27040 AO, 27141 AO, 27040 SE/SM, 27050 SM, color code: black
- 27040 OB For use with 27 Fr. sheaths 27040 AK, 27141 AK, color code: black.
- 27040 OC For use with 24/26 Fr. sheaths 27040 BO, 27141 BO, 27040 SD/SL, 27050 SL, color code: yellow.
- 27040 OD For use with 24 Fr. sheaths 27040 BK, 27141 BK, color code: yellow.



**Deflecting Obturators**

- 27048 AK For use with 27 Fr. sheaths 27040 AK, 27141 AK, color code: black.
- 27048 CO For use with 27/28 Fr. sheaths 27040 AO, 27141 AO, 27040 SE/SM, 27050 SM, color code: black.
- 27048 BK For use with 24 Fr. sheaths 27040 BK, 27141 BK, color code: yellow.
- 27048 CK For use with 24/26 Fr. sheaths 27040 BO, 27141 BO, 27040 SD/SL, 27050 SL, color code: yellow.



**LEUSCH Atraumatic Distending Obturators**



- 27049 AK For use with 27 Fr. sheaths 27040 AK, 27141 AK, color code: black.
- 27049 AO For use with 27/28 Fr. sheaths 27040 AO, 27141 AO, 27040 SE/SM, 27050 SM, color code: black.
- 27049 BK For use with 24 Fr. sheaths 27040 BK, 27141 BK, color code: yellow.
- 27049 BO For use with 24/26 Fr. sheaths 27040 BO, 27141 BO, 27040 SD/SL, 27050 SL, color code: yellow.

The above described distending obturators allow for atraumatic insertion of the sheaths.