

Resectoscope Sheaths

Metal-Sheaths with LUER-lock stopcock including connecting tube for irrigation.




- 27040 AO 27 Fr., oblique beak, with obturator 27040 OA, color code: black.
27040 BO 24 Fr., oblique beak, with obturator 27040 OC, color code: yellow.




- 27040 AK 27 Fr., short beak, with obturator 27040 OB, color code: black.
27040 BK 24 Fr., short beak, with obturator 27040 OD, color code: yellow.


Metal-Sheaths with Central Valve including connecting tubes for in- and outflow.



- 27041 AO 27 Fr., oblique beak, with obturator 27040 OA, color code: black.
27041 BO 24 Fr., oblique beak, with obturator 27040 OC, color code: yellow.



- 27041 AK 27 Fr., short beak, with obturator 27040 OB, color code: black.
27041 BK 24 Fr., short beak, with obturator 27040 OD, color code: yellow.



- 27040 AG MATOUSCHEK **Lubricant Applicator**, for use with
resectoscope sheaths 27040 AK/AO and 27041 AK/AO, 27 Fr.
27040 BG **Same**, for use with resectoscope sheaths 27040 BK/BO
and 27041 BK/BO, 24 Fr.

The MAUERMAYER Stone Punch 27043 M (see page RES-SC 2) may be used
with all **KARL STORZ** resectoscope sheaths.