


Continuous-Flow LASER-Cystoscope, 21 Fr.




27026 VT

27026 VL


27026 VT **LASER-Cystoscope-Urethroscope Sheath** 21 Fr.,
with continuous flow irrigation and suction,
with obturator 27026 VO and 2 LUER-lock adaptors.

27026 VL **Telescope Bridge**,
with 8 Fr. working channel, with stopcock.



27005 BA/FA

HOPKINS II® telescopes, autoclavable,
for use with LASER Cystoscope-Urethroscope Sheath 21 Fr.:



27005 FA **HOPKINS II**® Telescope 12°,
enlarged view, diameter 4 mm, **autoclavable**.
Fiber optic light transmission incorporated.
Color code: black.



27005 BA **HOPKINS II**® Forward-Oblique Telescope 30°,
enlarged view, diameter 4 mm, **autoclavable**.
Fiber optic light transmission incorporated.
Color code: red.