


Neonate Optical Urethrotome

for use with **HOPKINS**® Straight Forward Telescope 27017 A,
for transurethral treatment of strictures in children, size 10 Fr.


Working length 11 cm.



27045 EK **Working Element**, Set consisting of:
27045 E 1 Working Element.
27047 2 Knives.

27047 B **Urethrotome Sheath** with LUER-lock stopcock, 10 Fr.,
with obturator 27047 BO and 2 LUER-lock adaptors.

27047 F **Telescope Bridge** with one instrument channel.



Cutting Knives:



27047 TT **Cold Knife**, straight.



27047 V **Cold Knife**, round.



27047 W **Cold Knife**, sickle-shaped.



27047 X **Cold Knife**, hook-shaped.



280 **Protection Tube**, for sterilization and storage
of electrodes, curettes and knives.