

LASER CYSTOSCOPE



23Fr Continuous Flow Laser Cystoscope System

For use with the Green Light Laser™ for PVP Procedures and other side fire laser systems

- 502-532-030** 4.0 mm 30° Autoclavable Cystoscope with Integrated Green Light Filter
- 502-880-300** Continuous Flow Laser Cysto Sheath Kit (includes 502-880-299 and 502-880-303)
- 502-880-299** Continuous Flow Laser Cysto Sheath (Outer)
- 502-880-303** Continuous Flow Laser Cysto Sheath (Inner)
- 502-880-301** Visual Obturator, Laser Cystoscope
- 502-880-302** Sweep Guide, Laser Cystoscope

URETHROTOME



Visual Urethrotome

For use with 0° or 12° Cystoscopes to cut strictures in the urethra.



- 502-880-401** Working Element — Passive (Iglesias)
- 502-880-521** 21Fr Visual Urethrotome Sheath
- 502-880-525** 25Fr Visual Urethrotome Outer Sheath
- 502-880-524** 21Fr Visual Urethrotome Working/Optical Bridge
- 502-880-522** 21Fr Urethrotome Standard Obturator
- 502-880-523** 21Fr Urethrotome Visual Obturator
- 502-880-526** 25Fr Balloon Catheter Sheath

Urethrotome Cold Knives

For use with 0° or 12° Cystoscopes




- 502-880-501** Cold Knife, Half Moon
- 502-880-502** Cold Knife, Straight
- 502-880-503** Cold Knife, Serrated
- 502-880-504** Cold Knife, Hook
- 502-880-505** Cold Knife, Loop

SEMI-RIGID URETEROSCOPES

SRU-6/SRU-6X Ureteroscopes

- » Autoclavable, high-density fused quartz bundle for high-resolution images
- » 6.9Fr atraumatic SuperSlim tip for easy insertion
- » 3.4Fr and 2.5Fr dual working channels for consistent, maximized flow rate



SRU-6 Series Autoclavable Semi-Rigid Ureteroscopes


- 502-880-330** SRU-6 Autoclavable Semi-Rigid Ureteroscope, 6.9/10Fr, 33 cm
- 502-880-430** SRU-6X Autoclavable Semi-Rigid Ureteroscope, 6.9/10Fr, 43 cm

Parts 502-880-330 and 502-880-430 include the following products:

- 504-999-001** Cleaning Brush for Semi-Rigid Ureteroscope
- 233-032-890** Sterilization/Storage Tray for Semi-Rigid Ureteroscope

Disposable Baskets

For use with SRU-6, SRU-6X



- 504-803-070** 3Fr Flat Wire Basket, 70 cm
- 504-803-090** 3Fr Flat Wire Basket, 90 cm
- 504-803-115** 3Fr Flat Wire Basket, 115 cm
- 504-882-115** 2.4Fr Nitinol Basket, 115 cm
- 504-833-070** 3Fr 3-Wire Grasping Forcep, 70 cm
- 504-833-090** 3Fr 3-Wire Grasping Forcep, 90 cm
- 504-833-115** 3Fr 3-Wire Grasping Forcep, 115 cm
- 504-883-070** 3Fr Nitinol Basket, 70 cm
- 504-883-090** 3Fr Nitinol Basket, 90 cm
- 504-883-115** 3Fr Nitinol Basket, 115 cm