

ARTHROSCOPES



IDEAL EYES HD Arthroscopes

Improved Center To Edge Resolution

- » Resolution and brightness are enhanced across the entire field of view for optimal image quality
- » Enhanced outer radius line pairs per millimeter enhanced from 11.4 horizontal lines on current scopes to 12.8 horizontal lines on IDEAL EYES scopes

New Magnetic Drive Design For The Focus Ring

- » Optimal moisture resistance
- » Improved reliability during all sterilization methods
- » Precise focusing of the C-Mount scope with one hand

IDEAL EYES Autoclavable Arthroscopes C-Mount

Part Number	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Locking Mechanism
502-704-010	4.0 mm	0°	105°	140 mm	Standard	Speed-Lock
502-704-030	4.0 mm	30°	105°	140 mm	Standard	Speed-Lock
502-704-070	4.0 mm	70°	105°	140 mm	Standard	Speed-Lock
502-704-130	4.0 mm	30°	105°	140 mm	Reverse Cant	Speed-Lock
502-804-030	4.0 mm	30°	105°	165 mm	Standard	Speed-Lock
502-804-070	4.0 mm	70°	105°	165 mm	Standard	Speed-Lock
502-804-130	4.0 mm	30°	105°	165 mm	Reverse Cant	Speed-Lock
502-727-010	2.7 mm	0°	105°	120 mm	Standard	Speed-Lock
502-727-030	2.7 mm	30°	105°	120 mm	Standard	Speed-Lock
502-727-070	2.7 mm	70°	105°	120 mm	Standard	Speed-Lock
502-627-010	2.7 mm	0°	80°	75 mm	Standard	J-Lock
502-627-030	2.7 mm	30°	80°	75 mm	Standard	J-Lock
502-627-070	2.7 mm	70°	80°	75 mm	Standard	J-Lock
502-627-130	2.7 mm	30°	80°	75 mm	Reverse Cant	J-Lock
502-723-010	2.3 mm	0°	80°	72 mm	Standard	J-Lock
502-723-030	2.3 mm	30°	80°	72 mm	Standard	J-Lock
502-723-070	2.3 mm	70°	80°	72 mm	Standard	J-Lock
502-723-130	2.3 mm	30°	80°	72 mm	Reverse Cant	J-Lock
502-719-010	1.9 mm	0°	65°	58 mm	Standard	J-Lock
502-719-030	1.9 mm	30°	65°	58 mm	Standard	J-Lock