

Autoclavable C-Mount Arthroscopes ARTHROSCOPY



Autoclavable C-Mount Arthroscopes

- » Direct coupler mounted arthroscopes
- » Unique coupler mount design protects against moisture and fogging
- » Triple tube design for resistance against bending and lens cracking

Autoclavable Arthroscopes C-Mount


Part Number	Outer Diameter	Direction of View	Field of View	Working Length	Light Post Direction	Locking Mechanism
502-527-010	4.0 mm	0°	105°	140 mm	Standard	Speed-Lock
502-527-030	4.0 mm	30°	105°	140 mm	Standard	Speed-Lock
502-527-070	4.0 mm	70°	105°	140 mm	Standard	Speed-Lock
502-527-130	4.0 mm	30°	105°	140 mm	Reverse Cant	Speed-Lock
502-528-010	2.7 mm	0°	105°	120 mm	Standard	Speed-Lock
502-528-030	2.7 mm	30°	105°	120 mm	Standard	Speed-Lock
502-528-070	2.7 mm	70°	105°	120 mm	Standard	Speed-Lock
502-526-010	2.7 mm	0°	80°	75 mm	Standard	Speed-Lock
502-526-030	2.7 mm	30°	80°	75 mm	Standard	Speed-Lock
502-526-070	2.7 mm	70°	80°	75 mm	Standard	Speed-Lock
502-526-130	2.7 mm	30°	80°	75 mm	Reverse Cant	Speed-Lock
502-524-030	1.9 mm	30°	65°	58 mm	Standard	Speed-Lock

Hardware for 2.7 mm Arthroscopes **ARTHROSCOPY**

2.7 mm Hardware—Speed-Lock

4.0 mm Series Speed-Lock Hardware

For 2.7 mm, 120 mm Arthroscopes




- 343-331-430 4.0 mm Cannula, 1 Rotating Stopcock
- 343-331-450 4.0 mm Cannula, 2 Rotating Stopcocks
- 343-331-410 Pencil Tip Obturator for 4.0 mm Cannula
- 343-331-411 Blunt Tip Obturator for 4.0 mm Cannula
- 343-331-420 Trocar for 4.0 mm Cannula

2.7 mm Stubby Hardware—Speed-Lock

3.2 mm Series Speed-Lock Hardware


For 2.7 mm, 75 mm Arthroscopes
(Low-Flow)



- 502-344-530 3.2 mm Cannula, 1 Rotating Stopcock
- 502-344-550 3.2 mm Cannula, 1 Fixed Stopcock
- 502-344-540 3.2 mm Cannula, 2 Fixed Stopcocks
- 502-344-560 3.2 mm Cannula, 2 Rotating Stopcocks
- 502-344-500 Pencil Tip Obturator for 3.2 mm Cannula
- 502-344-510 Blunt Tip Obturator for 3.2 mm Cannula
- 502-344-520 Trocar for 3.2 mm Cannula

4.0 mm Series Speed-Lock Hardware

For 2.7 mm, 75 mm Arthroscopes
(High-Flow)




- 502-244-530 4.0 mm Cannula, 1 Rotating Stopcock
- 502-244-550 4.0 mm Cannula, 1 Fixed Stopcock
- 502-244-540 4.0 mm Cannula, 2 Fixed Stopcocks
- 502-244-560 4.0 mm Cannula, 2 Rotating Stopcocks
- 502-244-500 Pencil Tip Obturator for 4.0 mm Cannula
- 502-244-510 Blunt Tip Obturator for 4.0 mm Cannula
- 502-244-520 Trocar for 4.0 mm Cannula

2.7 mm Stubby Hardware—J-Lock

3.2 mm Series J-Lock Hardware

For 2.7 mm, 75 mm Arthroscopes
(Low-Flow)



- 502-713-530 3.2 mm Cannula, 1 Rotating Stopcock
- 502-713-550 3.2 mm Cannula, 1 Fixed Stopcock
- 502-713-560 3.2 mm Cannula, 2 Rotating Stopcocks
- 502-713-540 3.2 mm Cannula, 2 Fixed Stopcocks
- 502-713-500 Pencil Tip Obturator for 3.2 mm Cannula
- 502-713-510 Blunt Tip Obturator for 3.2 mm Cannula
- 502-713-520 Trocar for 3.2 mm Cannula

TRAYS

Arthroscopic Trays



Autoclavable Arthroscope Sterilization Tray



Autoclavable Arthroscope and Hardware Sterilization Tray



Arthroscope and Camera Sterilization Tray



Endoscope and Camera Sterilization Tray

Autoclavable Arthroscope Sterilization Tray	233-032-104	Holds two eyepiece arthroscopes L= 19" W= 8-1/2" H= 2-3/4"
Autoclavable Arthroscope and Hardware Sterilization Tray	233-032-116	Holds two C-mount arthroscopes and hardware L= 10-9/16" W= 6-1/4" H= 1-5/16"
Arthroscope and Camera Sterilization Tray	233-032-105	Holds three arthroscopes, one camera, and one light cable L= 19-1/2" W= 8-11/16" H= 2-15/16"
Endoscope and Camera Sterilization Tray	233-032-107	Holds two arthroscopes, two laparoscopes, one camera, and one light cable L= 19" W= 8-1/2" H= 2-3/4"

Instrument Trays



Manual Instruments Tray



Shaver Handpiece Tray



Cutter Accessory Tray



Glick Hip Kit Tray

Sterilization Trays for Arthroscopy Manual Instruments	242-000-012	12 Instrument Sterilization Tray
	242-000-024	23 Instrument Sterilization Tray
Shaver Handpiece Tray	272-700-000	Holds one handpiece L= 14" W= 9" H= 2"
Cutter Accessory Tray	272-707-000	Accommodates all cutter and bur hardware
Glick Hip Kit Tray	242-720-025	Sterilization Tray