

POWERED INSTRUMENTS



CROSSFIRE Integrated Arthroscopy Resection System

- » Increased cutting efficiency with all current Formula cutter and bur disposables
- » Improved RF ablation
- » The two-in-one console and footswitch saves space in the OR
- » iSWITCH wireless platform compatible
- » Completely customizable interface eliminates the need to switch individual settings between cases
- » Large, bright and easy to use LCD screen
- » 3 connection ports for Formula Shaver, SERFAS Energy probe, and CROSSFIRE footswitch
- » RFID tag recognition for all cutters, burs, and SERFAS Energy probes



CROSSFIRE

475-000-000 CROSSFIRE Console
475-000-100 CROSSFIRE Footswitch

iSWITCH



iSWITCH Universal Wireless Foot Control

Stryker developed a wireless universal footswitch that can eliminate the inconveniences and hazards standard footswitches and cords may create. The footswitch uses a wireless USB connection and has the ability to control the Stryker CORE and SERFAS Energy footswitch functions.

- » Completely sealed off from liquids
- » Long lasting battery life

iSWITCH Universal Wireless Receiver and Footswitch

277-200-100 iSWITCH Receiver
277-100-100 iSWITCH Universal Wireless Footswitch

iSWITCH Accessories

277-600-100 iSWITCH Valleylab Integration Kit
6000-001-020 2 Meter Firewire Cable
277-400-100 TPS Connector
277-500-100 iSWITCH Clean-and-Dry Bags
277-300-100 iSWITCH Battery Pack

RF ABLATION SYSTEM

CROSSFIRE Multiple Options

- » The hand-controlled SERFAS Energy Probes offer an ergonomic hand piece with multiple button functions. Each hand piece allows the surgeon to operate the cut and coagulation modes as well as change the power settings. The hand-controlled feature can now fully replace the need for a footswitch.
- » The SERFAS Energy System offers both suction and non-suction probes.
- » Each disposable SERFAS Energy Probe is designed with an integrated cable, allowing for easy case preparation and clean-up.
- » The SERFAS Energy Probes come equipped with a memory chip that contains specific optimized settings for each probe.



CROSSFIRE

- 475-000-000 CROSSFIRE Console
- 475-000-100 CROSSFIRE Footswitch



SERFAS Energy

- 279-000-000 SERFAS Energy Console
- 279-000-010 SERFAS Energy Footswitch
- 278-400-000 SERFAS Energy Probe Bender