

## Cannulas

### Color-Coded Cannulas

Use with Formula Series  
Cutters and Burs

**3.5 mm**

(Locks in all 375 Series Small-Joint Cutters and Burs—*except 3.0 mm Hooded Abrasion Bur*)

**4.5 mm**

(Fits 3.5 mm and 4.0 mm Cutters)

**5.0 mm**

(Fits all Cutters and Blue 5.0 mm Burs)

**5.6 mm**

(Fits all Cutters and Red 4.0 mm Burs)

**6.0 mm**

(Fits all Cutters and Burs)



**272-707-000** Cutter Accessory Tray  
Accommodates all cutter and bur hardware

**275-722-000** 4.5 mm Cutter Cannula  
» Fits all Cutters up to 4.0 mm

**275-722-300** 4.5 mm Cutter Cannula with Irrigation Holes at Distal Tip  
» Fits all Cutters up to 4.0 mm

**275-722-100** 4.5 mm Trocar

**275-722-200** 4.5 mm Obturator

**275-720-000** 5.0 mm-Bur Cannula  
» Fits all Cutters  
» Fits Burs up to 5.0 mm

**275-720-100** 5.0 mm-Bur Trocar

**275-720-200** 5.0 mm-Bur Obturator

**275-719-000** 5.6 mm Cutter Cannula  
» Fits all Cutters  
» Fits Burs up to 4.0 mm

**275-719-300** 5.6 mm Cutter Cannula with Irrigation Holes at Distal Tip  
» Fits all Cutters  
» Fits Burs up to 4.0 mm

**275-719-100** 5.6 mm Trocar

**275-719-200** 5.6 mm Obturator

**275-721-000** 6.0 mm-Bur Cannula  
» Fits all Cutters and Burs

**275-721-100** 6.0 mm-Bur Trocar

**275-721-200** 6.0 mm-Bur Obturator

**375-035-302** 3.5 mm Formula Small-Joint Cutter Cannula  
» Fits all Small-Joint Cutters  
» Fits all Burs, *except 3.0 mm Hooded Abrasion Bur*

**375-035-100** 3.5 mm Formula Trocar

**375-035-200** 3.5 mm Formula Obturator

**375-023-027** 2.7 mm Scope Bridge  
» For use with 2.7 mm Regular Length Speed-Lock Autoclavable Arthroscopes

### Non Color-Coded Cannulas

**270-713-300** 5.0 mm Cutter Cannula with Holes

**270-713-100** 5.0 mm Trocar

**270-713-200** 5.0 mm Obturator