

SERFAS Energy Suction Probes



Super 90-S
SERFAS Energy



90-S Max
SERFAS Energy



90-S
SERFAS Energy



Aardvark
SERFAS Energy



50-S (Angled)
SERFAS Energy



30-S (Angled)
SERFAS Energy

SERFAS Energy Suction Probe (135 mm)

All probes come with integrated cable and hand controls

- | | |
|--------------------|---------------------------------------------------------------------------------------------------------------|
| 279-351-300 | 3.5 mm Super 90-S SERFAS Energy Suction Probe
Most complete 90° suction probe offered |
| 279-401-100 | 4.0 mm 90-S MAX SERFAS Energy Suction Probe
Largest 90° suction probe for increased suction capabilities |
| 279-351-100 | 3.5 mm 90-S SERFAS Energy Suction Probe
Aggressive, all-purpose suction probe |
| 279-251-101 | 2.5 mm Aardvark SERFAS Energy Suction Probe
Bendable suction probe for tight spaces |
| 279-351-250 | 3.5 mm 50-S SERFAS Energy Suction Probe
50° all-purpose suction probe for use in posterior knee procedures |
| 279-351-230 | 3.5 mm 30-S SERFAS Energy Suction Probe
30° all-purpose suction probe |



SERFAS Energy Probes



90°
SERFAS Energy



Direct
SERFAS Energy



Lat
SERFAS Energy



Contour
SERFAS Energy



Hook
SERFAS Energy

Standard SERFAS Energy Probes (135 mm)

All probes come with integrated cable and hand controls

- | | |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 279-350-101 | 3.5 mm 90° SERFAS Energy Probe
An effective all-purpose probe; the 90° electrode orientation at the distal tip provides increased visibility during use |
| 279-350-401 | 3.5 mm Direct SERFAS Energy Probe
Provides excellent access for posterior areas of the knee and shoulder |
| 279-350-201 | 3.5 mm Lat SERFAS Energy Probe
Effective for ablating larger areas of soft tissue |
| 279-350-301 | 3.5 mm Contour SERFAS Energy Probe
An effective all-purpose probe |
| 279-350-501 | 3.5 mm Hook SERFAS Energy Probe
Specifically designed for a lateral release and precision cutting |

Hip Length SERFAS Energy Probes (180 mm)

All probes come with integrated cable and hand controls

- | | |
|--------------------|---------------|
| 279-351-103 | 90-S XL Probe |
| 279-351-253 | 50-S XL Probe |
| 279-350-502 | Hook XL Probe |