

# *157 Exam Light*



### So versatile

It can be a floor lamp—with casters or a stationary base. Or you can mount it on a wall or cabinet.

It's flexible, too. The neck twists and bends. Goes from a height of two to seven feet. Puts light wherever you need it. And it stays put. There's no drifting with this light. A conical reflector inside the head focuses a bright beam on the task at hand.

### Cool to the touch

Ritter's special dual ventilation design circulates air between the inner reflector and the outer shade. No burns or hot surprises from accidental contact.

### Long lasting

Everything about the 157 is built for maximum durability: the tough Lexan impact-resistant shade, the heavy-duty vinyl-covered flexarm, the advanced design aluminum heat sink that dissipates heat and contributes to longer switch and socket life, and a grounded 12', three-conductor cord set with molded plug for safety.

The 157 Exam Light is shown on the front of this page as a floor lamp with casters and in the wall-mounted and casework-mounted versions. Also pictured is the 129 Side Chair, 1000 Series Casework and 1705 Mobile Treatment Cart.

### Specifications

#### Lamp head

Uses standard incandescent bulb 100 watt capacity  
Shade construction Heavy-duty socket w/ aluminum heat sink  
Heavy-duty on/off switch  
Shade dimensions

#### Lamp neck

Dimensions

Neck material High carbon steel  
Neck covering Heavy duty vinyl

Color Black and chrome

Black Lexan polycarbonate

4½" (11.4 cm) diameter, 5¾" (14.6 cm) length

¾" (1.9 cm) diameter, 30" (76.2 cm) length  
High carbon steel  
Heavy duty vinyl

#### Extension tube

Adjustable From 31" (78.7 cm) to 59" (149.9 cm)  
Tube material Chrome plated carbon steel  
Inner tube size ⅝" (1.6 cm)  
Outer tube size ⅞" (2.2 cm)

#### Caster base (Model 157-001)

Construction Cast iron  
Size 18" (45.7 cm) diameter, 4" (10.2 cm) height  
Weight 17 lbs. (7.7 kg)  
Caster material Strong injection molded engineering grade nylon  
Caster wheel diameter 2⅝" (5.6 cm)

#### Standard base (Model 157-002)

Construction Cast iron  
Size 18" (45.7 cm) diameter, 3" (7.6 cm) height

Weight 16 lbs. (7.3 kg)  
Molded rubber floor cushions  
UL Listed, CSA Listed

### Options

Casework and/or wall mounting brackets

