



Typical Drop Ceiling Install Galaxy Gen 3 Dual 4-Pod Light



Touch Control Pad
5 Level Dimming

Galaxy Gen3 Dual 4-Pod Specifications

Input Voltage	100-240 VAC
Input Current	1.0 Amps @ 120 VAC
Input Wattage	120 Watts
Light Output per Head	55,000 Lux @ 1 Meter
Total Light Output	110,000 Lux @ 1 Meter
Moment Load	130 ft-lbs
Weight	80 Lbs.
Min. Ceiling Hgt.	96" @ 80" Pivot Hgt
Max. Ceiling Hgt.	126" @ 80" Pivot Hgt

Notes:

1. Customer is responsible for consulting with structural engineer to determine design of support structure.
2. Electrical connections to be completed by licensed electrical contractor in compliance with all local codes.

Mounting Info

