



Typical Drop Ceiling Install Galaxy Gen 3 Single 4-Pod Light




Touch Control Pad
5 Level Dimming


Galaxy Gen3 4-Pod Single Specifications

Input Voltage	100-240 VAC
Input Current	0.5 Amps @ 120 VAC
Input Wattage	60 Watts
Light Output	55,000 Lux @ 1 Meter
Moment Load	80 ft-lbs
Weight	51 Lbs.
Min. Ceiling Height	93" @ 80" Pivot Hgt.
Max. Ceiling Height	126" @ 80" Pivot Hgt.

Mounting Info



Front View



Notes:

1. Customer is responsible for consulting with structural engineer to determine design of support structure.
2. Electrical connections to be completed by licensed electrical contractor in compliance with all local codes.

